

AccuPrep® Universal RNA Extraction Kit (K-3140G, K-3141G)**I Before You Begin**

- 1) Add 10 μ l of β -mercaptoethanol per 1 ml RB Buffer.
- 2) Prepare additional ethanol (80% and 100%) that is not included.
- 3) Add correct amount of absolute ethanol to RWB1 Buffer and RWB2 Buffer, respectively (see bottle label).

II Cultured Cell Collection

- 1) Suspension cell culture:
Harvest cultured cells (10^4 - 10^6 cells) by centrifugation at 300 x g for 5 min to pellet cells. Discard supernatant carefully without disturbing a cell pellet and go to step III.
- 2) Monolayer cell culture: There are 2 different ways to collect cells grown in monolayer culture.
 - a. Direct cell harvesting on the culture dishes: Completely remove cell culture medium and go to step III.
 - b. Cell harvesting with trypsin: Remove cell culture medium and wash the cell monolayer with DPBS. Add 0.1-0.25% trypsin to the washed cell monolayer to detach the cells. When the cells are detached, add cell culture medium to inactivate trypsin. Transfer the cells into a RNase-free tube and centrifuge at 300 x g for 5 min. Discard supernatant carefully and go to step III.

III RNA Extraction from Cultured Cell

- 1) Add 400 μ l of RB Buffer to the cell pellet and mix by vortex mixer.
- 2) Add 300 μ l of 80% ethanol and mix immediately by using pipette.
- 3) Transfer the sample to a Binding column fitted in a 2 ml collection tube.
- 4) Close the lid and centrifuge at $\geq 14,000$ rpm for 20 sec.
- 5) Discard the flow-through from the collection tube and reuse the collection tube.
- 6) Add 700 μ l of RWB1 Buffer without wetting the rim, close the tube, and centrifuge at 14,000 rpm for 20 sec.
- 7) Discard the solution from the collection tube and reuse the collection tube.
- 8) Add 500 μ l of RWB2 Buffer without wetting the rim, close the tube, and centrifuge at 14,000 rpm for 20 sec.
- 9) Discard the solution from the collection tube and reuse the collection tube.
- 10) Add 500 μ l of RWB2 Buffer without wetting the rim, close the tube, and centrifuge at 14,000 rpm for 2 min.
- 11) Discard the solution from the collection tube and reuse the collection tube.

AccuPrep® Universal RNA Extraction Kit (K-3140G, K-3141G)

- 12) Centrifuge once more at **14,000 rpm** for **1 min** to completely remove ethanol, and check that there is no droplet clinging to the bottom of binding column tube.
- 13) Transfer the Binding column tube to a new 1.5 ml tube (supplied with a kit) for elution, add **50-200 µl** of **ER Buffer** onto Binding column tube, and wait for at least 1 min at RT (15-25°C).
- 14) Centrifuge at **10,000 rpm** for **1 min** to elute.

IV RNA Extraction from Plant Tissue

- 1) Grind Sample (up to 100 mg) using liquid nitrogen. Transfer the tissue powder to an appropriately sized tube and add **500 µl** of **RB Buffer** to a maximum of 100 mg tissue powder and vortex vigorously.
- 2) Incubate at **60°C** for **1-3 min**.
- 3) Centrifuge at full speed for **2 min**.
- 4) Transfer the supernatant to a new 1.5 ml tube.
- 5) Add **absolute ethanol (not provided)** by 0.5 volume of sample to the sample and mix immediately by using pipette.
- 6) Go to **step 3** of "RNA Extraction from Cultured Cell" in page 1 and follow the instructions accordingly.

V RNA Extraction from Animal Tissue

- 1) Homogenize the sample (20~30 mg) with a homogenizer, place them in a new 1.5 ml tube, and add **500 µl** of **RB Buffer**.
- 2) Centrifuge the lysate for 3 min at full speed, and transfer the supernatant to a new 1.5 ml tube.
- 3) Add **200 µl** of **absolute ethanol (not provided)** and mix immediately by using pipette.
- 4) Go to **step 3** of "RNA Extraction from Cultured Cell" in page 1 and follow the instructions accordingly.

※ For more information, please visit www.bioneer.com and refer to the User Guide of this kit.