

[Cat. No.] K-7110, K-7120, K-7130

## Introduction

AccuRapid™ Cloning Kit is designed for directional cloning of up to 3 PCR products into linearized vectors without the need for restriction enzymes, ligase, or phosphatase. This kit contains an optimized Enzyme Mix, which fuses PCR products and linearized vectors by annealing a 18-21 bp complementary sequence located at their ends. The 18-21 bp homologue sequence is added to the extension primers for the amplification of the inserts.




Figure 1. Workflow for cloning

## Features & Benefits

- **Convenient:** Directional cloning of inserts into a vector is possible without the use of restriction enzymes, a ligase, or a phosphatase.
- **Rapid:** Clones accurately in 30 min.
- **Accurate:** Cloning of up to 3 inserts is possible through precise design of extension primers.

## Components

Components	K-7110
AccuRapid™ Enzyme Mix	45 µl
2 kb pBHA Control Vector (25 ng/µl)	3 µl
750 bp Control Insert (50 ng/µl)	5 µl

\* **Note:** For research use only. Not for use in diagnostic or therapeutic procedures.

## Storage

- Store at -20°C.

## Online Resources



Korean



English

Visit our [product page](#) for additional information and protocols.

## Ordering Information

Description	Cat. No.
10 rxns x 1 ea	K-7110
AccuRapid™ Cloning Kit	10 rxns x 2 ea
	K-7120
	10 rxns x 5 ea
	K-7130

## Notice

BIONEER corporation reserves the right to make corrections, modifications, improvements and other changes to its products, services, specifications or product descriptions at any time without notice.

## Explanation of Symbols



Batch Code



Catalog Number



Caution



Consult Instructions For Use



Contains Sufficient for <n> tests



Manufacturer



Research Use Only









Temperature Limitation



Use-by Date

## Experimental Procedures

Steps		Procedure Details																								
1	 <p><b>Preparation of linearized vector</b></p>	<p>1. Prepare a linearized vector through PCR or treatment with restriction enzymes.</p>																								
2	 <p><b>Extension primer design</b></p>	 <p>2. Design extension primers which contain a 18-21 bp complementary sequence at the end of a linearized vector.</p> <p><b>* Note:</b> Refer to the User's guide for additional information.</p>																								
3	 <p><b>Insert amplification</b></p>	<p>3. Perform the PCR for the amplification of insert.</p> <p>4. Purify the samples using a PCR/Gel purification kit.</p>																								
4	 <p><b>Cloning</b></p>	<p>5. Prepare the reaction mixture.</p> <ul style="list-style-type: none"> <li>Preparation of reaction mixture</li> </ul> <table border="1"> <thead> <tr> <th>Components</th> <th>Negative</th> <th>Positive</th> <th>Sample</th> </tr> </thead> <tbody> <tr> <td>Linearized Vector (25-50 ng)</td> <td>1 µl</td> <td>1 µl</td> <td>1 µl</td> </tr> <tr> <td>Purified PCR products (70-150 ng)</td> <td>-</td> <td>1.5 µl</td> <td>Variable†</td> </tr> <tr> <td>Distilled water</td> <td>5 µl</td> <td>3.5 µl</td> <td>Variable</td> </tr> <tr> <td>AccuRapid™ Enzyme Mix</td> <td colspan="3">4 µl</td> </tr> <tr> <td>Total volume</td> <td colspan="3">10 µl</td> </tr> </tbody> </table> <p><b>* Note:</b> After cloning the reaction mixture should be stored either on ice or at -20°C until transformation.  † Purified PCR products need to be 70-150 ng/rxn totally, and 5 µl maximum in volume.  ex) 1 fragment: 5 µl, 2 fragments: 2.5 µl x 2, 3 fragments: 1.7 µl x 3</p> <p>6. Mix the reaction mixture by tapping the tube gently and briefly spin down.</p> <p>7. Incubate the reaction mixture for 30 min at 50°C.</p>	Components	Negative	Positive	Sample	Linearized Vector (25-50 ng)	1 µl	1 µl	1 µl	Purified PCR products (70-150 ng)	-	1.5 µl	Variable†	Distilled water	5 µl	3.5 µl	Variable	AccuRapid™ Enzyme Mix	4 µl			Total volume	10 µl		
Components	Negative	Positive	Sample																							
Linearized Vector (25-50 ng)	1 µl	1 µl	1 µl																							
Purified PCR products (70-150 ng)	-	1.5 µl	Variable†																							
Distilled water	5 µl	3.5 µl	Variable																							
AccuRapid™ Enzyme Mix	4 µl																									
Total volume	10 µl																									
5	 <p><b>Transformation</b></p>	<p>8. Add the reaction mixture to the <i>E. coli</i> cells to be transformed.</p>																								