AccuRapid™ Cloning Kit

Cat. No. K-7110
K-7120
K-7130

AccuRapid™ Cloning Kit

User Guide

K-7110
 **10**

K-7120
 **20**

K-7130
 **50**

Version No.: 2 (2022-02-23)

Please read all the information in booklet before using the unit

BIONEER Corporation

**8-11, Munpyeongseo-ro, Daedeok-gu, Daejeon
34302, Republic of Korea**

Tel: 1588-9788

Email: sales@bioneer.co.kr

www.bioneer.com

Intended Use

AccuRapid[™] Cloning Kit is developed and supplied for research purposes only. Certain applications possible with this kit may require special approval by appropriate local and/or national regulatory authorities in the country of use.

Safety Warning and Precaution

Wear appropriate protection when handling any irritant or harmful reagents. The use of a laboratory coat, protective gloves and safety goggles are highly recommended. For more information, please consult the appropriate Material Safety Data Sheet (MSDS).

Warranty and Liability

All BIONEER products undergo extensive Quality Control testing and validation. BIONEER guarantees quality during the warranty period as specified, when following the appropriate protocol as supplied with the product. It is the responsibility of the purchaser to determine the suitability of the product for its particular use. Liability is conditional upon the customer providing full details of the problem to BIONEER within 30 days.

Quality Management System ISO 9001 Certified

Every aspect of our quality management system from product development, production to quality assurance and supplier qualification meets the world-class standards.

Trademark

AccuRapid[™] is a trademark of BIONEER Corporation.

Copyright

Copyright 2022 BIONEER Corporation. All Rights Reserved.

Notice

BIONEER Corporation reserves the right to make corrections, modifications, improvements and other changes to its products, services, specifications or product descriptions at any time without notice.

Contents

Product Information	1
Components	1
Storage	1
Specifications	1
Introduction	2
Product Description	2
Experimental Procedures	3
Preparation of Linearized Vector	3
Primer Design	4
Insert Amplification	5
Cloning Reaction	7
Transformation	8
Results Analysis	9
Troubleshooting	10
Ordering Information	12
Related Products	12
Explanation of Symbols	13

Product Information

Components

Components	Amount		
	K-7110 (10 rxn)	K-7120 (20 rxn)	K-7130 (50 rxn)
<i>AccuRapid</i> [™] Enzyme Mix	45 µl	45 µl x 2 ea	45 µl x 5 ea
2 kb pBHA Control Vector (Amp ^R , 25 ng/µl)	3 µl	3 µl x 2 ea	3 µl x 5 ea
750 bp Control Insert (50 ng/µl)	5 µl	5 µl x 2 ea	5 µl x 5 ea

Storage

All components of *AccuRapid*[™] Cloning Kit should be stored at -20°C.

Specifications

	<i>AccuRapid</i> [™] Cloning Kit
Reaction time	30 min
Reaction temperature	50°C
Need for PCR/Gel purification	Recommendation

Introduction

Product Description

BIONEER's *AccuRapid*™ Cloning Kit allows an accurate and a rapid cloning of 1-3 inserts (PCR products) into a linearized vector. This product connects the linearized vector and the PCR-amplified inserts by recognizing 18-21 bp complementary sequences at each end of the linearized vector and the amplified inserts. Unlike the conventional method which performs restriction enzyme digestion and ligation, this product does not require a treatment of restriction enzyme with inserts. Instead, inserts to be cloned can be prepared rapidly through PCR and accurately positioned within vector for cloning. The PCR primers can be easily designed by adding 18-21 bp complementary sequences found at both ends of the vector to the 5' end of the primer.


Figure 1. Schematic protocol for *AccuRapid*™ Cloning Kit.

Experimental Procedures

Preparation of Linearized Vector

In order to use the *AccuRapid*[™] Cloning Kit, it is important to completely linearize the vector using restriction enzyme digestion or PCR.

Experiments with completely linearized vectors will increase the efficiency of cloning, forming many colonies with inserts as well as reducing colony formation of circular vectors without inserts.

*** Note:** In order to increase cloning efficiency, it is recommended to check the linearized vector through electrophoresis after restriction enzyme treatment and proceed to the next experiment.


1. The vector to be used for cloning is extracted from *E. coli* using the *AccuPrep*[®] Nano-Plus Plasmid Mini Extraction Kit (Cat. No. K-3111).
2. Linearize the vector by restriction enzyme treatment or PCR.
3. After the linearized vector is clearly separated from the circular vector through electrophoresis, the linearized vector is purified using *AccuPrep*[®] PCR/Gel Purification Kit (Cat. No. K-3038).
4. Prepare DNA at a concentration of 25-50 ng/μl.

Primer Design

Primers to amplify the insert should be designed to contain 18-21 bp complementary sequences at both ends of the linearized vector at the 5' end of the primers.

* You can easily design extension primers for insert amplification on the *AccuRapid™* Cloning Kit product homepage.

* BIONEER's oligo synthesis service allows you to receive high-quality primers quickly.

ex) When cloning into *Nde*I (CATATG), *Xho*I (CTCGAG) enzyme site.

Vector sequences:

TCTCGATCCCGCGAAATTAATACGACTCACTATAGGGGAATTGTGAGCGGATAACAATTCCCC
 TCTAGAAATAATTTTGTCTAACTTTAAGAAGGAGATATACCATGGGCAGCAGCCATCATCATCA
 TCATCACAGCAGCGGCCTG**GTGCCGCGCGGCAGCCATATG**(*Nde*I)GCTAGCATGACTGGTG
 GACAGCAAATGGGTCCGCGGATCCGAATTCGAGCTCCGTCGACAAGCTTGC GGCCGCA(*Xho*
I)CTCGAGCACCACCACCACCACCACTGAGATCCGGCTGCTAACAAAGCCCGAAAGGAAGC
 TGAGTTGGCTGCTGCCACCGCTGAGCAATAACTAGCATAACCCCTTGGGGCCTCTAAACGG
 GTCCTGAGGGGTTTTTGTCTGAAAGGAGGA ACTATATCCGGAT

3' end vector sequences (Forward primer)


GTGCCGCGCGGCAGCCATATGNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNN
 5' end insert sequences

5' end vector sequences (Reverse primer)

GTGGTGGTGGTGGTGCTCGAGNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNN
 3' end insert sequences (Reverse complement)

Insert Amplification

After primers design, PCR is performed with the synthesized extension primers to amplify the insert (Using BIONEER's PCR equipment, *AllInOneCycler*[™], you can amplify rapidly). After that, perform PCR purification or electrophoresis followed by gel purification to obtain a purified insert.

1. Using an existing template

- 1) Use the *AllInOneCycler*[™] PCR system (Cat. No. A-2041) and *AccuPower*[®] *ProFi Taq* PCR PreMix (Cat. No. K-2631) to amplify the insert from the template DNA.
- 2) Add template DNA, primers, and nuclease-free water into *AccuPower*[®] *ProFi Taq* PCR PreMix tubes to make a total volume of 20 μ l or 50 μ l. Do not calculate the dried pellet.

- Preparation of reaction mixture

Components	20 μ l reaction	50 μ l reaction
Template DNA (1-500 ng)	Variable (1-10 μ l)	Variable (1-25 μ l)
Forward primer (10 pmol/ μ l)	0.5-2 μ l	1-5 μ l
Reverse primer (10 pmol/ μ l)	0.5-2 μ l	1-5 μ l
Nuclease-free water	Variable	Variable
Total volume	20 μ l	50 μ l

- 3) Dissolve the vacuum-dried pellet by vortexing, and briefly spin down.
- 4) Perform the reaction under the following conditions.

Step	Temperature	Time	Cycles
Pre-denaturation	95°C	5 min	1 cycle
Denaturation	95°C	15~20 sec	30 cycles
Annealing	45~65°C	15~30 sec	
Extension	72°C	1 min/kb	
Final extension	72°C	3~5 min	1 cycle

5) After electrophoresis with the amplified insert, gel eluted with *AccuPrep*® PCR/Gel Purification Kit (Cat. No. K-3038).

2. When using AccuGeneBlock Service

When synthesizing an insert using that AccuGeneBlock service, it is recommended to add vector complementary sequences (18-21 bp) at both ends.

- 1) Access the BIONEER website (<http://www.bioneer.com>).
- 2) Click the banner in order of Gene/mRNA Synthesis and Gene Synthesis Custom Order.
- 3) Select AccuGeneBlock among service types.
- 4) Enter the gene name and sequences to be synthesized.
- 5) After selecting the Scale and End type, confirm the quotation.
- 6) Click Preview Invoice to confirm your order history and place your order.

(To complete your order, you must click Place Order in the shopping cart.)

Cloning Reaction

For the cloning reaction, mix the reagents as described below, spin down, and incubate at 50°C for 30 min using an incubator or thermal cycler. After the reaction, keep it in ice or at -20°C until use for transformation.

1. Prepare the reaction mixture under the following tables.

Components	Negative	Positive	Sample
Linearized Vector (25-50 ng)	1 µl (Sample Vector)	1 µl (2 kb pBHA Control Vector, Amp ^R)	1 µl (Sample Vector)
Purified PCR products (70-150 ng)	-	1.5 µl (750 bp Control Insert)	Variable*
Nuclease-free water	5 µl	3.5 µl	Variable
<i>AccuRapid</i> TM Enzyme Mix		4 µl	
Total volume		10 µl	

* Purified PCR products need to be 70-150 ng/rxn totally, and 5 µl maximum in volume.

ex) 1 fragment: 5 µl, 2 fragments: 2.5 µl x 2, 3 fragments: 1.7 µl x 3

2. Mix the reaction mixture by tapping the tube gently and briefly spin down.

3. Incubate the reaction mixture at 50°C for 30 min.

4. After the reaction, reaction mixture should be stored either on ice or at -20°C until transformation.

Transformation

1. Take competent cells from deep freezer and slowly thaw on ice.
2. Add 10 µl of cloning reaction mixture into 50-100 µl of competent cells and incubate on ice for 20 min.
3. Heat shock the cells at 42°C for 90 sec. Then, immediately place tubes on ice and incubate 2 min.
4. Add 1 ml of LB broth (not included antibiotics) and incubate at 37°C shaking incubator (200 rpm) for 1 hr.
5. After the incubation, centrifuge at 3,000 rpm for 5 min, remove the supernatant, add 100 µl of LB broth, and gently resuspend the cell pellet by tapping.
6. Transfer to LB agar plates (including antibiotics) and spread completely.
7. Incubate overnight at 37°C in an incubator.

Results Analysis

1. After transformation, check colonies in the cultured plate.
 - * **Note:** Colony numbers can be different depending on the size and number of inserts, conditions of competent cells, transformation conditions, etc.
2. Single colonies are inoculated into LB broth containing antibiotics and incubated at 37°C shaking incubator (200 rpm) for overnight. Then, extract the plasmid using *AccuPrep*[®] Nano-Plus Plasmid Mini Extraction Kit (Cat. No. K-3111) and check whether your gene of interest has been successfully cloned through restriction enzyme treatment, PCR, sequencing, etc.

Troubleshooting

1. Little or no colonies produced

Checking way	Cause	Solution
Colonies are normally formed in positive control.	Incorrect design of extension primers for insert amplification	Confirm vector sequences whether containing complementary sequences and design extension primers again.
	Incorrect vector sequences	Proceed the experiment after confirming the correct vector sequences through sequencing.
	The amount of vector or insert used in the reaction is small.	Even if the band is well identified on electrophoresis after restriction enzyme treatment or PCR, the insert may be contaminated or lost during gel purification or PCR purification. Therefore, after gel extraction or PCR purification, a small amount of insert or vector should be run through electrophoresis to confirm the appropriate concentration for the reaction.
Colonies are not formed in positive control	Decreased transformation efficiency of competent cells	Check the expiration date and efficiency of competent cells.
	More than 30 min of cloning reaction	Since the cloning efficiency may be reduced, the reaction proceeds for 30 min. If

	transformation is not possible immediately, store the reaction at -20°C until use.
Using large amounts of <i>AccuRapid</i> [™] Enzyme Mix	If you add <i>AccuRapid</i> [™] Enzyme Mix more than the appropriate amount, the cloning efficiency may decrease, so use only 4 µl.

2. A large number of colonies are generated, but the exact cloning product is not confirmed.

Checking way	Cause	Solution
Many colonies are formed in negative control.	Non-linearized vector	After restriction enzyme treatment or PCR, check the linearized vector through electrophoresis and proceed.
When PCR is performed with candidates, the correct band is not confirmed.	Contamination of insert (PCR fragment)	The insert may be contaminated during PCR or purification, so a small amount of insert should be reconfirmed by electrophoresis and used.
When confirming the size or sequences of the PCR product, sequences that is different in size or similar to the vector is found.	The same or similar sequences to the 18-21 bp complementary sequences of the primers also exist at the other sites within the vector	Proceed with a different cloning design, or use a ligation method after restriction enzyme treatment.

Ordering Information

Description		Cat. No
AccuRapid™ Cloning Kit	10 rxns	K-7110
	20 rxns	K-7120
	50 rxns	K-7130

Related Products

Description	Cat. No
Oligo Synthesis (Primer) Service	S-1001
AccuGeneBlock Service	S-2043
Gene Synthesis Service	S-2041
AccuPower® ProFi Taq PCR PreMix	K-2631
AccuPrep® Nano-Plus Plasmid Mini Extraction Kit	K-3111
AccuPrep® PCR/Gel Purification Kit	K-3038
AllInOneCycler™ PCR System	A-2041
Agarose, 100 g	C-9100
50X TAE ,500 ml	C-9004
Agaro-Power™ System	A-7020
DUALED Blue/White Transilluminator	A-6020
GreenStar™ Nucleic Acid Staining Solution I	C-9036

Explanation of Symbols

 <p>Batch Code</p>	 <p>Catalog Number</p>	 <p>Caution</p>	 <p>Consult Instructions For Use</p>
 <p>Contains Sufficient for <n> tests</p>	 <p>Manufacturer</p>	 <p>Research Use Only</p>	 <p>Temperature Limitation</p>
 <p>Use-by Date</p>			

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BIONEER Corporation - HQ
BIONEER Inc. - USA Branch
BIONEER Corp. - European Branch