

[Cat. No.] **K-6714, K-6715**

Introduction

AccuPower® Dual-HotStart™ RT-PCR PreMix (with UDG) is a one-step RT-PCR product that applied Pyro-HotStart RT technology and HotStart PCR technology to improve problems of non-specific reverse transcription reaction. It can effectively synthesize cDNA with small amounts of template RNA and complex secondary RNA structures. Moreover, application of uracil DNA glycosylase (UDG) system prevents carryover contaminations. UDG hydrolyzes the N-glycosidic bonds linking uracil and deoxyribose, breaking down the leftover templates inserted with uracil. Before the PCR, contaminants are removed by activating enzymes at 37°C for 2 min and enzymes are deactivated at the higher temperature during the cycles.




Figure 1. The 1st HotStart reaction. BIONEER's *RocketScript*™ Reverse Transcriptase is completely inhibited by pyrophosphate at temperatures below 50°C. However, *RocketScript*™ RTase becomes fully active at temperatures above 50°C via pyrophosphate hydrolysis with a thermostable pyrophosphatase. This prevents the formation of mis-primed products and primer dimers during the reaction set up process resulting in improved specificity of cDNA synthesis.




Figure 2. The 2nd HotStart reaction. BIONEER's HotStart *Taq* DNA Polymerase provides superior priming accuracy and specificity that cannot be achieved with other enzymes.

Applications

- Low copy viral/bacterial pathogen load determination in an earlier stage
- Low copy mRNA amplification
- Low copy target RNA quantification
- RNA amplification for microarray and NGS

Features & Benefits

- Carryover contamination prevention: Minimized false positives caused by a carryover contamination through application of uracil DNA glycosylase system.
- Specificity: Our premier *Dual-HotStart*™ RT-PCR reaction provides accurate results of target gene amplification by using Pyro-HotStart RT reaction and HotStart PCR.

Composition

Composition	Concentration
Uracil DNA glycosylase	1 U
<i>RocketScript</i> ™ Reverse Transcriptase	200 U

5X Reaction buffer	1X
DTT	0.25 mM
dNTPs with dUTP	1.2 mM
RNase inhibitor	1 U
HotStart <i>Taq</i> DNA Polymerase	1 U
Stabilizer and tracking dye	1X

Specifications

HotStart <i>Taq</i> DNA Polymerase	
5' to 3' exonuclease activity	Yes
3' to 5' exonuclease activity	No
3'-A overhang	Yes
Fragment size	Up to 3 kb

Enzyme Inactivation

UDG is inactivated by heating at 95°C for 5 min.

Storage

Store at -20°C. If stored in the recommended temperature, this product will be stable until the expiration date printed out on the label.

Precautions

UDG activities can be remained after finishing reactions, if it is kept on below 50°C. Therefore, reaction mixture is recommended to freeze immediately after the reaction.

Online Resources



Korean



English

Visit our [product page](#) for additional information and protocols

Ordering Information

Description			Cat. No.
0.2 ml thin-wall 8-tube	96 tubes	20 µl/rxn	K-6714
strips with attached cap		50 µl/rxn	K-6715




Notice

BIONEER corporation reserves the right to make corrections, modifications, improvements and other changes to its products, services, specifications or product descriptions at any time without notice.

Explanation of Symbols

Batch Code	Biological Risks	Catalog Number	Caution
Consult Instructions For Use	Contains Sufficient for <n> tests	Do not Re-use	Manufacturer
Research Use Only	Temperature Limitation	Use-by Date	

Experimental Procedures

Steps		Procedure Details																																
1	 Preparation of reaction mixture	<p>1. Add template RNA, primers and nuclease-free water into <i>AccuPower® Dual-HotStart™ RT-PCR PreMix (with UDG)</i> tubes to make a total volume of 20 µl or 50 µl. Do not calculate the dried pellet.</p> <ul style="list-style-type: none"> Amount of template RNA and primers <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th colspan="2">Components</th> <th>20 µl reaction</th> <th>50 µl reaction</th> </tr> </thead> <tbody> <tr> <td rowspan="2">Template RNA</td> <td>Total RNA</td> <td>10 pg-5 µg</td> <td>10 pg-5 µg</td> </tr> <tr> <td>Poly(A) RNA</td> <td>10 pg-5 µg</td> <td>10 pg-5 µg</td> </tr> <tr> <td>Primers</td> <td>Gene specific primer</td> <td>10-30 pmol</td> <td>10-30 pmol</td> </tr> </tbody> </table> <p>2. Dissolve the vacuum-dried pellet by pipetting or vortexing, and briefly spin down.</p>	Components		20 µl reaction	50 µl reaction	Template RNA	Total RNA	10 pg-5 µg	10 pg-5 µg	Poly(A) RNA	10 pg-5 µg	10 pg-5 µg	Primers	Gene specific primer	10-30 pmol	10-30 pmol																	
Components		20 µl reaction	50 µl reaction																															
Template RNA	Total RNA	10 pg-5 µg	10 pg-5 µg																															
	Poly(A) RNA	10 pg-5 µg	10 pg-5 µg																															
Primers	Gene specific primer	10-30 pmol	10-30 pmol																															
2	 RT-PCR	<p>3. Perform the reaction under the following conditions.</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Step</th> <th>Temperature</th> <th>Time</th> <th>Cycles</th> </tr> </thead> <tbody> <tr> <td>UDG activation</td> <td>37°C</td> <td>2 min</td> <td>1 cycle</td> </tr> <tr> <td>cDNA synthesis</td> <td>42-70°C</td> <td>10-60 min</td> <td>1 cycle</td> </tr> <tr> <td>Pre-denaturation</td> <td>95°C</td> <td>10 min</td> <td>1 cycle</td> </tr> <tr> <td>Denaturation</td> <td>95°C</td> <td>10-30 sec</td> <td></td> </tr> <tr> <td>Annealing</td> <td>50-65°C</td> <td>10-30 sec</td> <td>30-40 cycles</td> </tr> <tr> <td>Extension</td> <td>72°C</td> <td>1 kb/min</td> <td></td> </tr> <tr> <td>Final extension</td> <td>72°C</td> <td>5 min</td> <td>1 cycle</td> </tr> </tbody> </table> <p>* Note: Reaction temperature should be optimized according to T_m value of primers.</p>	Step	Temperature	Time	Cycles	UDG activation	37°C	2 min	1 cycle	cDNA synthesis	42-70°C	10-60 min	1 cycle	Pre-denaturation	95°C	10 min	1 cycle	Denaturation	95°C	10-30 sec		Annealing	50-65°C	10-30 sec	30-40 cycles	Extension	72°C	1 kb/min		Final extension	72°C	5 min	1 cycle
Step	Temperature	Time	Cycles																															
UDG activation	37°C	2 min	1 cycle																															
cDNA synthesis	42-70°C	10-60 min	1 cycle																															
Pre-denaturation	95°C	10 min	1 cycle																															
Denaturation	95°C	10-30 sec																																
Annealing	50-65°C	10-30 sec	30-40 cycles																															
Extension	72°C	1 kb/min																																
Final extension	72°C	5 min	1 cycle																															
3	 Analyze with gel electrophoresis	<p>4. After the reaction, maintain the reaction mixture at 4°C. The samples can be stored at -20°C until use.</p> <p>5. Load samples on agarose gel without adding a loading-dye mixture, and perform gel electrophoresis for analysis.</p>																																