



[Cat. No.] Please refer to the **Ordering Information**

### Introduction

AccuPower® Dual-HotStart™ RT-PCR PreMix is a one-step RT-PCR product that applied Pyro-HotStart RT technology and HotStart PCR technology to improve problems of non-specific reverse transcription reaction. It can effectively synthesize cDNA with small amounts of template RNA and complex secondary RNA structures. Moreover, it can be applied to various types of samples and provides accurate RT-PCR results with high specificity. This product contains vacuum-dried components including *RocketScript*™ Reverse Transcriptase, reaction buffer, DTT, dNTPs, RNase inhibitor, and HotStart *Taq* DNA Polymerase. It simplifies preparation of reverse transcription reaction mixture and PCR mixture by adding template RNA and primers without any extra process.



**Figure 1. The 1<sup>st</sup> HotStart reaction.** BIONEER's *RocketScript*™ Reverse Transcriptase is completely inhibited by pyrophosphate at temperatures below 50°C. However, *RocketScript*™ RTase becomes fully active at temperatures above 50°C via pyrophosphate hydrolysis with a thermostable pyrophosphatase. This prevents the formation of mis-primed products and primer dimers during the reaction set up process resulting in improved specificity of cDNA synthesis.



**Figure 2. The 2<sup>nd</sup> HotStart reaction.** BIONEER's HotStart *Taq* DNA Polymerase provides superior priming accuracy and specificity that cannot be achieved with other enzymes.

### Applications

- Low copy viral/bacterial pathogen load determination in an earlier stage
- Low copy mRNA amplification
- Low copy target RNA quantification
- RNA amplification for microarray and NGS

### Features & Benefits

- **Sensitivity:** Efficient detection even with small amounts of template RNA up to 10 pg of total RNA extracted from cells and tissues.
- **Comprehensive applications:** Blood and soil samples containing various types of PCR inhibitors can be applied due to its excellent reactivity.
- **Specificity:** Our premier *Dual-HotStart*™ RT-PCR reaction provides accurate results of target gene amplification by using Pyro-HotStart RT reaction and HotStart PCR.
- **Ease-of-use:** Reactants are individually packaged in each of the PCR tubes, it allows any user simply perform cDNA synthesis and PCR by adding template RNA and primers.
- **Reproducibility:** Mass production under ISO 9001 quality system allows minimized deviation between lots and reproducible results in replicated tests performed under same conditions and variation.

### Composition

Composition	Concentration
<i>RocketScript</i> ™ Reverse Transcriptase	200 U
5X Reaction buffer	1X
DTT	0.25 mM
dNTPs (dATP, dCTP, dGTP, dTTP)	1.2 mM
RNase inhibitor	1 U
HotStart <i>Taq</i> DNA Polymerase	1 U
Stabilizer and tracking dye	1X

### Specifications

HotStart <i>Taq</i> DNA Polymerase	
5' to 3' exonuclease activity	Yes
3' to 5' exonuclease activity	No
3'-A overhang	Yes
Fragment size	Up to 3 kb

### Storage

Store at -20°C. If stored in the recommended temperature, this product will be stable until the expiration date printed out on the label.

### Online Resources



Korean



English

Visit our **product page** for additional information and protocols

### Ordering Information

Description	Cat. No.
0.2 ml thin-wall 96 tubes	K-6710
8-tube strips	K-6711
with attached cap	K-6712
480 tubes	K-6713




### Notice

BIONEER corporation reserves the right to make corrections, modifications, improvements and other changes to its products, services, specifications or product descriptions at any time without notice.

### Explanation of Symbols

Batch Code	Biological Risks	Catalog Number	Caution
Consult Instructions For Use	Contains Sufficient for <n> tests	Do not Re-use	Manufacturer
Research Use Only	Temperature Limitation	Use-by Date	

**Experimental Procedures**

Steps		Procedure Details																												
1	 <b>Preparation of reaction mixture</b>	<p>1. Add template RNA, primers and nuclease-free water into <i>AccuPower® Dual-HotStart™</i> RT-PCR PreMix tubes to make a total volume of 20 µl or 50 µl. Do not include the dried pellet.</p> <ul style="list-style-type: none"> <li>Amount of template RNA and primers</li> </ul> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th colspan="2">Components</th> <th>20 µl reaction</th> <th>50 µl reaction</th> </tr> </thead> <tbody> <tr> <td rowspan="2">Template RNA</td> <td>Total RNA</td> <td>10 pg-5 µg</td> <td>10 pg-5 µg</td> </tr> <tr> <td>Poly(A) RNA</td> <td>10 pg-5 µg</td> <td>10 pg-5 µg</td> </tr> <tr> <td>Primers</td> <td>Gene specific primer</td> <td>10-30 pmol</td> <td>10-30 pmol</td> </tr> </tbody> </table> <p>2. Dissolve the vacuum-dried pellet by pipetting or vortexing, and briefly spin down.</p>	Components		20 µl reaction	50 µl reaction	Template RNA	Total RNA	10 pg-5 µg	10 pg-5 µg	Poly(A) RNA	10 pg-5 µg	10 pg-5 µg	Primers	Gene specific primer	10-30 pmol	10-30 pmol													
Components		20 µl reaction	50 µl reaction																											
Template RNA	Total RNA	10 pg-5 µg	10 pg-5 µg																											
	Poly(A) RNA	10 pg-5 µg	10 pg-5 µg																											
Primers	Gene specific primer	10-30 pmol	10-30 pmol																											
2	 <b>RT-PCR</b>	<p>3. Perform the reaction under the following conditions.</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Step</th> <th>Temperature</th> <th>Time</th> <th>Cycles</th> </tr> </thead> <tbody> <tr> <td>cDNA synthesis</td> <td>42-70°C</td> <td>10-60 min</td> <td>1 cycle</td> </tr> <tr> <td>Pre-denaturation</td> <td>95°C</td> <td>10 min</td> <td>1 cycle</td> </tr> <tr> <td>Denaturation</td> <td>95°C</td> <td>10-30 sec</td> <td></td> </tr> <tr> <td>Annealing</td> <td>50-65°C</td> <td>10-30 sec</td> <td>30 cycles</td> </tr> <tr> <td>Extension</td> <td>72°C</td> <td>1 kb/min</td> <td></td> </tr> <tr> <td>Final extension</td> <td>72°C</td> <td>5 min</td> <td>1 cycle</td> </tr> </tbody> </table> <p>* <b>Note:</b> Reaction temperature should be optimized T<sub>m</sub> value of primers.</p>	Step	Temperature	Time	Cycles	cDNA synthesis	42-70°C	10-60 min	1 cycle	Pre-denaturation	95°C	10 min	1 cycle	Denaturation	95°C	10-30 sec		Annealing	50-65°C	10-30 sec	30 cycles	Extension	72°C	1 kb/min		Final extension	72°C	5 min	1 cycle
Step	Temperature	Time	Cycles																											
cDNA synthesis	42-70°C	10-60 min	1 cycle																											
Pre-denaturation	95°C	10 min	1 cycle																											
Denaturation	95°C	10-30 sec																												
Annealing	50-65°C	10-30 sec	30 cycles																											
Extension	72°C	1 kb/min																												
Final extension	72°C	5 min	1 cycle																											
3	 <b>Analyze with gel electrophoresis</b>	<p>4. After the reaction, maintain the reaction mixture at 4°C. The samples can be stored at -20°C until use.</p> <p>5. Load samples on agarose gel without adding a loading-dye mixture, and perform gel electrophoresis for analysis.</p>																												