

[Cat. No.] Please refer to the **Ordering Information**

### Introduction

AccuPower® DualStar™ qPCR PreMix enables accurate and rapid quantification of target genes in a wide range of samples through real-time PCR with hydrolysis probe method. By applying BIONEER's patented enzyme-mediate HotStart technology, non-specific reactions are reduced during zero cycles and amplification efficiency is improved even with a trace amount of template DNA. This product contains vacuum-dried components for real-time PCR, except for template DNA, target-specific primers and fluorogenic probes. By just adding template DNA, target-specific primers and probes, reproducible results with high sensitivity and specificity can be obtained. It can be applied in hydrolysis probe-based real-time PCR experiments for the amplification and detection of genomic DNA and cDNA targets, differential gene expression profiling, Single Nucleotide Polymorphism (SNP) analysis, and evaluation of RNAi products.

### Applications

- Gene expression profiling
- Target DNA quantification
- Microbial detection
- Viral/bacterial pathogen load determination
- Evaluation of primer pair performance for probe-based real-time PCR

### Features & Benefits

- Specificity: Minimized non-specific amplification and maximized PCR efficiency by using BIONEER's differentiated PyroHotstart technology.
- Compatibility: Wide choice of real-time PCR instruments for optimal results.
- Comprehensiveness: Effective real-time PCR regardless of gene types, including DNA, cDNA and high GC templates.
- Convenience: Reactants are individually packaged in each of the PCR tubes, it allows any user simply perform real-time PCR by adding template DNA, target-specific primers, and probes.
- Stability: Included stabilizer enables enzyme activity to be stable for up to 2 years at -20°C.
- Reproducibility: Mass production under ISO 9001 quality system allows minimized deviation between lots and reproducible results in replicated tests performed under same conditions and variation.

### Components

Components	Tube/Plate	50X ROX dye	DEPC-D.W.
K-6100	96 tubes		1.2 ml x 2 ea
K-6110	96 tubes		1.2 ml x 4 ea
K-6103	96-well plate		1.2 ml x 2 ea
K-6113	96-well plate		1.2 ml x 4 ea
K-6101	96 tubes	0.2 ml	1.2 ml x 2 ea
K-6111	96 tubes	0.2 ml	1.2 ml x 4 ea
K-6104	96-well plate	0.2 ml	1.2 ml x 2 ea
K-6114	96-well plate	0.2 ml	1.2 ml x 4 ea
K-6102	96 tubes		1.2 ml x 2 ea
K-6112	96 tubes		1.2 ml x 4 ea

\* **Note:** ROX dye is used for normalization of intensity by background subtraction. The use of ROX dye is recommended for Applied Biosystems 7500 Real-Time PCR System (Applied Biosystems), but not required for *Exicycler*™ 96 Real-Time PCR System (BIONEER) and CFX96 Real-Time PCR System (Bio-Rad).

### Composition

Composition	Concentration
Taq DNA Polymerase	1 U
dNTPs (dATP, dCTP, dGTP, dTTP)	1.2 mM
HotStart reaction buffer with 1.5 mM MgCl <sub>2</sub>	1X
Stabilizer	1X

### Specifications

Taq DNA Polymerase	
5' to 3' exonuclease activity	Yes
3' to 5' exonuclease activity	No
3'-A overhang	Yes

### Storage

Store at -20°C. If stored in the recommended temperature, this product will be stable until the expiration date printed out on the label.

### Online Resources



Korean



English

Visit our **product page** for additional information and protocols

### Ordering Information

Description		Cat. No.
Exicycler	8-tube strips 20 µl	96 rxn K-6100
	50 µl optical film included	96 rxn K-6110
	96-well plate 20 µl	96 rxn K-6103
	50 µl	96 rxn K-6113
ABI7500	8-tube strips 20 µl	96 rxn K-6101
	50 µl optical film included	96 rxn K-6111
	96-well plate 20 µl	96 rxn K-6104
	50 µl	96 rxn K-6114
CFX96	8-tube strips 20 µl optical film included	96 rxn K-6102
	50 µl	96 rxn K-6112

### Notice

BIONEER corporation reserves the right to make corrections, modifications, improvements and other changes to its products, services, specifications or product descriptions at any time without notice.

### Explanation of Symbols



Batch Code



Biological Risks



Catalog Number



Caution



Consult Instructions For Use



Contains Sufficient for <n> tests



Do not Re-use



Manufacturer



Research Use Only





Temperature Limitation



Use-by Date

**Experimental Procedures**

Steps	Procedure Details																																			
<p>Recommended protocol for <i>Exicycler™</i> 96 (BIONEER), Applied Biosystems 7500 Real-Time PCR System (Applied Biosystems), and CFX96 Real-Time PCR System (Bio-Rad).</p>																																				
<p style="text-align: center;">1</p> <div style="text-align: center;">  <p><b>Preparation of reaction mixture</b></p> </div>	<p>1. Add template DNA, target-specific primers, hydrolysis probe (not provided), 50X ROX dye (optional), and DEPC-D.W. into <i>AccuPower® DualStar™</i> qPCR PreMix tubes to make a total volume of 20 µl or 50 µl. Do not include dried pellet.</p> <ul style="list-style-type: none"> <li>Amount of template</li> </ul> <table border="1" data-bbox="528 629 1465 792"> <thead> <tr> <th rowspan="2">Template DNA</th> <th colspan="2">Amount of template</th> </tr> <tr> <th>20 µl reaction</th> <th>50 µl reaction</th> </tr> </thead> <tbody> <tr> <td>Total genomic DNA</td> <td>10 pg-1 µg</td> <td>10 pg-1 µg</td> </tr> <tr> <td>cDNA</td> <td>10 pg-1 µg</td> <td>10 pg-1 µg</td> </tr> </tbody> </table> <ul style="list-style-type: none"> <li>Preparation of reaction mixture</li> </ul> <table border="1" data-bbox="528 860 1465 1173"> <thead> <tr> <th>Components</th> <th>20 µl reaction</th> <th>50 µl reaction</th> </tr> </thead> <tbody> <tr> <td>Template DNA</td> <td>Variable</td> <td>Variable</td> </tr> <tr> <td>Forward primer (10 pmol/µl)</td> <td>1-2 µl</td> <td>2-5 µl</td> </tr> <tr> <td>Reverse primer (10 pmol/µl)</td> <td>1-2 µl</td> <td>2-5 µl</td> </tr> <tr> <td>Hydrolysis probe (10 pmol/µl)</td> <td>1-2 µl</td> <td>2-5 µl</td> </tr> <tr> <td>(Optional) 50X ROX dye</td> <td>0.4 µl</td> <td>1 µl</td> </tr> <tr> <td>DEPC-D.W.</td> <td>Variable</td> <td>Variable</td> </tr> <tr> <td>Total volume</td> <td>20 µl</td> <td>50 µl</td> </tr> </tbody> </table> <p>* <b>Note:</b> This protocol was validated with the TaqMan® probe as a hydrolysis probe.</p> <p>2. Seal real-time PCR tubes or plate with adhesive optical sealing film (Cat. No. 3111-4110, provided).</p> <p>3. Dissolve the vacuum-dried pellet by vortexing, and briefly spin down.</p>	Template DNA	Amount of template		20 µl reaction	50 µl reaction	Total genomic DNA	10 pg-1 µg	10 pg-1 µg	cDNA	10 pg-1 µg	10 pg-1 µg	Components	20 µl reaction	50 µl reaction	Template DNA	Variable	Variable	Forward primer (10 pmol/µl)	1-2 µl	2-5 µl	Reverse primer (10 pmol/µl)	1-2 µl	2-5 µl	Hydrolysis probe (10 pmol/µl)	1-2 µl	2-5 µl	(Optional) 50X ROX dye	0.4 µl	1 µl	DEPC-D.W.	Variable	Variable	Total volume	20 µl	50 µl
Template DNA	Amount of template																																			
	20 µl reaction	50 µl reaction																																		
Total genomic DNA	10 pg-1 µg	10 pg-1 µg																																		
cDNA	10 pg-1 µg	10 pg-1 µg																																		
Components	20 µl reaction	50 µl reaction																																		
Template DNA	Variable	Variable																																		
Forward primer (10 pmol/µl)	1-2 µl	2-5 µl																																		
Reverse primer (10 pmol/µl)	1-2 µl	2-5 µl																																		
Hydrolysis probe (10 pmol/µl)	1-2 µl	2-5 µl																																		
(Optional) 50X ROX dye	0.4 µl	1 µl																																		
DEPC-D.W.	Variable	Variable																																		
Total volume	20 µl	50 µl																																		
<p style="text-align: center;">2</p> <div style="text-align: center;">  <p><b>Real-time PCR</b></p> </div>	<p>4. Perform the reaction under the following conditions.</p> <table border="1" data-bbox="528 1464 1465 1621"> <thead> <tr> <th>Step</th> <th>Temperature</th> <th>Time</th> <th>Cycles</th> </tr> </thead> <tbody> <tr> <td>Pre-denaturation</td> <td>95°C</td> <td>3-5 min</td> <td>1 cycle</td> </tr> <tr> <td>Denaturation</td> <td>95°C</td> <td>5-30 sec</td> <td rowspan="2">40-45 cycles</td> </tr> <tr> <td>Annealing &amp; Extension</td> <td>55-60°C</td> <td>30-35 sec</td> </tr> </tbody> </table> <p>* <b>Note:</b> Users can adjust the protocol according to their instrument and template DNA sequences to get optimal results.</p> <p>5. After the reaction is completed, analyze the results.</p>	Step	Temperature	Time	Cycles	Pre-denaturation	95°C	3-5 min	1 cycle	Denaturation	95°C	5-30 sec	40-45 cycles	Annealing & Extension	55-60°C	30-35 sec																				
Step	Temperature	Time	Cycles																																	
Pre-denaturation	95°C	3-5 min	1 cycle																																	
Denaturation	95°C	5-30 sec	40-45 cycles																																	
Annealing & Extension	55-60°C	30-35 sec																																		