

[Cat. No.] **K-6707**

### Introduction

AccuPower® Dual-HotStart™ RT-qPCR Master Mix is a one-step RT-qPCR product based on hydrolysis probe method. It is applied Pyro-HotStart RT technology and HotStart PCR technology to substantially improve the problems of non-specific reactions and enhance the sensitivity. In reverse transcription, enzyme-mediated HotStart method that provides robust, sensitive, and reliable cDNA synthesis results by using BIONEER's RocketScript™ reverse transcriptase. In real-time PCR, the second HotStart reaction using antibody-based HotStart Taq DNA Polymerase provides superior priming accuracy and reduced non-specific reactions such as mis-priming and primer dimer during PCR at a low temperature. This product is ready-to-use mixture containing all components for real-time PCR, except for template RNA, target-specific primers, and probes. By just adding template RNA, target-specific primers, and probes, reproducible results with high sensitivity and specificity can be obtained. Moreover, this product is a 2X Master Mix type, so it is easy to compatible with other company's equipment.

### Applications

- Gene expression profiling
- Target RNA quantification
- Microbial detection
- Viral/bacterial pathogen load determination

### Features & Benefits

- **Specificity:** The world's first Dual-HotStart™ RT-qPCR reaction provides accurate results of target gene amplification by using Pyro-HotStart RT reaction and HotStart PCR.
- **Sensitivity:** It can be detected from even a trace amount of template RNA with a wide dynamic range of 10 logs up to 10-10<sup>10</sup> copies.
- **Multiplex:** Compatible with many kinds of fluorescent dye (probes) to detect various kinds of target genes.
- **Comprehensive template RNA detection:** Included RocketScript™ RTase can perform reverse transcription at high temperature and even with secondary RNA structures.
- **Compatibility:** Wide range of samples such as blood and soil containing various PCR inhibitors can be applied due to its excellent reactivity.
- **Convenience:** Reactants are individually packaged in a tube, it allows any user simply perform cDNA synthesis and real-time PCR by adding template RNA and target-specific primers and probe.
- **Reproducibility:** Mass production under ISO 9001 quality system allows minimized deviation between lots and reproducible results in replicated tests performed under same conditions and variation.

### Components

Components	K-6707
2X Master Mix	1.25 ml x 2 ea
10X HotStart buffer	0.5 ml
50X ROX dye	0.1 ml
DEPC-D.W.	1.2 ml

### Composition

2X Master Mix	Concentration
HotStart Taq DNA Polymerase	1 U
RocketScript Reverse Transcriptase	200 U
10X Reaction buffer with 2 mM MgCl <sub>2</sub>	1X
10X HotStart buffer	1X
dNTPs (dATP, dCTP, dGTP, dTTP)	1.2 mM

### Specifications

HotStart Taq DNA Polymerase	
5' to 3' exonuclease activity	Yes
3' to 5' exonuclease activity	No
3'-A overhang	Yes

### Storage

Store at -20°C. If stored in the recommended temperature, this product will be stable until the expiration date printed out on the label.

### Online Resources



Korean



English

Visit our **product page** for additional information and protocols


### Ordering Information

Description	Cat. No.
2.5 ml of 2X Master Mix solution	1.25 ml x 2 ea K-6707




### Notice

BIONEER corporation reserves the right to make corrections, modifications, improvements and other changes to its products, services, specifications or product descriptions at any time without notice.

### Explanation of Symbols



**Experimental Procedures**

Steps		Procedure Details																				
1	 <b>Thaw reagents</b>	<ol style="list-style-type: none"> <li>1. Thaw <i>AccuPower® Dual-HotStart™</i> RT-qPCR Master Mix on ice and mix thoroughly before use. Then, briefly spin down components.</li> <li>2. Dispense appropriate volumes of <i>AccuPower® Dual-HotStart™</i> RT-qPCR Master Mix into PCR tubes (not provided). Use 25 µl of 2X Master Mix for 50 µl reaction.</li> </ol>																				
2	 <b>Preparation of reaction mixture</b>	<ol style="list-style-type: none"> <li>3. Add template RNA, target-specific primers, hydrolysis probe (not provided), 10X HotStart buffer, 50X ROX dye (optional), and DEPC-D.W. into PCR tubes containing <i>AccuPower® Dual-HotStart™</i> RT-qPCR Master Mix. <ul style="list-style-type: none"> <li>• Preparation of reaction mixture</li> </ul> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="text-align: left;">Components</th> <th style="text-align: left;">50 µl reaction</th> </tr> </thead> <tbody> <tr> <td>2X Master Mix</td> <td>25 µl</td> </tr> <tr> <td>Template RNA (10 pg-100 ng)</td> <td>Variable</td> </tr> <tr> <td>Forward primer (10 pmol/µl)</td> <td>0.5-5 µl</td> </tr> <tr> <td>Reverse primer (10 pmol/µl)</td> <td>0.5-5 µl</td> </tr> <tr> <td>Hydrolysis probe (10 pmol/µl)</td> <td>0.5-5 µl</td> </tr> <tr> <td>10X HotStart buffer</td> <td>5 µl</td> </tr> <tr> <td>(Optional) 50X ROX dye</td> <td>1 µl</td> </tr> <tr> <td>DEPC-D.W.</td> <td>Variable</td> </tr> <tr> <td>Total volume</td> <td>50 µl</td> </tr> </tbody> </table> </li> </ol> <p><b>* Note:</b> This protocol was validated with the TaqMan® probe as a hydrolysis probe</p> <ol style="list-style-type: none"> <li>4. Seal real-time PCR tubes with adhesive optical sealing film (Cat. No. 3111-4110, not provided).</li> <li>5. Mix the reaction mixture by vortexing, and briefly spin down.</li> </ol>	Components	50 µl reaction	2X Master Mix	25 µl	Template RNA (10 pg-100 ng)	Variable	Forward primer (10 pmol/µl)	0.5-5 µl	Reverse primer (10 pmol/µl)	0.5-5 µl	Hydrolysis probe (10 pmol/µl)	0.5-5 µl	10X HotStart buffer	5 µl	(Optional) 50X ROX dye	1 µl	DEPC-D.W.	Variable	Total volume	50 µl
Components	50 µl reaction																					
2X Master Mix	25 µl																					
Template RNA (10 pg-100 ng)	Variable																					
Forward primer (10 pmol/µl)	0.5-5 µl																					
Reverse primer (10 pmol/µl)	0.5-5 µl																					
Hydrolysis probe (10 pmol/µl)	0.5-5 µl																					
10X HotStart buffer	5 µl																					
(Optional) 50X ROX dye	1 µl																					
DEPC-D.W.	Variable																					
Total volume	50 µl																					
3	 <b>RT-qPCR</b>	<ol style="list-style-type: none"> <li>6. Perform the reaction under the following conditions. <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="text-align: left;">Step</th> <th style="text-align: left;">Temperature</th> <th style="text-align: left;">Time</th> <th style="text-align: left;">Cycles</th> </tr> </thead> <tbody> <tr> <td>cDNA synthesis</td> <td>50-70°C</td> <td>15 min</td> <td>1 cycle</td> </tr> <tr> <td>Pre-denaturation</td> <td>95°C</td> <td>3-5 min</td> <td>1 cycle</td> </tr> <tr> <td>Denaturation</td> <td>95°C</td> <td>5-30 sec</td> <td rowspan="2">40-45 cycles</td> </tr> <tr> <td>Annealing &amp; Extension</td> <td>55-60°C</td> <td>30-35 sec</td> </tr> </tbody> </table> </li> </ol> <p><b>* Note:</b> Users can adjust the protocol according to their instrument and template sequences to get optimal results.</p> <ol style="list-style-type: none"> <li>7. After the reaction is completed, analyze the results.</li> </ol>	Step	Temperature	Time	Cycles	cDNA synthesis	50-70°C	15 min	1 cycle	Pre-denaturation	95°C	3-5 min	1 cycle	Denaturation	95°C	5-30 sec	40-45 cycles	Annealing & Extension	55-60°C	30-35 sec	
Step	Temperature	Time	Cycles																			
cDNA synthesis	50-70°C	15 min	1 cycle																			
Pre-denaturation	95°C	3-5 min	1 cycle																			
Denaturation	95°C	5-30 sec	40-45 cycles																			
Annealing & Extension	55-60°C	30-35 sec																				