

[Cat. No.] Please refer to the **Ordering Information**

### Introduction

AccuPower® HotStart ProFi Taq PCR PreMix uses HotStart ProFi Taq DNA Polymerase maximizes PCR efficiency by removing PCR inhibitor, and it allows amplification of long template DNA with high fidelity. Moreover, by applying HotStart technology, this prevents the formation of mis-primed products and primer-dimers during the reaction set up process resulting in improved PCR specificity. This product contains vacuum-dried components for PCR including, HotStart ProFi Taq DNA Polymerase, dNTPs, reaction buffer, stabilizer, and tracking dye. It simplifies preparation of reaction mixture by adding template DNA and primers without any extra process. After the reaction, samples can be applied directly on agarose gel for analysis.

### Applications

- High specificity PCR
- Long-range amplification from genomic DNA
- High yield and high sensitivity PCR

### Features & Benefits

- **Specificity:** Minimized amplification of non-specific PCR products, maximized reaction efficiency, and effective amplification even with a small amount of template DNA through the application of HotStart technology.
- **Long-range PCR:** Effective amplification of large genomic DNA fragments up to 20 kb of human DNA.
- **User-friendly:** Reactants are individually packaged in each of the PCR tubes, it allows any user simply perform PCR by adding template DNA and primers. The tube also has tracking dye and sedimentation agents for electrophoresis, so sample loading buffer is not required.
- **Stability:** Included stabilizer provides prolonged storage and increased stability compared to solution-type products.
- **Reproducibility:** Reproducible experimental results with minimized lot-to-lot variation under ISO 9001 quality system.

### Composition

Composition	Concentration
HotStart ProFi Taq DNA Polymerase	1 U
dNTPs (dATP, dCTP, dGTP, dTTP)	Each 250 µM
Reaction buffer with 1.5 mM MgCl <sub>2</sub>	1X
Stabilizer and tracking dye	0

\* **Note:** For research use only. Not for use in diagnostic or therapeutic procedures.

### Specifications

HotStart ProFi Taq DNA Polymerase	
5' to 3' exonuclease activity	Yes
3' to 5' exonuclease activity	Yes
3'-A overhang	Yes
Fragment size	Up to 20 kb

### Storage

Store at -20°C. This product can be stable up to 2 years in a freezer.

### Online Resources



Korean



English

Visit our **product page** for additional information and protocols.


### Ordering Information

Description	Cat. No.
96 tubes 20 µl/rxn	K-2640
0.2 ml thin-wall 8-tube strips with attached cap	480 tubes 20 µl/rxn K-2641
96 tubes 50 µl/rxn	K-2642
480 tubes 50 µl/rxn	K-2643




### Notice

BIONEER corporation reserves the right to make corrections, modifications, improvements and other changes to its products, services, specifications or product descriptions at any time without notice.

### Explanation of Symbols



**Experimental Procedures**

Steps		Procedure Details																																							
1	 <b>Preparation of reaction mixture</b>	<p>1. Add template DNA, primers and nuclease-free water into <i>AccuPower® HotStart ProFi Taq PCR PreMix</i> tubes to a total volume of 20 µl or 50 µl. Do not calculate the dried pellet.</p> <ul style="list-style-type: none"> <li>Preparation of reaction mixture</li> </ul> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="text-align: center;">Components</th> <th style="text-align: center;">20 µl reaction</th> <th colspan="2" style="text-align: center;">50 µl reaction</th> </tr> </thead> <tbody> <tr> <td>Template DNA</td> <td style="text-align: center;">1-500 ng</td> <td colspan="2" style="text-align: center;">1-500 ng</td> </tr> <tr> <td>Forward primer (10 pmol/µl)</td> <td style="text-align: center;">0.5-2 µl</td> <td colspan="2" style="text-align: center;">0.5-2 µl</td> </tr> <tr> <td>Reverse primer (10 pmol/µl)</td> <td style="text-align: center;">0.5-2 µl</td> <td colspan="2" style="text-align: center;">0.5-2 µl</td> </tr> <tr> <td>Nuclease-free water</td> <td style="text-align: center;">Variable</td> <td colspan="2" style="text-align: center;">Variable</td> </tr> <tr> <td>Total volume</td> <td style="text-align: center;">20 µl</td> <td colspan="2" style="text-align: center;">50 µl</td> </tr> </tbody> </table> <p>2. Dissolve the vacuum-dried pellet completely by tapping or pipetting, and briefly spin down. If you are using the BIONEER's <i>ExiSpin™</i> Vortex/Centrifuge, follow the recommended settings below.</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="text-align: center;">Step</th> <th colspan="3" style="text-align: center;">Setting</th> </tr> </thead> <tbody> <tr> <td style="text-align: center;">Vortex</td> <td style="text-align: center;">High</td> <td style="text-align: center;">15 sec</td> <td style="text-align: center;">4 cycles</td> </tr> <tr> <td style="text-align: center;">Spin</td> <td style="text-align: center;">1,500 rpm</td> <td style="text-align: center;">5 sec</td> <td></td> </tr> </tbody> </table>				Components	20 µl reaction	50 µl reaction		Template DNA	1-500 ng	1-500 ng		Forward primer (10 pmol/µl)	0.5-2 µl	0.5-2 µl		Reverse primer (10 pmol/µl)	0.5-2 µl	0.5-2 µl		Nuclease-free water	Variable	Variable		Total volume	20 µl	50 µl		Step	Setting			Vortex	High	15 sec	4 cycles	Spin	1,500 rpm	5 sec	
Components	20 µl reaction	50 µl reaction																																							
Template DNA	1-500 ng	1-500 ng																																							
Forward primer (10 pmol/µl)	0.5-2 µl	0.5-2 µl																																							
Reverse primer (10 pmol/µl)	0.5-2 µl	0.5-2 µl																																							
Nuclease-free water	Variable	Variable																																							
Total volume	20 µl	50 µl																																							
Step	Setting																																								
Vortex	High	15 sec	4 cycles																																						
Spin	1,500 rpm	5 sec																																							
2	 <b>Perform the reaction in a thermal cycler</b>	<p>3. Perform the reaction under the following conditions.</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="text-align: center;">Step</th> <th style="text-align: center;">Temperature</th> <th style="text-align: center;">Time</th> <th style="text-align: center;">Cycles</th> </tr> </thead> <tbody> <tr> <td>Pre-denaturation</td> <td style="text-align: center;">95°C</td> <td style="text-align: center;">5 min</td> <td style="text-align: center;">1 cycle</td> </tr> <tr> <td>Denaturation</td> <td style="text-align: center;">95°C</td> <td style="text-align: center;">15-30 sec</td> <td></td> </tr> <tr> <td>Annealing</td> <td style="text-align: center;">45-65°C</td> <td style="text-align: center;">15-30 sec</td> <td style="text-align: center;">25-35 cycles</td> </tr> <tr> <td>Extension</td> <td style="text-align: center;">72°C</td> <td style="text-align: center;">1 min/kb</td> <td></td> </tr> <tr> <td>Final extension</td> <td style="text-align: center;">72°C</td> <td style="text-align: center;">3-5 min</td> <td style="text-align: center;">1 cycle</td> </tr> </tbody> </table> <p>* <b>Note:</b> The annealing temperature depends on the melting point of the primers.</p>				Step	Temperature	Time	Cycles	Pre-denaturation	95°C	5 min	1 cycle	Denaturation	95°C	15-30 sec		Annealing	45-65°C	15-30 sec	25-35 cycles	Extension	72°C	1 min/kb		Final extension	72°C	3-5 min	1 cycle												
Step	Temperature	Time	Cycles																																						
Pre-denaturation	95°C	5 min	1 cycle																																						
Denaturation	95°C	15-30 sec																																							
Annealing	45-65°C	15-30 sec	25-35 cycles																																						
Extension	72°C	1 min/kb																																							
Final extension	72°C	3-5 min	1 cycle																																						
3	 <b>Analyze with gel electrophoresis</b>	<p>4. After the reaction, maintain the reaction mixture at 4°C. The samples can be stored at -20°C until use.</p> <p>5. Load 5 µl of PCR product samples on agarose gel without adding a loading-dye mixture, and perform gel electrophoresis for analysis.</p>																																							