

[Cat. No.] **K-2940**

Introduction

AccuPower® rPLU/rVIV PCR Kit is a ready-to-use premix for nested PCR that can be used to detect malaria infection.

Plasmodium spp. is the causative organism of malaria, which is an epidemic that has a high risk and a mortality rate. Each year, it infects more than 200 million people worldwide. When bitten by an infected female mosquito, people may show various symptoms, such as headache, splenomegaly, jaundice, and anemia. Patients may also show signs of chill and high fever above 39-41°C. In the case of South Korea, *Plasmodium vivax* is the major cause of malaria. Polymerase chain reaction (PCR) is used for an accurate diagnosis of the infectious agent.

This product contains vacuum-dried components specific to *Plasmodium vivax* including DNA polymerase, primers, dNTPs, and reaction buffer required for PCR. This ready-to-use kit simplifies preparation of PCR mixture as the user only has to add template DNA and nuclease-free water. After the reaction, since tracking dye is included, the samples can be applied directly on agarose gel for analysis without adding extra solution.

Features & Benefits

- **Nested PCR:** Detect *Plasmodium vivax* with nested PCR consisting of rPLU PCR Kit and rVIV PCR Kit.
- **Convenience & Reproducibility:** All reactants necessary for PCR including primers are lyophilized in each PCR tube, providing reproducible results in a convenient way.
- **Sensitivity:** By applying the patented PyroHotStart (Enzyme-mediated HotStart) technology that minimizes non-specific reactions and maximizes reaction efficiency, only the target gene can be effectively amplified even with a trace amount of template DNA.
- **Stability:** Included stabilizer in the PCR reaction mixture provides increased stability compared to solution-type products.

Composition

Composition	20 µl reaction	
	1 st PCR	2 nd PCR
Top DNA Polymerase	1 U	1 U
dNTPs (dATP, dCTP, dGTP, dTTP)	Each 250 µM	Each 250 µM
Reaction buffer with 1.5 mM MgCl ₂	1X	1X
Stabilizer and tracking dye	O	O
rPLU 5 Forward primer	0.45 µM	X
rPLU 6 Reverse primer	0.45 µM	X
rVIV Forward primer	X	0.55 µM
rVIV Reverse primer	X	0.55 µM

* **Note:** For research use only. Not for use in diagnostic or therapeutic procedures.

Specifications

Top DNA Polymerase	
5'→3' exonuclease activity	No
3'→5' exonuclease activity	No
3'-A overhang	Yes

1 st fragment size	1,200 bp
2 nd fragment size	120 bp

Storage

Store at -20°C. If stored in the recommended temperature, this product will be stable until the expiration date printed out on the label.

Online Resources



English

Visit our **product page** for additional information and protocols

Ordering Information

Description	Cat. No.
AccuPower® rPLU PCR Kit, 0.2 ml thin-wall 8-tube strips with attached cap / 96 tubes	K-2940
AccuPower® rVIV PCR Kit, 0.2 ml thin-wall 8-tube strips with attached cap / 96 tubes	






Notice

BIONEER corporation reserves the right to make corrections, modifications, improvements and other changes to its products, services, specifications or product descriptions at any time without notice.

Explanation of Symbols

Batch Code	Biological Risks	Catalog Number	Caution
Consult Instructions For Use	Contains Sufficient for <n> tests	Do not Re-use	Manufacturer
Research Use Only	Temperature Limitation	Use-by Date	

Experimental Procedures

Steps		Procedure Details																						
1	 Preparation of 1st reaction mixture	<ol style="list-style-type: none"> After preparing the template DNA and nuclease-free water, add the template DNA to the <i>AccuPower®</i> rPLU PCR Kit. Add nuclease-free water into PCR tubes to make a total volume of 20 µl. (Do not include the volume of the dried premix in the PCR tubes.) Completely dissolve the vacuum-dried pellet by vortexing, and briefly spin down. 																						
2	 Incubate reactions in a thermal cycler	<ol style="list-style-type: none"> Place PCR tubes on the thermal cycler and perform the reaction under the following conditions. <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Step</th> <th>Temperature</th> <th>Time</th> <th>Cycles</th> </tr> </thead> <tbody> <tr> <td>Pre-denaturation</td> <td>95°C</td> <td>5 min</td> <td>1 cycle</td> </tr> <tr> <td>Denaturation</td> <td>92°C</td> <td>30 sec</td> <td rowspan="2">30 cycles</td> </tr> <tr> <td>Annealing</td> <td>60°C</td> <td>60 sec</td> </tr> <tr> <td>Extension</td> <td>72°C</td> <td>60 sec</td> <td rowspan="2">1 cycle</td> </tr> <tr> <td>Final extension</td> <td>72°C</td> <td>5 min</td> </tr> </tbody> </table> <p>* Note: Users can adjust the protocol according to their instrument and template sequences to get optimal results.</p>	Step	Temperature	Time	Cycles	Pre-denaturation	95°C	5 min	1 cycle	Denaturation	92°C	30 sec	30 cycles	Annealing	60°C	60 sec	Extension	72°C	60 sec	1 cycle	Final extension	72°C	5 min
Step	Temperature	Time	Cycles																					
Pre-denaturation	95°C	5 min	1 cycle																					
Denaturation	92°C	30 sec	30 cycles																					
Annealing	60°C	60 sec																						
Extension	72°C	60 sec	1 cycle																					
Final extension	72°C	5 min																						
3	 Preparation of 2nd reaction mixture	<ol style="list-style-type: none"> After the 1st PCR reaction is completed, add 1-5 µl of the reaction solution into the <i>AccuPower®</i> rVIV PCR Kit. Add nuclease-free water into PCR tubes to make a total volume of 20 µl. (Do not include the volume of the dried premix in the PCR tubes.) Completely dissolve the vacuum-dried pellet by vortexing, and briefly spin down. 																						
4	 Incubate reactions in a thermal cycler	<ol style="list-style-type: none"> Place PCR tubes on the thermal cycler and perform the reaction under the following conditions. <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Step</th> <th>Temperature</th> <th>Time</th> <th>Cycles</th> </tr> </thead> <tbody> <tr> <td>Pre-denaturation</td> <td>95°C</td> <td>5 min</td> <td>1 cycle</td> </tr> <tr> <td>Denaturation</td> <td>92°C</td> <td>20 sec</td> <td rowspan="2">30 cycles</td> </tr> <tr> <td>Annealing</td> <td>58°C</td> <td>20 sec</td> </tr> <tr> <td>Extension</td> <td>72°C</td> <td>20 sec</td> <td rowspan="2">1 cycle</td> </tr> <tr> <td>Final extension</td> <td>72°C</td> <td>5 min</td> </tr> </tbody> </table> <p>* Note: Users can adjust the protocol according to their instrument and template sequences to get optimal results.</p>	Step	Temperature	Time	Cycles	Pre-denaturation	95°C	5 min	1 cycle	Denaturation	92°C	20 sec	30 cycles	Annealing	58°C	20 sec	Extension	72°C	20 sec	1 cycle	Final extension	72°C	5 min
Step	Temperature	Time	Cycles																					
Pre-denaturation	95°C	5 min	1 cycle																					
Denaturation	92°C	20 sec	30 cycles																					
Annealing	58°C	20 sec																						
Extension	72°C	20 sec	1 cycle																					
Final extension	72°C	5 min																						
5	 Analyze with gel electrophoresis	<ol style="list-style-type: none"> After the reaction, maintain the reaction mixture at 4-8°C. Load samples on agarose gel without adding a loading-dye mixture, and perform gel electrophoresis for analysis. 																						