

[Cat. No.] **K-2992**

Introduction

AccuPower® PRV Master Mix is a product that can detect Piscine orthoreovirus, which causes the disease heart and skeletal muscle inflammation (HSMI) in farmed Atlantic salmon, by real-time PCR.

Piscine orthoreovirus (PRV) is a serious virus in salmon aquaculture belonging to the family Reoviridae. PRV has been found present at higher concentration in fish with various diseases. This disease includes HSMI, jaundice syndrome, proliferative darkening syndrome and erythrocytic body inclusion syndrome. PRV is thought to mainly affect aquacultured and maricultured fish stocks.

This product contains all real-time PCR components specific to Piscine orthoreovirus, including RTase, DNA polymerase, primers, dNTPs, Reaction buffer. The users can easily prepare reaction mixture simply by adding template RNA, internal positive control (IPC), oligo mix, and DEPC-D.W.

Features & Benefits

- Convenience: All necessary reactants for real-time PCR are included in a tube (i.e., Master Mix Type), allowing the users to perform reaction simply by adding template RNA, oligo mix, and DEPC-D.W.
- Sensitivity: By using BIONEER's PyroHotStart RT reaction and HotStart Taq DNA Polymerase that minimizes non-specific reactions and maximizes reaction efficiency, only the target gene can be effectively amplified even with a trace amount of template RNA.

Components

Components	Amount
2X Master Mix	625 µl x 2ea
Oligo Mix	500 µl
Positive control (1*10 ⁸ copies/µl)	50 µl
DEPC-D.W.	1.3 ml

* **Note:** For research use only. Not for use in diagnostic or therapeutic procedures.

Composition

Composition	25 µl reaction
RocketScript™ Reverse transcriptase	0.5 U
2X Master Mix	Taq DNA Polymerase 3.5 U
	dNTPs (dATP, dCTP, dGTP, dTTP) Each 300 µM
	Reaction buffer with 2 mM MgCl ₂ 1X
	PRV Forward primer 0.4 µM
Oligo Mix	PRV Reverse primer 0.4 µM
	PRV Probe (FAM) 0.4 µM

Specifications

Taq DNA Polymerase	
5'→3' exonuclease activity	Yes
3'→5' exonuclease activity	No
3'-A overhang	Yes

Storage

Store at -20°C. If stored in the recommended temperature, this product will be stable until the expiration date printed out on the label.

Online Resources



Visit our **product page** for additional information and protocols.


Ordering Information

Description	Cat. No.
AccuPower® PRV Master Mix 1.25 ml of 2X Master Mix solution, 100 tests	K-2992




Notice

BIONEER corporation reserves the right to make corrections, modifications, improvements and other changes to its products, services, specifications or product descriptions at any time without notice.

Explanation of Symbols



Experimental Procedures

Steps		Procedure Details																							
1	 Preparation of reaction mixture	<p>1. Before use, thaw all components of <i>AccuPower®</i> PRV Master Mix on ice and mix them thoroughly. Then, briefly spin down components.</p>																							
2	 Composition of reaction mixture	<p>2. Add all components into PCR tubes or a plate referring to the following list of components. (based on 1 test)</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="text-align: left;">Components</th> <th style="text-align: right;">Volume (ul)</th> </tr> </thead> <tbody> <tr> <td>2X Mater Mix</td> <td style="text-align: right;">12.5</td> </tr> <tr> <td>Oligo Mix</td> <td style="text-align: right;">5</td> </tr> <tr> <td>Template RNA (Positive Control)</td> <td style="text-align: right;">1~5</td> </tr> <tr> <td>DEPC-DW</td> <td style="text-align: right;">Up to 25</td> </tr> <tr> <td>Total reaction volume</td> <td style="text-align: right;">25</td> </tr> </tbody> </table>	Components	Volume (ul)	2X Mater Mix	12.5	Oligo Mix	5	Template RNA (Positive Control)	1~5	DEPC-DW	Up to 25	Total reaction volume	25											
Components	Volume (ul)																								
2X Mater Mix	12.5																								
Oligo Mix	5																								
Template RNA (Positive Control)	1~5																								
DEPC-DW	Up to 25																								
Total reaction volume	25																								
3	 Real-time PCR	<p>3. Place PCR tubes or plate on the real-time quantitative thermal cycler.</p> <p>4. Perform the reaction under the following conditions.</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="text-align: left;">Step</th> <th style="text-align: left;">Temperature</th> <th style="text-align: left;">Time</th> <th style="text-align: left;">Cycles</th> </tr> </thead> <tbody> <tr> <td>Reverse Transcription</td> <td>50 °C</td> <td>15 min</td> <td>1 cycle</td> </tr> <tr> <td>Pre-denaturation</td> <td>95 °C</td> <td>5 min</td> <td>1 cycle</td> </tr> <tr> <td>Denaturation</td> <td>95 °C</td> <td>10 sec</td> <td rowspan="2">45 cycles</td> </tr> <tr> <td>Annealing& Extension</td> <td>60 °C</td> <td>20 sec</td> </tr> <tr> <td>Scan</td> <td></td> <td></td> <td></td> </tr> </tbody> </table> <p>* Note: Users can adjust the protocol according to their instrument and template sequences to get optimal results.</p> <p>5. After the reaction is completed, analyze the results.</p>	Step	Temperature	Time	Cycles	Reverse Transcription	50 °C	15 min	1 cycle	Pre-denaturation	95 °C	5 min	1 cycle	Denaturation	95 °C	10 sec	45 cycles	Annealing& Extension	60 °C	20 sec	Scan			
Step	Temperature	Time	Cycles																						
Reverse Transcription	50 °C	15 min	1 cycle																						
Pre-denaturation	95 °C	5 min	1 cycle																						
Denaturation	95 °C	10 sec	45 cycles																						
Annealing& Extension	60 °C	20 sec																							
Scan																									