

[Cat. No.] **E-3121, E-3122**

Introduction

Moloney Murine Leukemia Virus (*M-MLV*) Reverse Transcriptase (RTase) is an RNA-dependent DNA polymerase. It can synthesize complementary DNA by using RNA molecules as a template. *M-MLV* RTase is isolated from a *E. coli* strain containing a recombinant clone and it can synthesize cDNA up to 9 kb. Furthermore, by using the RNase H+ of *M-MLV* RTase, template RNA is removed after the cDNA synthesis to minimize the effect of residual template RNA in PCR. This product can be applied for first-strand cDNA synthesis for RT-PCR and qRT-PCR.

Applications

- First-strand cDNA synthesis from RNA
- RT-PCR and qRT-PCR

Components

Components	E-3121	E-3122
<i>M-MLV</i> Reverse Transcriptase	10,000 U (50 µl)	50,000 U (50 µl x 5)
5X Reaction buffer	0.5 ml	0.5 ml x 5
10 mM dNTPs	0.2 ml	0.2 ml x 5
100 mM DTT	0.3 ml	0.3 ml x 5

* **Note:** For research use only. Not for use in diagnostic or therapeutic procedures.

Specifications

<i>M-MLV</i> Reverse Transcriptase	
5' to 3' exonuclease activity	No
3' to 5' exonuclease activity	No
3'-A overhang	No
Fragment size	Up to 9 kb

Buffer Composition

5X Reaction buffer	150 mM Tris-HCl, 250 mM KCl, 10 mM MgCl ₂ , pH 8.1
--------------------	---

Storage Buffer

M-MLV Reverse Transcriptase is supplied in 50% (v/v) glycerol containing 20 mM Tris-HCl, 150 mM NaCl, 0.1 mM EDTA, 1 mM DTT, 50 mM (NH₄)₂SO₄, and stabilizer, pH 7.6.

Unit Definition

One unit is defined as the amount of enzyme that incorporates 1 nmole of dTTP into acid-insoluble products in 10 min at 37°C using poly(A)-oligo(dT) as template primer.

Quality Control

- Nuclease Contamination Assay: Nuclease activity is not detected after incubation of 1 µg of DNA and RNA with 200 U of *M-MLV* Reverse Transcriptase at 37-42°C for 3 hrs.

Storage

Store at -20°C. If stored in the recommended temperature, this product will be stable until the expiration date printed out on the label.

Online Resources



Korean



English

Visit our [product page](#) for additional information and protocols.

Ordering Information

Description	Cat. No
<i>M-MLV</i> Reverse Transcriptase	E-3121
10,000 U (50 rxn)	
50,000 U (250 rxn)	E-3122

Notice

BIONEER corporation reserves the right to make corrections, modifications, improvements and other changes to its products, services, specifications or product descriptions at any time without notice.

Explanation of Symbols



Batch Code



Catalog Number



Caution



Consult Instructions For Use



Contains Sufficient for <n> tests



Research Use Only






Temperature Limitation



Use-by Date

Experimental Procedures

Steps		Procedure Details																		
1	 Primer annealing	<p>1. Mix template RNA and primers in a sterile tube (not provided) indicated as below.</p> <ul style="list-style-type: none"> Amount of template RNA and primers <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th colspan="2">Components</th> <th>20 µl reaction</th> </tr> </thead> <tbody> <tr> <td rowspan="2" style="text-align: center;">Template RNA</td> <td style="text-align: center;">Total RNA</td> <td style="text-align: center;">1 µg</td> </tr> <tr> <td style="text-align: center;">RNA</td> <td style="text-align: center;">5-100 ng</td> </tr> <tr> <td rowspan="2" style="text-align: center;">Primers</td> <td style="text-align: center;">Oligo dT or Random primer</td> <td style="text-align: center;">10-100 pmol</td> </tr> <tr> <td style="text-align: center;">Sequence specific</td> <td style="text-align: center;">10-30 pmol</td> </tr> </tbody> </table> <p>2. Incubate the mixture at 65°C for 10 min and place it on ice directly.</p>	Components		20 µl reaction	Template RNA	Total RNA	1 µg	RNA	5-100 ng	Primers	Oligo dT or Random primer	10-100 pmol	Sequence specific	10-30 pmol					
Components		20 µl reaction																		
Template RNA	Total RNA	1 µg																		
	RNA	5-100 ng																		
Primers	Oligo dT or Random primer	10-100 pmol																		
	Sequence specific	10-30 pmol																		
2	 Preparation of reaction mixture	<p>3. Transfer the incubated mixture to other sterile tubes (not provided) and mix with 5X Reaction buffer, 10 mM dNTPs, 100 mM DTT, RNase inhibitor (not provided), M-MLV Reverse Transcriptase, and nuclease-free water to a total volume of 20 µl.</p> <ul style="list-style-type: none"> Preparation of reaction mixture <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Components</th> <th>20 µl reaction</th> </tr> </thead> <tbody> <tr> <td>Step 1 mixture</td> <td style="text-align: center;">Variable</td> </tr> <tr> <td>5X Reaction buffer</td> <td style="text-align: center;">4 µl</td> </tr> <tr> <td>10 mM dNTPs</td> <td style="text-align: center;">2 µl or Variable</td> </tr> <tr> <td>100 mM DTT</td> <td style="text-align: center;">2 µl</td> </tr> <tr> <td>RNase inhibitor (Not provided)</td> <td style="text-align: center;">Variable</td> </tr> <tr> <td>M-MLV Reverse Transcriptase (200 U/ µl)</td> <td style="text-align: center;">200 U</td> </tr> <tr> <td>Nuclease-free water</td> <td style="text-align: center;">Variable</td> </tr> <tr> <td>Total volume</td> <td style="text-align: center;">20 µl</td> </tr> </tbody> </table>	Components	20 µl reaction	Step 1 mixture	Variable	5X Reaction buffer	4 µl	10 mM dNTPs	2 µl or Variable	100 mM DTT	2 µl	RNase inhibitor (Not provided)	Variable	M-MLV Reverse Transcriptase (200 U/ µl)	200 U	Nuclease-free water	Variable	Total volume	20 µl
Components	20 µl reaction																			
Step 1 mixture	Variable																			
5X Reaction buffer	4 µl																			
10 mM dNTPs	2 µl or Variable																			
100 mM DTT	2 µl																			
RNase inhibitor (Not provided)	Variable																			
M-MLV Reverse Transcriptase (200 U/ µl)	200 U																			
Nuclease-free water	Variable																			
Total volume	20 µl																			
3	 cDNA synthesis	<p>4. Perform the reaction under the following conditions.</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Step</th> <th>Temperature</th> <th>Time</th> </tr> </thead> <tbody> <tr> <td>cDNA synthesis</td> <td style="text-align: center;">37-42°C</td> <td style="text-align: center;">1 hr</td> </tr> <tr> <td>Heat inactivation</td> <td style="text-align: center;">95°C</td> <td style="text-align: center;">5 min</td> </tr> </tbody> </table> <p>5. After the reaction, maintain the reaction mixture at 4°C. The samples can be stored at -20°C until use.</p>	Step	Temperature	Time	cDNA synthesis	37-42°C	1 hr	Heat inactivation	95°C	5 min									
Step	Temperature	Time																		
cDNA synthesis	37-42°C	1 hr																		
Heat inactivation	95°C	5 min																		