

[Cat. No.] **E-3141, E-3142**

## Introduction

RocketScript™ Reverse Transcriptase developed by BIONEER is originated from *M-MLV* Reverse Transcriptase and genetically engineered to provide increased thermal stability and outstanding processivity. Its enhanced thermal stability enables to effectively synthesize cDNA with complex secondary RNA structures. It can be used at various temperatures (42-70°C) allowing increased specificity, higher yields of cDNA, and more full-length product than other reverse transcriptase.

## Applications

- Gene synthesis
- First-strand cDNA synthesis from RNA molecules
- RT-PCR
- Random priming reactions
- Library construction
- Probe labeling
- mRNA 5'end mapping by primer extension analysis
- Real-time PCR

## Components

Components	E-3141	E-3142
RocketScript™ Reverse Transcriptase	10,000 U (50 µl)	50,000 U (50 x 5 µl)
5X Reaction buffer	0.5 ml	0.5 ml x 5
10 mM dNTPs	0.2 ml	0.2 ml x 5
100 mM DTT	0.3 ml	0.3 ml x 5
RNase inhibitor	50 µl	50 µl x 5

\* Note: For research use only. Not for use in diagnostic or therapeutic procedures.

## Specifications

RocketScript™ Reverse Transcriptase	
5' to 3' exonuclease activity	No
3' to 5' exonuclease activity	No
3'-A overhang	No
Fragment size	Up to 10 kb

## Buffer Composition

5X Reaction buffer	250 mM Tris-HCl, 375 mM KCl, 15 mM MgCl <sub>2</sub> , and stabilizer, pH 8.5
--------------------	---

## Storage Buffer

RocketScript™ Reverse Transcriptase is supplied in 50% (v/v) glycerol containing 20 mM Tris-HCl, 150 mM NaCl, 0.1 mM EDTA, 1 mM DTT, and stabilizer, pH 7.6.

## Unit Definition

One unit is defined as the amount of enzyme that incorporates 1 nmole of dTTP into acid-insoluble products in 10 min at 37°C using poly(A)-oligo(dT) as template primer.

## Quality Control

- Nuclease Contamination Assay: Nuclease activity is not detected after incubation of 1 µg of DNA and RNA with 200 U of RocketScript™ Reverse Transcriptase at 37-42°C for 4 hrs.

## Storage

Store at -20°C. If stored in the recommended temperature, this product will be stable until the expiration date printed out on the label.

## Online Resources



Korean



English

Visit our [product page](#) for additional information and protocols.

## Ordering Information

Description	Cat. No
RocketScript™ Reverse Transcriptase	10,000 U (50 rxn) E-3141
	50,000 U (250 rxn) E-3142

## Notice

BIONEER corporation reserves the right to make corrections, modifications, improvements and other changes to its products, services, specifications or product descriptions at any time without notice.

## Explanation of Symbols



Batch Code



Catalog Number



Caution



Consult Instructions For Use



Contains Sufficient for <n> tests



Research Use Only





Temperature Limitation



Use-by Date

**Experimental Procedures**

Steps		Procedure Details																																																																					
1	 <b>Preparation of reaction mixture</b>	<p>1. Mix template RNA and primers in a sterile tube (not provided) indicated as below.</p> <ul style="list-style-type: none"> <li>Amount of template RNA and primers</li> </ul> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th colspan="2">Components</th> <th>20 µl reaction</th> </tr> </thead> <tbody> <tr> <td rowspan="2" style="text-align: center;">Template RNA</td> <td style="text-align: center;">Total RNA</td> <td style="text-align: center;">1 µg</td> </tr> <tr> <td style="text-align: center;">RNA</td> <td style="text-align: center;">5-100 ng</td> </tr> <tr> <td rowspan="2" style="text-align: center;">Primers</td> <td style="text-align: center;">Oligo dT or Random primer</td> <td style="text-align: center;">10-100 pmol</td> </tr> <tr> <td style="text-align: center;">Sequence specific</td> <td style="text-align: center;">10-30 pmol</td> </tr> </tbody> </table> <p>2. Add all components for cDNA synthesis into PCR tubes (not provided) to a total volume of 20 µl.</p> <ul style="list-style-type: none"> <li>Preparation of reaction mixture</li> </ul> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th colspan="2">Components</th> <th>20 µl reaction</th> </tr> </thead> <tbody> <tr> <td colspan="2">Mixture of template RNA and primers</td> <td style="text-align: center;">Variable</td> </tr> <tr> <td colspan="2">5X Reaction buffer</td> <td style="text-align: center;">4 µl</td> </tr> <tr> <td colspan="2">10 mM dNTPs</td> <td style="text-align: center;">2 µl or Variable</td> </tr> <tr> <td colspan="2">100 mM DTT</td> <td style="text-align: center;">2 µl</td> </tr> <tr> <td colspan="2">RNase inhibitor (100 ng/µl)</td> <td style="text-align: center;">0.5-1 µl</td> </tr> <tr> <td colspan="2">RocketScript™ Reverse Transcriptase (200 U/µl)</td> <td style="text-align: center;">200 U</td> </tr> <tr> <td colspan="2">Nuclease-free water</td> <td style="text-align: center;">Variable</td> </tr> <tr> <td colspan="2">Total volume</td> <td style="text-align: center;">20 µl</td> </tr> </tbody> </table>	Components		20 µl reaction	Template RNA	Total RNA	1 µg	RNA	5-100 ng	Primers	Oligo dT or Random primer	10-100 pmol	Sequence specific	10-30 pmol	Components		20 µl reaction	Mixture of template RNA and primers		Variable	5X Reaction buffer		4 µl	10 mM dNTPs		2 µl or Variable	100 mM DTT		2 µl	RNase inhibitor (100 ng/µl)		0.5-1 µl	RocketScript™ Reverse Transcriptase (200 U/µl)		200 U	Nuclease-free water		Variable	Total volume		20 µl																													
Components		20 µl reaction																																																																					
Template RNA	Total RNA	1 µg																																																																					
	RNA	5-100 ng																																																																					
Primers	Oligo dT or Random primer	10-100 pmol																																																																					
	Sequence specific	10-30 pmol																																																																					
Components		20 µl reaction																																																																					
Mixture of template RNA and primers		Variable																																																																					
5X Reaction buffer		4 µl																																																																					
10 mM dNTPs		2 µl or Variable																																																																					
100 mM DTT		2 µl																																																																					
RNase inhibitor (100 ng/µl)		0.5-1 µl																																																																					
RocketScript™ Reverse Transcriptase (200 U/µl)		200 U																																																																					
Nuclease-free water		Variable																																																																					
Total volume		20 µl																																																																					
2	 <b>cDNA synthesis</b>	<p>3. Perform the reaction under the following conditions.</p> <p>3-1. Cyclic reverse transcription (example 1)</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th rowspan="2">Step</th> <th colspan="4">Temperature</th> <th rowspan="2">Time</th> <th rowspan="2">Cycles</th> </tr> <tr> <th>dN<sub>6</sub></th> <th>dN<sub>12</sub></th> <th>dT<sub>20</sub></th> <th>Sequence specific</th> </tr> </thead> <tbody> <tr> <td>Primer annealing</td> <td style="text-align: center;">15°C</td> <td style="text-align: center;">30°C</td> <td style="text-align: center;">37°C</td> <td style="text-align: center;">T<sub>m</sub> of primers</td> <td style="text-align: center;">10-30 sec</td> <td rowspan="3" style="text-align: center;">12 cycles or less</td> </tr> <tr> <td>cDNA synthesis</td> <td colspan="3"></td> <td style="text-align: center;">50°C</td> <td style="text-align: center;">4 min</td> </tr> <tr> <td>Melting secondary structure &amp; cDNA synthesis</td> <td colspan="3"></td> <td style="text-align: center;">55°C</td> <td style="text-align: center;">30 sec</td> </tr> <tr> <td>Heat inactivation</td> <td colspan="3"></td> <td style="text-align: center;">95°C</td> <td style="text-align: center;">5 min</td> <td style="text-align: center;">1 cycle</td> </tr> </tbody> </table> <p>3-2. Single temperature reaction (example 2)</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th rowspan="2">Step</th> <th colspan="4">Temperature</th> <th rowspan="2">Time</th> <th rowspan="2">Cycles</th> </tr> <tr> <th>dN<sub>6</sub></th> <th>dN<sub>12</sub></th> <th>dT<sub>20</sub></th> <th>Sequence specific</th> </tr> </thead> <tbody> <tr> <td>Primer annealing</td> <td style="text-align: center;">15°C</td> <td style="text-align: center;">30°C</td> <td style="text-align: center;">37°C</td> <td style="text-align: center;">T<sub>m</sub> of primers</td> <td style="text-align: center;">1 min</td> <td style="text-align: center;">1 cycle</td> </tr> <tr> <td>cDNA synthesis</td> <td colspan="3"></td> <td style="text-align: center;">42-70°C</td> <td style="text-align: center;">10-60 min</td> <td style="text-align: center;">1 cycle</td> </tr> <tr> <td>Heat inactivation</td> <td colspan="3"></td> <td style="text-align: center;">95°C</td> <td style="text-align: center;">5 min</td> <td style="text-align: center;">1 cycle</td> </tr> </tbody> </table> <p>* <b>Note:</b> Primer annealing step can be omissible. Perform cDNA synthesis reaction as follows: cDNA synthesis, 50°C, 60 min; Heat inactivation, 95°C, 5 min. RT reaction temperature should be selected to fit the T<sub>m</sub> value of primers.</p> <p>4. After the reaction, maintain the reaction mixture at 4°C. The samples can be stored at -20°C until use.</p>	Step	Temperature				Time	Cycles	dN <sub>6</sub>	dN <sub>12</sub>	dT <sub>20</sub>	Sequence specific	Primer annealing	15°C	30°C	37°C	T <sub>m</sub> of primers	10-30 sec	12 cycles or less	cDNA synthesis				50°C	4 min	Melting secondary structure & cDNA synthesis				55°C	30 sec	Heat inactivation				95°C	5 min	1 cycle	Step	Temperature				Time	Cycles	dN <sub>6</sub>	dN <sub>12</sub>	dT <sub>20</sub>	Sequence specific	Primer annealing	15°C	30°C	37°C	T <sub>m</sub> of primers	1 min	1 cycle	cDNA synthesis				42-70°C	10-60 min	1 cycle	Heat inactivation				95°C	5 min	1 cycle
Step	Temperature				Time	Cycles																																																																	
	dN <sub>6</sub>	dN <sub>12</sub>	dT <sub>20</sub>	Sequence specific																																																																			
Primer annealing	15°C	30°C	37°C	T <sub>m</sub> of primers	10-30 sec	12 cycles or less																																																																	
cDNA synthesis				50°C	4 min																																																																		
Melting secondary structure & cDNA synthesis				55°C	30 sec																																																																		
Heat inactivation				95°C	5 min	1 cycle																																																																	
Step	Temperature				Time	Cycles																																																																	
	dN <sub>6</sub>	dN <sub>12</sub>	dT <sub>20</sub>	Sequence specific																																																																			
Primer annealing	15°C	30°C	37°C	T <sub>m</sub> of primers	1 min	1 cycle																																																																	
cDNA synthesis				42-70°C	10-60 min	1 cycle																																																																	
Heat inactivation				95°C	5 min	1 cycle																																																																	