

[Cat. No.] Please refer to the **Ordering Information**

Introduction

HotStart Taq DNA Polymerase is designed to perform Hotstart PCR. The HotStart Taq DNA Polymerase is inhibited at lower temperature, but is fully activated at temperature above 70°C. This prevents the formation of mis-primed products and primer-dimers during the reaction setup process, resulting in improved PCR specificity.

Applications

- Real-time quantification of DNA and cDNA targets using dual probe and dsDNA binding dye
- HotStart PCR
- Multiplex PCR
- Automated PCR

Components

Components	E-2017	E-2017-2	E-2017-1
HotStart Taq DNA Polymerase	250 U (50 µl)	500 U (50 µl x 2)	1,000 U (50 µl x 4)
10X Reaction buffer	0.5 ml	0.5 ml x 2	0.5 ml x 4
10 mM dNTPs	0.5 ml	0.5 ml x 2	0.5 ml x 4
Dilution buffer	0.5 ml	0.5 ml x 2	0.5 ml x 4

Components	E-2017-3	E-2017-4
HotStart Taq DNA Polymerase	250 U (50 µl)	1,000 U (50 µl x 4)
10X Reaction buffer	0.5 ml	0.5 ml x 4
10 mM dNTPs	-	-
Dilution buffer	0.5 ml	0.5 ml x 4

* Note: For research use only. Not for use in diagnostic or therapeutic procedures.

Specifications

HotStart Taq DNA Polymerase	
5' to 3' exonuclease activity	Yes
3' to 5' exonuclease activity	No
3'-A overhang	Yes
Fragment size	Up to 10 kb

Buffer Composition

10X Reaction buffer	100 mM Tris-HCl, 450 mM KCl, 15 mM MgCl ₂ , pH 9.0
Dilution buffer	20 mM Tris-HCl, 0.5 mM EDTA, 1 mM DTT, 100 mM KCl, stabilizer, 50% glycerol, pH 8.0

Storage Buffer

HotStart Taq DNA Polymerase is supplied in 50% (v/v) glycerol containing 20 mM Tris-HCl, 0.5 mM EDTA, 1 mM DTT, 100 mM KCl, and stabilizer, pH 8.0.

Unit Definition

One unit is defined as the amount of enzyme that incorporates 10 nmole of dNTP into acid-insoluble products in 30 min at 72°C.

Quality Control

- Nuclease Contamination Assay: Nuclease activity is not detected after incubation of 1 µg of substrate DNA (supercoiled plasmid and Lambda/Hind III DNA) with 5 U of HotStart Taq DNA Polymerase in 50 µl reaction volume at 37°C and 70°C for 18 hrs.

Storage

Store at -20°C. If stored in the recommended temperature, this product will be stable until the expiration date printed out on the label.

Online Resources



Korean



English

Visit our **product page** for additional information and protocols.

Ordering Information

Description	Cat. No
250 U	10 mM dNTP E-2017
500 U	10 mM dNTP E-2017-2
1,000 U	10 mM dNTP E-2017-1
	E-2017-4





Notice

BIONEER corporation reserves the right to make corrections, modifications, improvements and other changes to its products, services, specifications or product descriptions at any time without notice.

Explanation of Symbols

LOT Batch Code	REF Catalog Number	Caution	Consult Instructions For Use
Contains Sufficient for <n> tests	Research Use Only	Temperature Limitation	Use-by Date

Experimental Procedures

Steps		Procedure Details																											
1	 Thaw reagents	<p>1. Thaw 10X Reaction buffer and 10 mM dNTPs on ice and mix thoroughly before use. Then, briefly spin down all components including template DNA, HotStart Taq DNA Polymerase, primers and nuclease-free water.</p>																											
2	 Preparation of reaction mixture	<p>2. Add all components into PCR tubes (not provided) to a total volume of 20 µl or 50 µl.</p> <ul style="list-style-type: none"> Preparation of reaction mixture <table border="1"> <thead> <tr> <th>Components</th> <th>20 µl reaction</th> <th>50 µl reaction</th> </tr> </thead> <tbody> <tr> <td>Template DNA*</td> <td>Variable</td> <td>Variable</td> </tr> <tr> <td>Forward primer (10 pmol/µl)</td> <td>1-2 µl</td> <td>2.5-5 µl</td> </tr> <tr> <td>Reverse primer (10 pmol/µl)</td> <td>1-2 µl</td> <td>2.5-5 µl</td> </tr> <tr> <td>10X Reaction buffer</td> <td>2 µl</td> <td>5 µl</td> </tr> <tr> <td>10 mM dNTPs</td> <td>2 µl or Variable</td> <td>5 µl or Variable</td> </tr> <tr> <td>HotStart Taq DNA Polymerase (5 U/µl)</td> <td>0.5-1 U</td> <td>1-2.5 U</td> </tr> <tr> <td>Nuclease-free water</td> <td>Variable</td> <td>Variable</td> </tr> <tr> <td>Total volume</td> <td>20 µl</td> <td>50 µl</td> </tr> </tbody> </table> <p>* Recommended amounts of template DNA is as follows; plasmid and lambda DNA, > 1 pg; bacterial genomic DNA, > 10 pg; human genomic DNA, > 1 ng.</p> <p>3. Mix the reaction mixture by tapping or pipetting, and briefly spin down.</p>	Components	20 µl reaction	50 µl reaction	Template DNA*	Variable	Variable	Forward primer (10 pmol/µl)	1-2 µl	2.5-5 µl	Reverse primer (10 pmol/µl)	1-2 µl	2.5-5 µl	10X Reaction buffer	2 µl	5 µl	10 mM dNTPs	2 µl or Variable	5 µl or Variable	HotStart Taq DNA Polymerase (5 U/µl)	0.5-1 U	1-2.5 U	Nuclease-free water	Variable	Variable	Total volume	20 µl	50 µl
Components	20 µl reaction	50 µl reaction																											
Template DNA*	Variable	Variable																											
Forward primer (10 pmol/µl)	1-2 µl	2.5-5 µl																											
Reverse primer (10 pmol/µl)	1-2 µl	2.5-5 µl																											
10X Reaction buffer	2 µl	5 µl																											
10 mM dNTPs	2 µl or Variable	5 µl or Variable																											
HotStart Taq DNA Polymerase (5 U/µl)	0.5-1 U	1-2.5 U																											
Nuclease-free water	Variable	Variable																											
Total volume	20 µl	50 µl																											
3	 Incubate reactions in a thermal cycler	<p>4. Perform the reaction under the following conditions.</p> <table border="1"> <thead> <tr> <th>Step</th> <th>Temperature</th> <th>Time</th> <th>Cycles</th> </tr> </thead> <tbody> <tr> <td>Pre-denaturation</td> <td>94°C</td> <td>15 min</td> <td>1 cycle</td> </tr> <tr> <td>Denaturation</td> <td>94°C</td> <td>0.5-1 min</td> <td></td> </tr> <tr> <td>Annealing</td> <td>45-65°C†</td> <td>0.5-1 min</td> <td>25-35 cycles</td> </tr> <tr> <td>Extension</td> <td>72°C</td> <td>1 min/kb</td> <td></td> </tr> <tr> <td>Final extension</td> <td>72°C</td> <td>5-10 min</td> <td>1 cycle</td> </tr> </tbody> </table> <p>* Note: For maximum yield and specificity, temperatures and cycling times should be optimized for each new template DNA or primers. † Set the annealing temperature to 3-5°C lower than the T_m of the primers.</p>	Step	Temperature	Time	Cycles	Pre-denaturation	94°C	15 min	1 cycle	Denaturation	94°C	0.5-1 min		Annealing	45-65°C†	0.5-1 min	25-35 cycles	Extension	72°C	1 min/kb		Final extension	72°C	5-10 min	1 cycle			
Step	Temperature	Time	Cycles																										
Pre-denaturation	94°C	15 min	1 cycle																										
Denaturation	94°C	0.5-1 min																											
Annealing	45-65°C†	0.5-1 min	25-35 cycles																										
Extension	72°C	1 min/kb																											
Final extension	72°C	5-10 min	1 cycle																										
4	 Analyze with gel electrophoresis	<p>5. After the reaction, maintain the reaction mixture at 4-8°C. The samples can be stored at -20°C until use.</p> <p>6. Load samples on agarose gel with adding a loading-dye mixture, and perform gel electrophoresis for analysis.</p>																											