ExiProgen[™] **EC-Tagfree Protein Synthesis Kit**

Cat. No. K-7320

***ExiProgen*[™] EC-Tagfree Protein Synthesis Kit**

User Guide

K-7320

Version No.: 3 (2022-06-23)

Please read all the information in booklet before using the unit

BIONEER Corporation

**Bioneer Global Center, 71, Techno-2-ro,
Yuseong-gu, Daejeon, 34013, Republic of Korea**

Tel: 1588-9788

Email: sales@bioneer.co.kr

www.bioneer.com

Intended Use

ExiProgen[™] EC-Tagfree Protein Synthesis Kit is developed and supplied for research purposes only. Certain applications possible with this kit may require special approval by appropriate local and/or national regulatory authorities in the country of use.

Safety Warning and Precaution

Wear appropriate protection when handling any irritant or harmful reagents. The use of a laboratory coat, protective gloves and safety goggles are highly recommended. For more information, please consult the appropriate Material Safety Data Sheet (MSDS).

Warranty and Liability

All BIONEER products undergo extensive Quality Control testing and validation. BIONEER guarantees quality during the warranty period as specified, when following the appropriate protocol as supplied with the product. It is the responsibility of the purchaser to determine the suitability of the product for its particular use. Liability is conditional upon the customer providing full details of the problem to Bioneer within 30 days.

Quality Management System ISO 9001 Certified

Every aspect of our quality management system from product development, production to quality assurance and supplier qualification meets the world-class standards.

Patent

ExiProgen[™] and its kits are protected by the patents KR10-2011-0085824, PCT/KR2012/006715, KR10-2012-0138335, and KR10-2013-0022324.

Trademark

ExiProgen[™] is a trademark of BIONEER Corporation.

Copyright

Copyright 2022. BIONEER Corporation. All Rights Reserved.

Notice

BIONEER corporation reserves the right to make corrections, modifications, improvements and other changes to its products, services, specifications or product descriptions at any time without notice.

Contents

Product Information	1
Components	1
Storage	1
Specifications.....	2
Precautions	2
Introduction	3
Cell-free Protein Synthesis System	3
The <i>ExiProgen</i> [™] Protein Synthesis System	5
Product Description	7
Experimental Procedures	9
Components Details	9
Experiments Preparation	12
Protocols	16
Analysis of Sample	21
Maintenance	24
Troubleshooting	25
Appendix A: pBIVT Vectors	27
Appendix B: Screening of Template DNA Concentration	28
References	29
Ordering Information	31
Related Products	31
Explanation of Symbols	32

Product Information

Components

ExiProgen™ EC-Tagfree Protein Synthesis Kit ①		
Components	K-7320 (8 rxn)	Storage
Cartridge ①	96-well x 1 ea	4°C
Dialysis tube	1 pack (16 ea/pack)	
Disposable filter tip	1 pack (8 ea/pack)	-
Protection Cover	1 ea	-

ExiProgen™ EC-Tagfree Protein Synthesis Kit ②		
Components	K-7320 (8 rxn)	Storage
Cartridge ②	96-well x 1 ea	-20°C (<i>E. coli</i> extract recommended at -70°C or less)
<i>E. coli</i> extract	8-tube strip (Yellow) x 1 ea	
Master mix	8-tube strip (Violet) x 1 ea	
DEPC DW	8-tube strip (White) x 1 ea	
Storage buffer	35 ml x 2 bottle	
TEV protease	1.5 ml tube x 1 ea	
Positive control DNA	1.5 ml tube x 1 ea	

Storage

The Cartridge ① and Dialysis tube included in Kit ① should be stored at 4°C. The Cartridge ②, Master mix, DEPC DW, Storage buffer, TEV protease, and Positive control DNA included in Kit ② should be stored at -20°C. In particular, *E. coli* extract contains T7 RNA polymerase and ribosome, so it is recommended to store at a temperature below -70°C immediately upon receipt, and if the extract repeatedly frozen and thawed, the efficiency may decrease.

Specifications

ExiProgen™ EC-Tagfree Protein Synthesis Kit	
	K-7320 (8 rxn)
Number of reactions (Up to 8 rxn/ 1 times)	8 rxn
Operating time	26 hours
Required sample DNA	1-6 µg [Plasmid DNA: 0.6 µg/kb (Vector)]
Yield	Up to 200 µg*

* **Note:** The protein yield can be varied depending on the type of target protein.

Precautions

The Cartridge ① and ② of this Synthesis Kit are each covered with sealing film in order to prevent cross-contamination, evaporation, or leakage of solutions inside. All of the plastic products and buffers in this kit are provided under nuclease- and protease-free condition, hence, please be careful not to contaminate any part of the kit with nuclease or protease.

Introduction

Cell-free Protein Synthesis System

Since protein composes essential components in biological reactions, such as enzymes, hormones, and structural proteins, researches on roles and structures of protein are being actively conducted in the post-genomic era. These studies start with producing a specific type of protein.

The various protein expression systems are based on different cell, such as *E. coli*, yeast, animal, and plants cells, which are used for production of recombinant proteins. In order to express proteins based on these systems, a vector containing target genes is transformed into cells and the cell line is cultured to express a large amount of protein. After that, cells are usually disrupted and proteins are purified from cell lysates. However, the general protein expression systems require much time and labor, as these have to go through a series of steps such as selection of: 1) strains, 2) cell lines stably expressing recombinant proteins, 3) cell culture, 4) cell disruption, and 5) protein purification. In addition, if the target protein is toxic to strains and cell lines, it takes several weeks or up to several months to synthesize a single protein because it is difficult to express a toxic protein and requires optimization of various conditions for its expression.

To overcome the limitations of cell-based system, cell-free protein synthesis method and related products have been developed (Keum et al., 2009; Kim et al., 2007; Kim et al., 2011; Park & Hamad-Schifferli, 2010). Cell-free protein synthesis system is a technology for protein synthesis based on cell extract by transcribing and translating *in vitro*. To carry out protein synthesis *in vitro*, the following should be prepared: 1) cell extract including T7 RNA polymerase, ribosome, tRNA, and enzyme, and 2) protein expression solution containing a substrate and a substance that performs transcription and translation such as amino acids, rNTPs, and energy sources. After that, a gene encoding the target protein is added, and the reaction is conducted at an appropriate temperature and time to express the recombinant protein (Figure 1). In case of cell-free protein synthesis systems, the reaction duration is at least 1 to 3 hours, which can dramatically reduce the time taken for protein expression. In addition, it is easier to express proteins that show toxicity in cells because the expression is based on cell extracts, not living cells, and since it is an open system, conditions for protein expression can be easily adjusted (Hino et al., 2008).

Cell-free protein synthesis systems are widely used in enzyme engineering, protein labeling, research on protein-protein interaction mechanism, and research on protein active site (Forster et

al., 2004; Josephson et al., 2005; Keum et al., 2006; Kigawa et al., 2002; Ohashi et al., 2007; Villemagne et al., 2006).


Figure 1. Principles of cell-free protein synthesis system.

The *ExiProgen*[™] Protein Synthesis System

BIONEER has developed *ExiProgen*[™], a fully automatic protein synthesis instrument that applies cell-free protein synthesis technology and affinity purification technology using magnetic nanoparticles. Therefore, high-purity proteins can be obtained quickly and fully automatically by using our protein expression and purification kits. Furthermore, the instrument *ExiProgen*[™] is highly versatile as it can fully automatically extract the required DNA and RNA from various samples using nucleic acid extraction kit series.

Our protein expression and purification kits use T7 expression system of *E. coli* and can be used immediately with template DNA (Ahn et al., 2008; Rungpragayphan et al., 2003). Protein synthesis using *ExiProgen*[™] allows the addition of template DNA to perform a series of processes, including protein expression and purification, to finally obtain samples containing target proteins.

Protein expression and purification kits are divided into three product groups (Figure 2): 1) Template preparation kits that produce template DNA used for cell-free protein expression, 2) *AccuRapid*[™] Kit Series that is responsible for easy manual expression or synthesis of recombinant proteins, and 3) *ExiProgen*[™] Kit Series that is responsible for an automatic expression or synthesis of recombinant proteins. Each protein expression and purification kit can express and synthesize proteins from µg up to mg scale.

In addition, there are other products that can be selected depending on the purpose of the experiment, such as a kit for improving protein expression efficiency, a kit solely for protein purification or buffer exchange, etc.

Please refer to the BIONEER's website (www.bioneer.com) for a detailed description of each kit.


Figure 2. Product information related to protein expression and synthesis.

Product Description

ExiProgen[™] EC-Tagfree Protein Synthesis Kit can be applied to *ExiProgen*[™], a fully automated protein synthesis and nucleic acid extraction system developed by BIONEER, and up to 200 µg of 8 different histidine tag-free proteins can be obtained with high purity in a single reaction.

This product uses Stepwise Exchange Cell-Free (SECF) protein synthesis technology to create an optimal environmental for protein expression. During dialysis using a membrane, by-products of protein expression are removed, and at the same time, energy source and amino acids are continuously supplied to the reaction mixture (Figure 3) (Kim & Choi, 1996; Lim & Kim, 2019). The optimized protein expression environment leads to high protein synthesis yield.

The expressed protein is purified by affinity method using histidine tag. After purification, TEV protease acts on the TEV cleavage site that connected the target protein and histidine tag, finally obtaining histidine tag-free protein (Figure 4) (Li et al., 2013).

All reactions including protein expression/purification and buffer exchange through dialysis are carried out fully automated. Therefore, the user only needs to add the target DNA, and 8 different types of final target proteins can be recovered in a dissolved state in a storage buffer.


Figure 3. SECF Principles of *ExiProgen*[™] EC-Tagfree Protein Synthesis Kit.


Figure 4. Principles of *ExiProgen*™ EC-Tagfree Protein Synthesis Kit.

Experimental Procedures

Components Details

1. Cartridges

F2					0.9 ml
F1					0.3 ml
E2					0.9 ml
E1					0.9 ml
D2					0.9 ml
D1					0.9 ml
C2					0.9 ml
C1					0.9 ml
B2					0.9 ml
B1					0.9 ml
A2					0.25 ml
A1					1 ml

Ni-NTA magnetic bead

Binding/Washing buffer

Buffer to set equilibrium so that expressed protein sample can be bound onto Ni-NTA magnetic bead, and after binding of the target protein, used for washing process to remove impurities.

Sterile distilled water

Used for tip cleaning before elution.

L2					1.2 ml
L1					1 ml
K2					1 ml
K1					1 ml
J2					1 ml
J1					1 ml
I2					1 ml
I1					1 ml
H2					1 ml
H1					1 ml
G2					1 ml
G1					0.3 ml

Feeding buffer

Protein expression solution for SECF reaction and contains NTP, amino acids, and energy sources.

Storage buffer

After protein purification, the buffer for storing the final target proteins.

TEV protease and reaction buffer

(When adding 100 µl of TEV protein, final volume)

Used for removal of 6x histidine-tag.

2. *E. coli* extract

E. coli extract included in Kit ② is a cell lysate and supplies T7 RNA polymerase, ribosome, tRNA, etc., required for protein expression. It is provided as an 8-tube strip and **130 µl** is dispensed in each tube.

3. Master mix

Master mix, a components of Kit ②, contains NTPs, amino acids, energy source and salt necessary for protein expression. It is supplied as an 8-tube strip, and **220 µl** is dispensed in each tube.

4. DEPC DW

DEPC DW is used to adjust the final volume of the protein expression solution, and it is supplied as an 8-tube strip, and **130 µl** is dispensed in each tube.

5. Dialysis tube

Dialysis tube, which is specifically manufactured by BIONEER, is used to continuously supply the energy source from the feeding buffer to the reaction solution during protein expression, and to exchange the purification buffer into storage buffer after purification of a target protein. It is provided in 20% ethanol and each tube is individually packaged, so as many can be used as needed according to the need. Before use, it is required to remove 20% ethanol and rinse it out with sterile distilled water for once.

6. Storage buffer

The composition of the storage buffer provided in the kit is as follows.

- Buffer composition: 50 mM Tris-HCl, 100 mM NaCl, 1 mM DTT, 0.1 mM EDTA, 0.05% (v/v) NaN₃, 50% (v/v) glycerol, pH 7.6


Storage buffer is provided in the bottle, and just before use, dispense 1 ml to each of the column of rows G2-J1 of Cartridge ②. If you want to use a buffer with a different composition, you can manufacture it according to the required composition. However, 10% or higher concentration of glycerol must be added to the buffer. If not, a large amount of storage buffer will flow into the Dialysis tube due to the concentration difference between the solutions, leading to the overflow of the solution containing the purified protein out of the Reaction block.

7. TEV protease

TEV protease is used to remove 6x histidine-tag, and just before use, **dispense 100 µl to each column of the row G1 of Cartridge ②.**

3. Note

Structure of *ExiProgen*™ Baseplate

Experiments Preparation

Preparation of Template DNA

The template DNA must have “T7 promoter - Ribosome binding site (RBS) – 6x histidine tag – **TEV cleavage site** - Target gene - T7 terminator” structure to apply for ExiProgen™ EC-Tagfree Protein Synthesis Kit (Figure 5). Proteins synthesized based on this sequence will have an additional Glycine (or Serine) amino acid in front of the target protein after the purification process. The position X can be any of 20 amino acids, while the DNA sequence provided as a positive control encodes Glu-Asn-Leu-Tyr-Phe-Gln-Gly.


Figure 5. Template DNA structures and TEV site sequences for ExiProgen™ EC-Tagfree Protein Synthesis Kit.

For cloning of a target gene, pBIVT Vector Set-1 (Cat. No. K-7350, see Appendix A.) can be used as an expression vector, and pK7, pIVEX, pET vector can be used elsewhere.

DNA sequences must be optimized in order for *E. coli* codon to maximize protein synthesis, and if it is difficult to optimize codon, template DNA (vector) can be synthesized through our Gene Synthesis Service (Trotta, E., 2011). In addition, since it is provided as cloned into our *in vitro* translation vector, pBT7-N-His vector and pBT7-C-His vector, it can be applied quickly and easily. For more information, visit the website (www.bioneer.com).

Example of Template DNA Preparation

1. When using user's expression vector

After constructing a PCR product containing TEV cleavage site* between the 6x histidine tag and the target gene, clone it into the desired expression vector.

* **Example:** GAAAACTTATATTTTCAGGGC, Glu-Asn-Leu-Tyr-Phe-Gln-Gly

2. When using Positive control DNA as an expression vector

Prepare the primers as follows to obtain the PCR product.

Forward primer:

5'- CGCCTACGGATCCGAAAACTTATATTTTCAGGGCXXXXXXXXXXXXXXXXXXXX-3'
*Bam*H I TEV site Overlapping sequences to target gene (18-mer)


Reverse primer:

5'- CGCCTACGCGGCCGCXXXXXXXXXXXXXXXXXXXX-3'
Not I Overlapping sequences to target gene (18-mer)

The prepared PCR product can be cloned into positive control DNA.

Before You Begin

1. Take out the Cartridge ② from the Kit ② box and thaw at room temperature.
2. Take out *E. coli* extract, Master mix, DEPC DW according to the number of reaction samples from the Kit ② box and thaw it on ice. If necessary, take out the Positive control DNA tube and thaw it together.
 - * **Note:** It takes about 2 hours for the solution to completely thaw. Before use, make sure that all the solutions are completely thawed.
 - * **Note:** The Positive control is provided with pBIVT-(TEV site)-AcGFP and is approximately 3.8 kb (28 kDa).
3. Prepare the Reaction block (ExiProgen™ accessory).
4. Take out twice as many Dialysis tube as the number of samples from Kit ① box. After taking out the Dialysis tube with tweezers, use a pipette or shake it up and down to remove the ethanol inside. Then, rinse inside and out of the tube using a squeeze bottle containing sterile distilled water.
 - * **Note:** Sterile distilled water is not provided, so please prepare it.
5. Completely remove the water inside the Dialysis tube using a pipette, and install it on the Reaction block as shown below. Fill the tube in row B with **500 µl of sterile distilled water**.

6. Take out Cartridge ①, Disposable filter tip, and Protection Cover from the Kit ① box and prepare it.

7. Preparation of protein expression solution

Before starting protein synthesis with *ExiProgen*[™], the protein expression solution containing the DNA to be expressed must be prepared.

- 1) Prepare protein expression solution using pre-thawed *E. coli* extract, Master mix, and DEPC DW. After spin down each tube, mix the solution with a pipette before use.

Components	Sample	Positive control DNA
Template DNA	X µl	6 µl
<i>E. coli</i> extract	120 µl	120 µl
Master mix	210 µl	210 µl
DEPC DW	(120-X) µl	114 µl
Total volume	450 µl	450 µl


* **Note:** Amount of template DNA – Plasmid DNA: 0.6 µg/kb (Vector)

Example) When using Plasmid DNA containing insert DNA with a size of 5 kb and a concentration of 100 ng/µl
→ Use 30 µl (= 3 µg)

* **Note:** For optimal protein synthesis, please refer to Appendix B. and conduct screening of template DNA concentration in advance.




- 2) Add each **protein expression solution** prepared above into the tubes in **row A** of the Reaction block.


Example) For 1 sample

- 3) Finish the preparation for the experiment by installing the Protection Cover on the Reaction block.


Protocols

Steps	Procedure Details
Example 1) For 1 sample  Example 2) For 8 samples 	1. Punch holes in the sealing films of Cartridge ① and ② using 6 Hole Punch (<i>ExiProgen</i>™'s accessory) according to the number of samples. * Note: Refer to the left pictures and punch holes depending on the required number of samples. 2. Add 1 ml of storage buffer to each column of the rows G2-J1 and 100 µl of TEV protease to the columns of the row G1 of Cartridge ②. * Note: If you want to use a buffer with a different composition, you can manufacture it according to the required composition. However, 10% or higher concentration of glycerol must be added to the buffer. If not, a large amount of storage buffer will flow into the Dialysis tube due to the concentration difference between the solutions, leading to the overflow of the solution containing the purified protein out of the Reaction block.
	3. Open the door of <i>ExiProgen</i>™ and pull out the baseplate completely.
	4. Load Cartridge ② in the position of ② on the baseplate. * Note: Place the row L of Cartridge ② in the direction of the Heating block first and then install it, and afterwards check if it is properly fixed.

5. Load Cartridge ① in the position of ① on the baseplate.

*** Note:** There are silicon rings embedded on both sides of the Cartridge ① installation position. Place the left side first, and then the right side.

6. Install Waste Tray between Cartridge ① and ②.

*** Note:** Follow the exact installation order of Cartridge ② → Cartridge ① → Waste Tray. Ensure that Cartridges and Waste Tray are firmly installed, without showing any movement.

7. Install Reaction block prepared in pp. 14-15 to the Heating and cooling block (Magnetic part) of the baseplate.

*** Note:** Ensure that the row A of the Reaction block is located towards the inside of the instrument and the row B towards the Disposable Tip Rack.


8. Place Disposable filter tips in **row B** of Disposable Tip Rack according to the number of samples.

*** Note:** Tips should be placed in the corresponding positions with the punched holes of the Cartridges. Do not insert tips in the corresponding positions of the cartridge that are not punched.


9. Finally, confirm holes in the cartridges and position of samples and tips (Refer to the pictures on page 18). Push the baseplate completely until you hear the click sound, then close the door.

*** Note**

Example 1) For 1 Sample

Example 2) For 8 Samples

10. Turn on the *ExiProgen™* and tap 'Press to start' button. When the 'Press to start' button is pressed, the *ExiProgen™* screen appears as shown in the lower left, and after the scroll bar moves, it moves to the next screen.

*** Note:** This process initializes the X, Y, Z axis values of the equipment. If the next screen does not proceed normally, turn off the power of the equipment and contact the A/S center.

11. In the MENU screen, press 'Start' button to select a proper protocol.


12. The PREP SETUP screen appears as shown in the left, and a screen to select the protocol number for each kit appears. Press '**904**' to confirm that the following information is displayed on the screen, and then press the 'Enter' button.

Prep Type: Protein

Sample SRC: Synthesis_TF

13. After selecting the protocol, the screen to select the Elution volume appears. It is used for nucleic acid extraction, and when using this kit, please press the 'ok' button to move next steps.

14. After the Elution volume selection screen, the screen to select protein synthesis temperature appears. When using this kit, please select **30°C** and press the 'ok' button.

15. When the 'CHECK LIST' screen appears, confirm once again that cartridges and other accessories are fixed to each location properly, and then click the 'ok' button.

16. When the 'Running Mode' screen appears as shown on the left, finally check the following information and click the 'RUN' button. After this, the instrument starts and it takes approximately **26 hours** to finish run.

Prep Type: Protein

Sample SRC: Synthesis_TF

17. 'Work Completion' screen appears when the protocol is completed. Remove all components used in the experiment, and choose a desired button.

* **Note:** If you want to quit and press the 'ok' button, the UV lamp will be turned on automatically.

Analysis of Sample

After protein synthesis with *ExiProgen*™, the final target proteins of about 180-190 µl is collected from Dialysis tube located at the **row B** of the Reaction block (Figure 6). Final protein samples may contain Ni-NTA magnetic beads, but this can be removed by centrifugation before use.

If necessary, reaction samples from protein expression and purification can be recovered from the rows K1, K2, and L2 of Cartridges ② as shown in Figure 6.


Figure 6. Samples for each row of Cartridge ②.

- Expression sample: Samples after expression, but without purification.
- Unbound sample: Samples not bound to Ni-NTA magnetic bead.
- Bead sample: Used bead samples for purification.

Each sample can be checked through SDS-PAGE whether target protein is synthesized as desired.

1. Prepare the loading mixture as shown in the table.

Components	Expression/ Unbound sample	Purified protein sample	Bead sample*
Sample	5 µl	5 µl	15 µl
4X loading dye	5 µl	5 µl	5 µl
Sterile distilled water	10 µl	10 µl	-
Total volume	20 µl	20 µl	20 µl

* **Note:** Add 360-400 µl of sterile distilled water to the well containing the beads for suspension and proceed with the sampling afterwards.

2. Incubate the prepared loading mixture at 95°C for 5-10 min.

3. Load the following amounts of each sample to the wells of 10% or 12% SDS-PAGE gel [10 x 8 (cm), 10-well] and run the SDS-PAGE.

- Expression/Unbound samples: 5 µl/well
- Purified protein/Bead sample: 10 µl/well

4. Confirm whether the target protein is synthesized through staining with Coomassie blue and de-staining procedures (Figure 7, 8).


Figure 7. Confirmation of SDS-PAGE results for protein synthesis using ExiProgen™ EC-Tagfree Protein Synthesis Kit. M, Protein Size Marker; Lane 1, Dual specificity protein phosphatase 3 (22 kDa); Lane 2, Human growth hormone (23 kDa); Lane 3, CAT (24 kDa); Lane 4, AcGFP (28 kDa).


Figure 8. SDS-PAGE data to confirm removal of 6x histidine tag and Western Blot result using His-tag antibody from proteins synthesized through *ExiProgen*[™] EC-Tagfree Protein Synthesis Kit. M, Protein Size Marker; Lane 1, Expression sample of pBIVT-TEV-AcGFP; Lane 2, Histidine tag-TEV-AcGFP (histidine tagged sample); Lane 3, AcGFP (histidine tag-free sample).

Maintenance

After protein synthesis, Cartridges and other accessories should be washed and stored as follows.

1. Reaction block

Rinse all the used wells of the Reaction block with sterile distilled water and store it in 70% ethanol at a room temperature. Before use, wash with sterile distilled water and dry.

2. Waste Tray and Disposable Tip Rack

Discard all the solution in the Waste Tray, wash it in running water and clean it with 20% ethanol. Also, if there are any impurities on the Disposable Tip Rack, clean it with 20% ethanol.

3. Cartridge ①, ②

After the reaction is completed, keep the remaining Cartridges with the lid covered, and store Cartridge ① at 4°C and Cartridge ② at -20°C at the appropriate temperature, respectively.

Troubleshooting

If there is a problem with protein synthesis, please refer to the following information (However, keep in mind that the following suggests solutions for general problems regarding protein synthesis, and may not correspond to solutions for all protein synthesis problems).

1. No protein synthesis, including positive control protein.

Cause	Solution
Nuclease (DNase, RNase) contamination	When performing the experiment, wear gloves and use DNase- and RNase-free pipette tips.
No addition of DNA or missed components	Make sure to add the correct amount of Positive control DNA and ensure that all components are installed correctly.
Storage condition of reagent	Store all reagents and components in the kit at the recommended temperatures. In particular, avoid repeated freeze-thaw cycles of <i>E. coli</i> extract.

2. Positive control protein is synthesized, but the target protein is not synthesized.

Cause	Solution
Sequences of template DNA	If there is a mutation in the ORF (Open Reading Frame) of the template DNA, translation may be interrupted, so please check the ORF sequence.
Structures of template DNA	Confirm that the structures of the T7 promoter, T7 terminator, and histidine tag are correctly positioned in the template DNA.
Contamination of template DNA	If DNA is degraded by nuclease contamination in the process of preparing template DNA, protein cannot be synthesized. It is recommended to use template DNA obtained using nuclease-free elution buffer.
Low expression efficiency	<i>ExiProgen</i> [™] Protein Expression Optimization Kit can be applied to increase protein expression efficiency.

3. Low synthesis amount of target protein.

Cause	Solution
Purity of template DNA	If the purity of template DNA is low, the protein may not be synthesized. It is recommended to use DNA as follows: $A_{260/280}$: 1.7-2.0, $A_{260/230} > 1.5$.
Amount of template DNA	There may be differences in the amount of protein synthesis depending on the amount of template DNA added. To increase the amount of protein synthesis, optimize the DNA concentration through expression level screening test prior to the experiment (See Appendix B).
Unoptimized template DNA for <i>E. coli</i> codons	If the sequences of the template DNA are not optimized for <i>E. coli</i> codons, protein synthesis may not occur or the amount of synthesis may be low. It is recommended to use codon-optimized DNA.
Position of histidine tag	Depending on the location of the histidine tag, the target protein expression may be inhibited. Or, it might not affect protein expression, but it may affect purification when the tag is not exposed to the outside during the formation of the tertiary structure. In this case, change the position of the tag.


4. The activity or solubility of the target protein is low.

Cause	Solution
Proteins that require post-translational modification	In the cell-free protein synthesis system using <i>E. coli</i> extract, it is impossible to synthesize proteins that require post-translational modifications such as glycosylation and phosphorylation.
When specific elements are required for protein activity	If the synthesized protein requires specific elements to exhibit activity, the activity should be confirmed after adding the elements to the final purified protein solution.
Aggregation due to low solubility of proteins	If the protein synthesis temperature is lowered or the chaperone protein is added to help fold the protein, the solubility may be improved (Frydman, J., 2001).

Appendix A: pBIVT Vectors

BIONEER's cell-free protein synthesis vector, pBIVT Vector Set-1, is compatible with all kits of *ExiProgen™* protein synthesis system. This vector set contains two types of plasmids, so 6x histidine tagging is possible at either the N-terminus or C-terminus. The structure of the vector is as shown below.

- pBIVT Vector Map

- pBIVT-1's MCS (Multiple Cloning Site) sequences

ATG	CATATG	CACCACCACCACCACCAC	GGATCC	GAGCTC	AAGCTT	GCGGCCGC	A	TAG	GTCGAC
Start	NdeI	6X His-tag	BamHI	SacI		NotI		Stop	SalI
codon								codon	

- pBIVT-2's MCS (Multiple Cloning Site) sequences

ATG	CATATG	GGATCC	GAGCTC	AAGCTT	GCGGCCGC	A	CACCACCACCACCACCAC	TAG	GTCGAC
Start	NdeI	BamHI	SacI		NotI		6X His-tag	Stop	SalI
codon								codon	

Appendix B: Screening of Template DNA Concentration

When synthesizing proteins using the *ExiProgen*™ protein synthesis system, there may be differences in the amount of protein synthesis depending on the amount of template DNA. Optimal protein synthesis results can be obtained by first performing expression screening according to the concentration of template DNA as shown below and then synthesizing the protein using *ExiProgen*™.

1. AccuRapid™ Cell-Free Protein Expression Kit (Cat. No. K-7250)

1) Prepare by measuring the concentration of template DNA.

Example) 100 ng/μl

2) Prepare the reaction mixture as follows.

DNA concentration screening						
Components	Volume	120 ng	240 ng	360 ng	480 ng	600 ng
<i>E. coli</i> extract	12 μl	12 μl	12 μl	12 μl	12 μl	12 μl
Master Mix	21 μl	21 μl	21 μl	21 μl	21 μl	21 μl
Template DNA	X μl	1.2 μl	2.4 μl	3.6 μl	4.8 μl	6 μl
DEPC DW	12-X μl	10.8 μl	9.6 μl	8.4 μl	7.2 μl	6 μl
Total volume	45 μl	45 μl	45 μl	45 μl	45 μl	45 μl

3) Incubate the reaction mixtures at 30°C for 3 hrs.

4) Confirm the optimal expression concentration through SDS-PAGE.

2. ExiProgen™ EC-Tagfree Protein Synthesis Kit (Cat. No. K-7320)

1) Confirm the DNA concentration equivalent to **10 times** the optimal DNA concentration identified in the previous experiment, add protein expression solution each Dialysis tube in the **row A** of the Reaction block.

Example) EC-Tagfree Kit expression volume (450 μl) ÷ screening expression volume (45 μl) ≈ 10

Optimal DNA concentration through screening (600 ng) x 10 ≈ 6 μg

Optimal DNA concentration through screening (120 ng) x 10 ≈ 1.2 μg

Reference

- Ahn, J. H., Kang, T. J., & Kim, D. M. (2008). Tuning the expression level of recombinant proteins by modulating mRNA stability in a cell-free protein synthesis system. *Biotechnology and bioengineering*, 101(2), 422-427.
- Forster, A. C., Cornish, V. W., & Blacklow, S. C. (2004). Pure translation display. *Analytical biochemistry*, 333(2), 358-364.
- Frydman, J. (2001). Folding of newly translated proteins in vivo: the role of molecular chaperones. *Annual review of biochemistry*, 70(1), 603-647.
- Gualerzi, C. O., Giuliodori, A. M., & Pon, C. L. (2003). Transcriptional and post-transcriptional control of cold-shock genes. *Journal of molecular biology*, 331(3), 527-539.
- Hino, M., Kataoka, M., Kajimoto, K., Yamamoto, T., Kido, J. I., Shinohara, Y., & Baba, Y. (2008). Efficiency of cell-free protein synthesis based on a crude cell extract from Escherichia coli, wheat germ, and rabbit reticulocytes. *Journal of biotechnology*, 133(2), 183-189.
- Josephson, K., Hartman, M. C., & Szostak, J. W. (2005). Ribosomal synthesis of unnatural peptides. *Journal of the American Chemical Society*, 127(33), 11727-11735.
- Keum, J. W., Ahn, J. H., Choi, C. Y., Lee, K. H., Kwon, Y. C., & Kim, D. M. (2006). The presence of a common downstream box enables the simultaneous expression of multiple proteins in an E. coli extract. *Biochemical and biophysical research communications*, 350(3), 562-567.
- Keum, J. W., Ahn, J. H., Kang, T. J., & Kim, D. M. (2009). Combinatorial, selective and reversible control of gene expression using oligodeoxynucleotides in a cell-free protein synthesis system. *Biotechnology and bioengineering*, 102(2), 577-582.
- Kigawa, T., Yamaguchi-Nunokawa, E., Kodama, K., Matsuda, T., Yabuki, T., Matsuda, N., ... & Yokoyama, S. (2002). Selenomethionine incorporation into a protein by cell-free synthesis. *Journal of structural and functional genomics*, 2(1), 29-35.
- Kim, D. M., & Choi, C. Y. (1996). A semicontinuous prokaryotic coupled transcription/translation system using a dialysis membrane. *Biotechnology progress*, 12(5), 645-649.

Kim, T. W., Oh, I. S., Keum, J. W., Kwon, Y. C., Byun, J. Y., Lee, K. H., ... & Kim, D. M. (2007). Prolonged cell-free protein synthesis using dual energy sources: Combined use of creatine phosphate and glucose for the efficient supply of ATP and retarded accumulation of phosphate. *Biotechnology and bioengineering*, 97(6), 1510-1515.

Kim, H. C., Kim, T. W., & Kim, D. M. (2011). Prolonged production of proteins in a cell-free protein synthesis system using polymeric carbohydrates as an energy source. *Process Biochemistry*, 46(6), 1366-1369.

Li, Q., Yi, L., Marek, P., & Iverson, B. L. (2013). Commercial proteases: present and future. *FEBS letters*, 587(8), 1155-1163.

Lim, H. J., & Kim, D. M. (2019). Cell-free metabolic engineering: recent developments and future prospects. *Methods and Protocols*, 2(2), 33.

Ohashi, H., Shimizu, Y., Ying, B. W., & Ueda, T. (2007). Efficient protein selection based on ribosome display system with purified components. *Biochemical and biophysical research communications*, 352(1), 270-276.

Park, S., & Hamad-Schifferli, K. (2010). Enhancement of in vitro translation by gold nanoparticle– DNA conjugates. *ACS nano*, 4(5), 2555-2560.

Rungpragayphan, S., Nakano, H., & Yamane, T. (2003). PCR-linked in vitro expression: a novel system for high-throughput construction and screening of protein libraries. *FEBS letters*, 540(1-3), 147-150.

Trotta, E. (2011). The 3-base periodicity and codon usage of coding sequences are correlated with gene expression at the level of transcription elongation. *PloS one*, 6(6), e21590.

Villemagne, D., Jackson, R., & Douthwaite, J. A. (2006). Highly efficient ribosome display selection by use of purified components for in vitro translation. *Journal of immunological methods*, 313(1-2), 140-148.

Ordering Information

Description		Cat. No
<i>ExiProgen</i> TM EC-Tagfree Protein Synthesis Kit	8 rxn	K-7320

Related Products

Description		Cat. No
<i>AccuRapid</i> TM Cell-Free Protein Expression Kit		K-7250
<i>AccuRapid</i> TM Midi Protein Expression Kit		K-7260
<i>AccuRapid</i> TM Maxi Protein Expression Kit		K-7270
<i>AccuRapid</i> TM Protein Synthesis Kit		K-7280
<i>ExiProgen</i> TM ProXpress PCR Template Kit		K-7400
First primer F/R sets (N terminus 6x His tag) (each 5 nmole)		N-8229
First primer F/R sets (C terminus 6x His tag) (each 5 nmole)		N-8230
<i>ExiProgen</i> TM Protein Expression Optimization Kit		K-7410
pBIVT Vector Set-1		K-7350
<i>ExiProgen</i> TM EC-Maxi Protein Synthesis Kit		K-7310
<i>ExiProgen</i> TM EC Protein Synthesis Kit		K-7300
<i>ExiProgen</i> TM EC-Disulfide Protein Synthesis Kit		K-7330
<i>ExiProgen</i> TM EC-Bulk Protein Synthesis Kit		K-7340
<i>ExiProgen</i> TM His-tagged Protein Purification Kit		K-7220
<i>ExiProgen</i> TM Dialysis Kit		K-7240
<i>ExiProgen</i> TM Consumable SET		KA-3001
Gene Synthesis Service		S-2041
Protein Synthesis Service		S-2500
<i>ExiProgen</i> TM		A-5041

Explanation of Symbols

 Batch Code	 Biological Risks	 Catalog Number	 Caution
 Consult Instructions For Use	 Contains Sufficient for <n> tests	 Do not Re-use	 Manufacturer
 Research Use Only	 Temperature Limitation	 Use-by date	

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1. *AccuRapid*[™] Cell-Free Protein Expression Kit (Cat. No. K-7250)

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2. *ExiProgen*[™] EC-Tagfree Protein Synthesis Kit (Cat. No. K-7320)

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제품명		Cat. No
AccuRapid™ Cell-Free Protein Expression Kit		K-7250
AccuRapid™ Midi Protein Expression Kit		K-7260
AccuRapid™ Maxi Protein Expression Kit		K-7270
AccuRapid™ Protein Synthesis Kit		K-7280
ExiProgen™ ProXpress PCR Template Kit		K-7400
First primer F/R sets (N terminus 6x His tag) (each 5 nmole)		N-8229
First primer F/R sets (C terminus 6x His tag) (each 5 nmole)		N-8230
ExiProgen™ Protein Expression Optimization Kit		K-7410
pBIVT Vector Set-1		K-7350
ExiProgen™ EC-Maxi Protein Synthesis Kit		K-7310
ExiProgen™ EC Protein Synthesis Kit		K-7300
ExiProgen™ EC-Disulfide Protein Synthesis Kit		K-7330
ExiProgen™ EC-Bulk Protein Synthesis Kit		K-7340
ExiProgen™ His-tagged Protein Purification Kit		K-7220
ExiProgen™ Dialysis Kit		K-7240
ExiProgen™ Consumable SET		KA-3001
Gene Synthesis Service		S-2041
Protein Synthesis Service		S-2500
ExiProgen™		A-5041

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BIONEER Global Center

Address 71, Techno 2-ro, Yuseong-gu, Daejeon, 34013, Republic of Korea
E-mail sales@bioneer.co.kr
Web www.bioneer.com

BIONEER R&D Center

Address Korea Bio Park BLDG #B-702, 700 Daewangpangyo-ro, Bundang-gu, Seongnam-si
Gyeonggi-do, 13488, Republic of Korea
E-mail sales@bioneer.co.kr
Web www.bioneer.com

BIONEER Inc. - USA Branch

Address 155 Filbert St. Suite 216 Oakland, CA 94607, USA
E-mail order.usa@bioneer.com
Web us.bioneer.com

BIONEER Corp. - European Branch

Address Ludwig-Erhard-Strasse 30-34, 65760 Eschborn, Germany
E-mail euinfo@bioneer.com
Web www.bioneer.com

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BIONEER Corporation - HQ