ExiProgen[™] **His-tagged Protein Purification Kit**

Cat. No. K-7220
K-7221

***ExiProgen*[™] His-tagged Protein Purification Kit**

User Guide

K-7220

K-7221

Version No.: 2 (2022-06-23)

Please read all the information in booklet before using the unit

BIONEER Corporation

**Bioneer Global Center, 71, Techno-2-ro,
Yuseong-gu, Daejeon, 34013, Republic of Korea**

Tel: 1588-9788

Email: sales@bioneer.co.kr

www.bioneer.com

Intended Use

ExiProgen[™] His-tagged Protein Purification Kit is developed and supplied for research purposes only. Certain applications possible with this kit may require special approval by appropriate local and/or national regulatory authorities in the country of use.

Safety Warning and Precaution

Wear appropriate protection when handling any irritant or harmful reagents. The use of a laboratory coat, protective gloves and safety goggles are highly recommended. For more information, please consult the appropriate Material Safety Data Sheet (MSDS).

Warranty and Liability

All BIONEER products undergo extensive Quality Control testing and validation. BIONEER guarantees quality during the warranty period as specified, when following the appropriate protocol as supplied with the product. It is the responsibility of the purchaser to determine the suitability of the product for its particular use. Liability is conditional upon the customer providing full details of the problem to Bioneer within 30 days.

Quality Management System ISO 9001 Certified

Every aspect of our quality management system from product development, production to quality assurance and supplier qualification meets the world-class standards.

Patent

ExiProgen[™] and its kits are protected by the patents KR10-2011-0085824, and PCT/KR2012/006715.

Trademark

ExiProgen[™] is a trademark of BIONEER Corporation.

Copyright

Copyright 2022. BIONEER Corporation. All Rights Reserved.

Notice

BIONEER corporation reserves the right to make corrections, modifications, improvements and other changes to its products, services, specifications or product descriptions at any time without notice.

Contents

Product Information	1
Components	1
Storage	1
Specifications.....	1
Precautions	1
Introduction	2
Cell-free Protein Synthesis System	2
The <i>ExiProgen</i> ™ Protein Synthesis System	4
Product Description	6
Experimental Procedures	7
Components Details	7
Experiments Preparation	9
Protocols	10
Analysis of Sample	16
Maintenance	19
Troubleshooting	20
References	22
Ordering Information	24
Related Products	24
Explanation of Symbols	25

Product Information

Components

<i>ExiProgen</i>[™] His-tagged Protein Purification Kit			
Components	K-7220 (16 rxn)	K-7221 (32 rxn)	Storage
Cartridge ①	96-well x 1 ea	96-well x 2 ea	4°C
Cartridge ②	96-well x 1 ea	96-well x 2 ea	
Disposable filter tip	2 pack (8 ea/pack)	1 pack (33 ea/pack)	-
Elution tubes and caps	8-tube strips x 2 ea	8-tube strips x 4 ea	-

Storage

The Cartridge ① and ② included in *ExiProgen*[™] His-tagged Protein Purification Kit contain magnetic Ni-NTA magnetic nanoparticles and buffers used for protein purification, so they should be stored at 4°C.

Specifications

<i>ExiProgen</i>[™] His-tagged Protein Purification Kit		
	K-7220 (16 rxn)	K-7221 (32 rxn)
Number of reactions (Up to 16 rxn/ 1 times)	16 rxn	32 rxn
Operating time	2 hours	
Sample	Cell lysate, Cell-free protein expression sample (<i>in vitro</i>)	
Binding capacity	3 mg of protein/ml of beads	

Precautions

The Cartridge ① and ② of this Synthesis Kit are each covered with sealing film in order to prevent cross-contamination, evaporation, or leakage of solutions inside. All of the plastic products and buffers in this kit are provided under nuclease- and protease-free condition, hence, please be careful not to contaminate any part of the kit with nuclease or protease.

Introduction

Cell-free Protein Synthesis System

Since protein composes essential components in biological reactions, such as enzymes, hormones, and structural proteins, researches on roles and structures of protein are being actively conducted in the post-genomic era. These studies start with producing a specific type of protein.

The various protein expression systems are based on different cell, such as *E. coli*, yeast, animal, and plants cells, which are used for production of recombinant proteins. In order to express proteins based on these systems, a vector containing target genes is transformed into cells and the cell line is cultured to express a large amount of protein. After that, cells are usually disrupted and proteins are purified from cell lysates. However, the general protein expression systems require much time and labor, as these have to go through a series of steps such as selection of: 1) strains, 2) cell lines stably expressing recombinant proteins, 3) cell culture, 4) cell disruption, and 5) protein purification. In addition, if the target protein is toxic to strains and cell lines, it takes several weeks or up to several months to synthesize a single protein because it is difficult to express a toxic protein and requires optimization of various conditions for its expression.

To overcome the limitations of cell-based system, cell-free protein synthesis method and related products have been developed (Keum et al., 2009; Kim et al., 2007; Kim et al., 2011; Park & Hamad-Schifferli, 2010). Cell-free protein synthesis system is a technology for protein synthesis based on cell extract by transcribing and translating *in vitro*. To carry out protein synthesis *in vitro*, the following should be prepared: 1) cell extract including T7 RNA polymerase, ribosome, tRNA, and enzyme, and 2) protein expression solution containing a substrate and a substance that performs transcription and translation such as amino acids, rNTPs, and energy sources. After that, a gene encoding the target protein is added, and the reaction is conducted at an appropriate temperature and time to express the recombinant protein (Figure 1). In case of cell-free protein synthesis systems, the reaction duration is at least 1 to 3 hours, which can dramatically reduce the time taken for protein expression. In addition, it is easier to express proteins that show toxicity in cells because the expression is based on cell extracts, not living cells, and since it is an open system, conditions for protein expression can be easily adjusted (Hino et al., 2008).

Cell-free protein synthesis systems are widely used in enzyme engineering, protein labeling, research on protein-protein interaction mechanism, and research on protein active site (Forster et

al., 2004; Josephson et al., 2005; Keum et al., 2006; Kigawa et al., 2002; Ohashi et al., 2007; Villemagne et al., 2006).


Figure 1. Principles of cell-free protein synthesis system.

The ExiProgen™ Protein Synthesis System

BIONEER has developed *ExiProgen*™, a fully automatic protein synthesis instrument that applies cell-free protein synthesis technology and affinity purification technology using magnetic nanoparticles. Therefore, high-purity proteins can be obtained quickly and fully automatically by using our protein expression and purification kits. Furthermore, the instrument *ExiProgen*™ is highly versatile as it can fully automatically extract the required DNA and RNA from various samples using nucleic acid extraction kit series.

Our protein expression and purification kits use T7 expression system of *E. coli* and can be used immediately with template DNA (Ahn et al., 2008; Rungpragayphan et al., 2003). Protein synthesis using *ExiProgen*™ allows the addition of template DNA to perform a series of processes, including protein expression and purification, to finally obtain samples containing target proteins.

Protein expression and purification kits are divided into three product groups (Figure 2): 1) Template DNA preparation kits that produce template DNA used for cell-free protein expression, 2) *AccuRapid*™ Kit Series that is responsible for easy manual expression or synthesis of recombinant proteins, and 3) *ExiProgen*™ Kit Series that is responsible for an automatic expression or synthesis of recombinant proteins. Each protein expression and purification kit can express and synthesize proteins from µg up to mg scale.

In addition, there are other products that can be selected depending on the purpose of the experiment, such as a kit for improving protein expression efficiency, a kit solely for protein purification or buffer exchange, etc.

Please refer to the BIONEER's website (www.bioneer.com) for a detailed description of each kit.


Figure 2. Product information related to protein expression and synthesis.

Product Description

ExiProgen™ His-tagged Protein Purification Kit can be applied to *ExiProgen™*, a fully automated protein synthesis and nucleic acid extraction system developed by BIONEER, and target proteins can be obtained from protein expression samples with high purity.

This product purifies proteins by affinity method using Ni-NTA beads and Histidine-tag from cell lysate or cell-free protein expression samples. The target protein is captured to the surface of the Ni-NTA magnetic bead equilibrated with Equilibrium buffer and washed with washing buffer. After that, the target protein is eluted from the beads by treating the elution buffer including high concentration of imidazole.

A series of protein purification processes can save a lot of time and labor since users do not have to perform the entire process except for sample addition by applying an automated system (Figure 3). In addition, it can have high reproducibility, unlike manual type experiments that can obtain different results for each user.



Figure 3. Principles of *ExiProgen™* His-tagged Protein Purification Kit.

Experimental Procedures


Components Details

1. Cartridges

2. Note

Structure of ExiProgen™ Baseplate

Experiments Preparation

Preparation of Sample

Prepare the cell lysate or cell-free protein expression sample to be applied to *ExiProgen™* His-tagged Protein Purification Kit.





1. Cell lysate (e.g., *E. coli* cell lysate)


- 1) Harvest cells of target protein that is overexpressed by centrifuging at 3,000 rpm for 10 min.
- 2) Add an appropriate amount of lysis buffer* (not provided) to the cell pellet and resuspend the cells completely.
* **Example of lysis buffer:** 20 mM Tris-HCl, 1 mM 2-Mercaptoethanol, pH 7.6
- 3) Lyse the cells on ice with a sonicator.
(Microtip, On/10 sec, Off/20 sec, Total time/5 min)
- 4) Transfer the cell lysate to the 1.5 ml tube and centrifuge at 13,000 rpm for 1 min to separate supernatant from pellet.
- 5) Confirm whether the target protein in the supernatant is expressed or not with SDS-PAGE.

2. Cell-free protein expression sample (*In vitro*)

- 1) Prepare the protein expression solution with a final volume of 750 µl in a 1.5 ml tube using the *AccuRapid™* Midi Protein Expression Kit (Cat. No. K-7260).
- 2) Place the prepared protein expression solution in an incubator and incubate at 30°C for 3 hrs.
- 3) Confirm whether the target protein in the supernatant is expressed or not with SDS-PAGE.


Protocols

Steps	Procedure Details
<p>Example 1) For 1 sample</p>  <p>Example 2) For 8 samples</p>  <p>Example 3) For 16 samples</p> 	<p>1. Punch holes in the sealing films of Cartridge ① and ② using 6 Hole Punch (<i>ExiProgen™</i> accessory) according to the number of samples.</p> <p>* Note: Refer to the left pictures and punch holes depending on the required number of samples.</p>
	<p>2. Add 700-750 µl of cell lysate or protein expression sample prepared on page 9 to each column in rows G of Cartridge ②.</p>
	<p>3. Open the door of <i>ExiProgen™</i> and pull out the baseplate completely.</p>


4. Load Cartridge ② in the position of ② on the baseplate.

*** Note:** Place the row L of Cartridge ② in the direction of the Heating block first and then install it, and afterwards check if it is properly fixed.

5. Load Cartridge ① in the position of ① on the baseplate.

*** Note:** There are silicon rings embedded on both sides of the Cartridge ① installation position. Place the left side first, and then the right side.

6. Install Waste Tray between Cartridge ① and ②.


*** Note:** Follow the exact installation order of Cartridge ② → Cartridge ① → Waste Tray. Ensure that Cartridges and Waste Tray are firmly installed, without showing any movement.


7. Load Elution Tube Rack (*ExiProgen™* accessory) with Elution tube on the Heating and cooling block (Magnetic part) of the baseplate.

*** Note:** Elution tube should be installed on the Elution Tube Rack in an open state.


*** Note:** Please refer to the illustration below when placing Elution tube to Elution Tube Rack and install it in the correct position in the rows in alphabet. Elution tube can be used according to the number of samples and they must be placed on the same row as punched cartridges.

8. Place the Protection Cover as shown in the picture on the left.

*** Note:** Press down Protection Cover firmly so that it covers Elution Tube Rack all the way down to the baseplate.

9. Place Disposable filter tips in Disposable Tip Rack according to the number of samples.

*** Note:** Tips should be placed in the corresponding positions with the punched holes of the Cartridges. Do not insert tips in the corresponding positions of the cartridge that are not punched.


10. Finally, confirm holes in the cartridges and position of samples and tips (Refer to the pictures on page 13).


Push the baseplate completely until you hear the click sound, then close the door.

*** Note**


Example 1) For 1 Sample

Example 2) For 8 Samples

Example 3) For 16 Samples

11. Turn on the *ExiProgen™* and tap 'Press to start' button. When the 'Press to start' button is pressed, the *ExiProgen™* screen appears as shown in the lower left, and after the scroll bar moves, it moves to the next screen.

*** Note:** This process initializes the X, Y, Z axis values of the equipment. If the next screen does not proceed normally, turn off the power of the equipment and contact the A/S center.

12. In the MENU screen, press 'Start' button to select a proper protocol.


13. The PREP SETUP screen appears as shown in the left, and a screen to select the protocol number for each kit appears. Press '901' to confirm that the following information is displayed on the screen, and then press the 'Enter' button.

Prep Type: Protein

Sample SRC: Protein_Purification

14. After selecting the protocol, the screen to select the Elution volume appears. It is used for nucleic acid extraction, and when using this kit, please press the 'ok' button to move next steps.

15. After the Elution volume selection screen, the screen to select protein purification temperature appears. When using this kit, please select 30°C and press the 'ok' button.

16. When the 'CHECK LIST' screen appears, confirm once again that cartridges and other accessories are fixed to each location properly, and then click the 'ok' button.

17. When the 'Running Mode' screen appears as shown on the left, finally check the following information and click the 'RUN' button. After this, the instrument starts and it takes approximately **2 hours** to finish run.


Prep Type: Protein

Sample SRC: Protein_Purification

18. 'Work Completion' screen appears when the protocol is completed. Remove all components used in the experiment, and choose a desired button.

*** Note:** If you want to quit and press the 'ok' button, the UV lamp will be turned on automatically.

19. After taking out the sample rack, Protection Cover can be removed by using Protection Cover Separation Tool (*ExiProgen*TM's accessory).


Take out Elution Tube Rack from *ExiProgen*TM, place it on top of the Protection Cover Separation Tool.

*** Note:** When placing Elution Tube Rack on Protection Cover Separation Tool, the lever must be facing left-hand side.


20. Fix the Protection Cover and Separation Tool so that it does not move with one hand, and turn the lever clockwise by 180° with the other hand.

*** Note:** Rotate the lever until Elution Tube Rack is firmly fixed to Protection Cover Separation Tool.

21. After fixing the Elution Tube Rack, press down both sides of Separation Tool. When the lever is pulled, the Protection Cover rises upwards and gets separated from the Elution Tube Rack for easy detachment.

*** Note:** When rotating lever, while holding the Protection Cover with one hand, pull the lever sequentially to prevent the solution from splashing in the tube.

Analysis of Sample

After protein purification with *ExiProgen*™, the final target proteins of about 250 µl is collected from Elution tube. Final protein samples may contain Ni-NTA magnetic beads, but this can be removed by centrifugation before use.

If necessary, reaction samples from protein purification can be recovered from the rows L, K, J of Cartridges ② as shown in Figure 4.


Figure 4. Samples for each row of Cartridge ②.

- Unbound sample: Samples not bound to Ni-NTA magnetic bead.
- 1st washing sample: Samples after 1st washing step in purification process.
- Bead sample: Used bead samples for purification.

Each sample can be checked through SDS-PAGE whether target protein is purified as desired.

1. Prepare the loading mixture as shown in the table.

Components	Unbound/1 st washing sample	Purified protein/Bead sample*
Sample	5 µl	15 µl
4X loading dye	5 µl	5 µl
Sterile distilled water	10 µl	-
Total volume	20 µl	20 µl

* **Note:** Add 200 µl of sterile distilled water to the well containing the beads for suspension and proceed with the sampling afterwards.

2. Incubate the prepared loading mixture at 95°C for 5-10 min.

3. Load the following amounts of each sample to the wells of 10% or 12% SDS-PAGE gel [10 x 8 (cm), 10-well] and run the SDS-PAGE.

- Unbound/ 1st washing samples: 5 µl/well
- Purified protein/Bead sample: 10 µl/well

4. Confirm whether the target protein is synthesized and purified through staining with Coomassie blue and de-staining procedures (Figure 5, 6).


Figure 5. Purification of target protein using *E. coli* cell lysate. (A) *E. coli* cell lysate overexpressed DUSP3. M, Protein Size Marker; BI, Pre-induction sample; I, Post-induction sample. **(B)** Purification of target protein using *ExiProgen*TM. M, Protein Size Marker; CL, Recombinant cell lysate; P, Purification sample.


Figure 6. Purification of the target protein using cell-free protein expression sample. (A) Protein expression sample. **(B)** Purification of target protein using *ExiProgen™*. M, Protein Size Marker; 1, DUSP3 (22 kDa); 2, CAT (24 kDa); 3, AcGFP (28 kDa); 4, ERFP (31 kDa); 5, EF-Ts (34 kDa); 6, VF (45 kDa).

Maintenance

After protein purification, Cartridges and other accessories should be washed and stored as follows.

1. Waste Tray and Disposable Tip Rack

Discard all the solution in the Waste Tray, wash it in running water and clean it with 20% ethanol. Also, if there are any impurities on the Disposable Tip Rack, clean it with 20% ethanol.

2. Cartridge ①, ②

After the reaction is completed, keep the remaining Cartridges with the lid covered, and store Cartridge ① and ② at 4°C.

Troubleshooting

If there is a problem with protein purification, please refer to the following information (However, keep in mind that the following suggests solutions for general problems regarding protein purification only, and may not correspond to solutions for all protein purification problems).

1. No protein purification or low amount of purified target protein.

Cause	Solution
Low amount of expression of target protein	<p>1. When using cell lysate sample. Optimize the protein expression level by changing the incubation temperature, inoculation concentration, time, etc.</p> <p>2. When using cell-free protein expression sample. 2-1) Optimize codons according to the cell line of cell extract with sequences of template DNA used for protein expression. 2-2) Use the manual <i>AccuRapid™</i> kit to increase the yield by screening for the optimal expression conditions, through template concentration optimization process on a small scale.</p>
Formation of inclusion body	If the expressed target protein forms an aggregate, it may not bind to the Ni-NTA magnetic bead. Optimize the expression process to minimize protein aggregation.
Low amount of loading sample	Reduce the amount of lysis buffer when preparing cell lysate. Then, high concentration samples can be prepared for protein purification.

No histidine tag	Check if the Histidine tag is attached to the N-terminus or C-terminus, and if not, prepare the tagged sample again.
Position of histidine tag	Depending on the location of the histidine tag, the target protein expression may be inhibited. Or, it might not affect protein expression, but it may affect purification when the tag is not exposed to the outside during the formation of the tertiary structure. In this case, change the position of the tag.
Degradation of target DNA	Add protease inhibitor during cell lysis.

2. When the target protein is identified in Unbound or 1st washing sample.

Cause	Solution
In case of overloading	This corresponds to the case where an amount of sample exceeding the binding capacity of the Ni-NTA magnetic bead is loaded. In this case, dilute the sample before use.
Protein precipitation	Solubilization reagent (e.g., 0.1% Triton X-100 or Tween 20) can be added to the buffer.

Reference

- Ahn, J. H., Kang, T. J., & Kim, D. M. (2008). Tuning the expression level of recombinant proteins by modulating mRNA stability in a cell-free protein synthesis system. *Biotechnology and bioengineering*, 101(2), 422-427.
- Forster, A. C., Cornish, V. W., & Blacklow, S. C. (2004). Pure translation display. *Analytical biochemistry*, 333(2), 358-364.
- Hino, M., Kataoka, M., Kajimoto, K., Yamamoto, T., Kido, J. I., Shinohara, Y., & Baba, Y. (2008). Efficiency of cell-free protein synthesis based on a crude cell extract from Escherichia coli, wheat germ, and rabbit reticulocytes. *Journal of biotechnology*, 133(2), 183-189.
- Josephson, K., Hartman, M. C., & Szostak, J. W. (2005). Ribosomal synthesis of unnatural peptides. *Journal of the American Chemical Society*, 127(33), 11727-11735.
- Keum, J. W., Ahn, J. H., Choi, C. Y., Lee, K. H., Kwon, Y. C., & Kim, D. M. (2006). The presence of a common downstream box enables the simultaneous expression of multiple proteins in an E. coli extract. *Biochemical and biophysical research communications*, 350(3), 562-567.
- Keum, J. W., Ahn, J. H., Kang, T. J., & Kim, D. M. (2009). Combinatorial, selective and reversible control of gene expression using oligodeoxynucleotides in a cell-free protein synthesis system. *Biotechnology and bioengineering*, 102(2), 577-582.
- Kim, T. W., Oh, I. S., Keum, J. W., Kwon, Y. C., Byun, J. Y., Lee, K. H., Choi, C. Y., & Kim, D. M. (2007). Prolonged cell-free protein synthesis using dual energy sources: Combined use of creatine phosphate and glucose for the efficient supply of ATP and retarded accumulation of phosphate. *Biotechnology and bioengineering*, 97(6), 1510-1515.
- Kim, H. C., Kim, T. W., & Kim, D. M. (2011). Prolonged production of proteins in a cell-free protein synthesis system using polymeric carbohydrates as an energy source. *Process Biochemistry*, 46(6), 1366-1369.
- Kigawa, T., Yamaguchi-Nunokawa, E., Kodama, K., Matsuda, T., Yabuki, T., Matsuda, N., Yabuki, T., Matsuda, N., Ishitani, R., Nureki, O., & Yokoyama, S. (2002). Selenomethionine incorporation into a protein by cell-free synthesis. *Journal of structural and functional genomics*, 2(1), 29-35.

Ohashi, H., Shimizu, Y., Ying, B. W., & Ueda, T. (2007). Efficient protein selection based on ribosome display system with purified components. *Biochemical and biophysical research communications*, 352(1), 270-276.

Park, S., & Hamad-Schifferli, K. (2010). Enhancement of in vitro translation by gold nanoparticle– DNA conjugates. *ACS nano*, 4(5), 2555-2560.

Rungpragayphan, S., Nakano, H., & Yamane, T. (2003). PCR-linked in vitro expression: a novel system for high-throughput construction and screening of protein libraries. *FEBS letters*, 540(1-3), 147-150.

Villemagne, D., Jackson, R., & Douthwaite, J. A. (2006). Highly efficient ribosome display selection by use of purified components for in vitro translation. *Journal of immunological methods*, 313(1-2), 140-148.

Ordering Information

Description		Cat. No
<i>ExiProgen™</i> His-tagged Protein Purification Kit	16 rxn	K-7220
	32 rxn	K-7221

Related Products

Description		Cat. No
<i>AccuRapid™</i> Cell-Free Protein Expression Kit		K-7250
<i>AccuRapid™</i> Midi Protein Expression Kit		K-7260
<i>AccuRapid™</i> Maxi Protein Expression Kit		K-7270
<i>AccuRapid™</i> Protein Synthesis Kit		K-7280
<i>ExiProgen™</i> ProXpress PCR Template Kit		K-7400
First primer F/R sets (N terminus 6x His tag) (each 5 nmole)		N-8229
First primer F/R sets (C terminus 6x His tag) (each 5 nmole)		N-8230
<i>ExiProgen™</i> Protein Expression Optimization Kit		K-7410
pBIVT Vector Set-1		K-7350
<i>ExiProgen™</i> EC Protein Synthesis Kit		K-7300
<i>ExiProgen™</i> EC-Maxi Protein Synthesis Kit		K-7310
<i>ExiProgen™</i> EC-Tagfree Protein Synthesis Kit		K-7320
<i>ExiProgen™</i> EC-Disulfide Protein Synthesis Kit		K-7330
<i>ExiProgen™</i> EC-Bulk Protein Synthesis Kit		K-7340
<i>ExiProgen™</i> Dialysis Kit		K-7240
<i>ExiProgen™</i> Consumable SET		KA-3001
Gene Synthesis Service		S-2041
Protein Synthesis Service		S-2500
<i>ExiProgen™</i>		A-5041

Explanation of Symbols

 LOT	Batch Code		Biological Risks	 REF	Catalog Number		Caution
	Consult Instructions For Use		Contains Sufficient for <n> tests		Do not Re-use		Manufacturer
 RUO	Research Use Only		Temperature Limitation		Use-by date		

ExiProgen™ His-tagged Protein Purification Kit

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Ahn, J. H., Kang, T. J., & Kim, D. M. (2008). Tuning the expression level of recombinant proteins by modulating mRNA stability in a cell-free protein synthesis system. *Biotechnology and bioengineering*, 101(2), 422-427.

Forster, A. C., Cornish, V. W., & Blacklow, S. C. (2004). Pure translation display. *Analytical biochemistry*, 333(2), 358-364.

Hino, M., Kataoka, M., Kajimoto, K., Yamamoto, T., Kido, J. I., Shinohara, Y., & Baba, Y. (2008). Efficiency of cell-free protein synthesis based on a crude cell extract from Escherichia coli, wheat germ, and rabbit reticulocytes. *Journal of biotechnology*, 133(2), 183-189.

Josephson, K., Hartman, M. C., & Szostak, J. W. (2005). Ribosomal synthesis of unnatural peptides. *Journal of the American Chemical Society*, 127(33), 11727-11735.

Keum, J. W., Ahn, J. H., Choi, C. Y., Lee, K. H., Kwon, Y. C., & Kim, D. M. (2006). The presence of a common downstream box enables the simultaneous expression of multiple proteins in an E. coli extract. *Biochemical and biophysical research communications*, 350(3), 562-567.

Keum, J. W., Ahn, J. H., Kang, T. J., & Kim, D. M. (2009). Combinatorial, selective and reversible control of gene expression using oligodeoxynucleotides in a cell-free protein synthesis system. *Biotechnology and bioengineering*, 102(2), 577-582.

Kim, T. W., Oh, I. S., Keum, J. W., Kwon, Y. C., Byun, J. Y., Lee, K. H., Choi, C. Y., & Kim, D. M. (2007). Prolonged cell-free protein synthesis using dual energy sources: Combined use of creatine phosphate and glucose for the efficient supply of ATP and retarded accumulation of phosphate. *Biotechnology and bioengineering*, 97(6), 1510-1515.

Kim, H. C., Kim, T. W., & Kim, D. M. (2011). Prolonged production of proteins in a cell-free protein synthesis system using polymeric carbohydrates as an energy source. *Process Biochemistry*, 46(6), 1366-1369.

Kigawa, T., Yamaguchi-Nunokawa, E., Kodama, K., Matsuda, T., Yabuki, T., Matsuda, N., Yabuki, T., Matsuda, N., Ishitani, R., Nureki, O., & Yokoyama, S. (2002). Selenomethionine incorporation into a protein by cell-free synthesis. *Journal of structural and functional genomics*, 2(1), 29-35.

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제품명	Cat. No
<i>AccuRapid</i> TM Cell-Free Protein Expression Kit	K-7250
<i>AccuRapid</i> TM Midi Protein Expression Kit	K-7260
<i>AccuRapid</i> TM Maxi Protein Expression Kit	K-7270
<i>AccuRapid</i> TM Protein Synthesis Kit	K-7280
<i>ExiProgen</i> TM ProXpress PCR Template Kit	K-7400
First primer F/R sets (N terminus 6x His tag) (each 5 nmole)	N-8229
First primer F/R sets (C terminus 6x His tag) (each 5 nmole)	N-8230
<i>ExiProgen</i> TM Protein Expression Optimization Kit	K-7410
pBIVT Vector Set-1	K-7350
<i>ExiProgen</i> TM EC Protein Synthesis Kit	K-7300
<i>ExiProgen</i> TM EC-Maxi Protein Synthesis Kit	K-7310
<i>ExiProgen</i> TM EC-Tagfree Protein Synthesis Kit	K-7320
<i>ExiProgen</i> TM EC-Disulfide Protein Synthesis Kit	K-7330
<i>ExiProgen</i> TM EC-Bulk Protein Synthesis Kit	K-7340
<i>ExiProgen</i> TM Dialysis Kit	K-7240
<i>ExiProgen</i> TM Consumable SET	KA-3001
Gene Synthesis Service	S-2041
Protein Synthesis Service	S-2500
<i>ExiProgen</i> TM	A-5041

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 Batch Code	 Biological Risks	 Catalog Number	 Caution
 Consult Instructions For Use	 Contains Sufficient for <n> tests	 Do not Re-use	 Manufacturer
 Research Use Only	 Temperature Limitation	 Use-by date	

BIONEER Corporation - HQ

Address 8-11 Munpyeongseo-ro, Daedeok-gu, Daejeon, 34302, Republic of Korea
E-mail sales@bioneer.co.kr
Web www.bioneer.com

BIONEER Global Center

Address 71, Techno 2-ro, Yuseong-gu, Daejeon, 34013, Republic of Korea
E-mail sales@bioneer.co.kr
Web www.bioneer.com

BIONEER R&D Center

Address Korea Bio Park BLDG #B-702, 700 Daewangpangyo-ro, Bundang-gu, Seongnam-si
Gyeonggi-do, 13488, Republic of Korea
E-mail sales@bioneer.co.kr
Web www.bioneer.com

BIONEER Inc. - USA Branch

Address 155 Filbert St. Suite 216 Oakland, CA 94607, USA
E-mail order.usa@bioneer.com
Web us.bioneer.com

BIONEER Corp. - European Branch

Address Ludwig-Erhard-Strasse 30-34, 65760 Eschborn, Germany
E-mail euinfo@bioneer.com
Web www.bioneer.com