

[Cat. No.] **S-6042-PH1**

## Introduction

BIONEER offers 96-well plates coated with 88 genes which are related with apoptosis invasion & metastasis, angiogenesis, oncogenes & tumor suppressor, signal transduction & transcription factors, and cell cycle & DNA damage repair, reference 5 genes, and 3 control primers. Human Cancer qPCR panel kit is an easy-to-use product as it simplifies preparation of real-time PCR mixture by making the user add the template DNA, 2X Master Mix (intercalating dye type), and nuclease-free water only. All primers are designed and validated in accordance with the MIQE (Minimum Information for Publication of Quantitative Real-Time PCR Experiments) guidelines\*.

\* Bustin, S.A., et al. 2009. The MIQE Guidelines: Minimum Information for Publication of Quantitative Real-Time PCR Experiments, Clinical Chemistry 55:4, 611-622.

## Features & Benefits

- Fidelity: qPCR primer with qPCR efficiency of 90-110%, with guaranteed detection limit of 100 copies.
- User-friendly: Simplified procedure starting just by adding the template and Master Mix you want to analyze.
- Economic: Reduced time and cost of primer design, synthesis and efficiency assays.

## Components

Components	Amount
Human Cancer qPCR panel kit (96 well plate)	1 plate
Adhesive Optical Sealing Film (Cat.No. 3111-4110)	1 sheet per plate

## Storage

- This product is shipped at ambient temperature.
- Store at room temperature. If stored in the recommended temperature, this product will be stable for 2 years after the delivery date.

## Plate Map

	1	2	3	4	5	6	7	8	9	10	11	12
A	ABL1	BAD	CASP8	CDKN2A	FGFR2	JAK2	MTA2	NTRK3	RB1	TFE3	TNF	ACTB
B	AKT1	BAX	CCND1	CFLAR	FGFR3	MAP2K1	MYC	PDGFRA	RET	TFEB	TNFRSF10B	B2M
C	ALK	BCL2	CCND2	CHEK2	FOS	MDM2	NF1	PDGFRB	RUNX1	TGFB1	TNFRSF1A	GAPDH
D	ANGPT1	BCL2L1	CCND3	CTNNB1	FOXO1A	MET	NF2	PLAU	SERPINB5	TGFBR1	TNFRSF25	HPRT1
E	ANGPT2	BCL6	CDK2	ERBB2	FOXO3A	KMT2A	NFKB1	PLAUR	SERPINE1	THBS1	TP53	RPLP0
F	APAF1	BRAF	CDK4	MECOM	GMPS	MMP2	NFKBIA	PTEN	SYK	TIMP1	TRADD	NTC
G	APC	BRCA1	CDK6	EWSR1	HMGA2	MMP9	NOTCH1	RAF1	TAL1	TIMP2	VEGFA	GDC
H	ATR	BRCA2	CDKN1A	FADD	IGF1	MTA1	NTRK1	RARA	TERT	TLX1	WT1	PPC

**Table 1. Layout of Human Cancer qPCR panel kit**

The panel is involved 88 target genes (A1 through H11), 5 reference genes (A12-E12), and 3 control primers (F12-H12).

## Online Resources



Korean



English

Visit our **product page** for additional information and protocols

## Ordering Information

Description		Cat. No.
AccuPower® qPCR Array System: Human Cancer qPCR panel kit	96 genes	S-6042-PH1
AccuPower® qPCR Array System: Single gene qPCR Primer Set	200 rxn	S-6042-S200

\* **Note:** The selected reference gene can be ordered in 200 rxn increments.

## Notice

BIONEER corporation reserves the right to make corrections, modifications, improvements and other changes to its products, services, specifications or product descriptions at any time without notice.

## Explanation of Symbols



Caution



Consult  
Instructions  
For Use







Do not  
Re-use



Use-by  
Date

## Experimental Procedures

Steps		Procedure Details																															
1	 <b>Preparation of reaction mixture</b>	<p>1. Prepare template DNA, AccuPower® 2X GreenStar™ qPCR Master Mix (K-6251, not provided), and nuclease-free water in a tube to make a total volume of 50 µl as described in following table.</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="text-align: left;">Components</th> <th style="text-align: left;">50 µl reaction</th> </tr> </thead> <tbody> <tr> <td>AccuPower® 2X GreenStar™ qPCR Master Mix</td> <td>25 µl</td> </tr> <tr> <td>Template DNA</td> <td>5 pg-100 ng</td> </tr> <tr> <td>Nuclease-free water</td> <td>Variable</td> </tr> <tr> <td>Total volume</td> <td>50 µl</td> </tr> </tbody> </table>	Components	50 µl reaction	AccuPower® 2X GreenStar™ qPCR Master Mix	25 µl	Template DNA	5 pg-100 ng	Nuclease-free water	Variable	Total volume	50 µl																					
Components	50 µl reaction																																
AccuPower® 2X GreenStar™ qPCR Master Mix	25 µl																																
Template DNA	5 pg-100 ng																																
Nuclease-free water	Variable																																
Total volume	50 µl																																
2	 <b>Resuspension of primers</b>	<p>2. Carefully remove the covered film of panel and dispense 50 µl of reaction mixture into each well of AccuPower® qPCR Array System: Human Cancer qPCR panel kit.  <b>* Note:</b> Change pipette tips following each pipetting step to avoid cross-contamination among the wells.</p> <p>3. Seal the plate with adhesive optical sealing film and briefly spin down.</p> <p>4. Then, completely mix by vortexing to resuspend lyophilized primers and spin down again.  <b>* Note:</b> Before start, check carefully if there are residues on the film.</p>																															
3	 <b>Real-time PCR</b>	<p>5. Perform the reaction under the following conditions.</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="text-align: left;">Step</th> <th style="text-align: left;">Temperature</th> <th style="text-align: left;">Time</th> <th style="text-align: left;">Cycles</th> </tr> </thead> <tbody> <tr> <td>Pre-denaturation</td> <td>95°C</td> <td>10 min</td> <td>1 cycle</td> </tr> <tr> <td>Denaturation</td> <td>95°C</td> <td>5 sec</td> <td></td> </tr> <tr> <td>Annealing</td> <td>58°C</td> <td>25 sec</td> <td rowspan="2">40 cycles</td> </tr> <tr> <td>Extension</td> <td>72°C</td> <td>30 sec</td> </tr> <tr> <td>Detection</td> <td></td> <td>Scan</td> <td></td> </tr> <tr> <td>Final extension</td> <td>65°C</td> <td>5 min</td> <td>1 cycle</td> </tr> <tr> <td>Melting</td> <td>65-95°C</td> <td>1 sec</td> <td>-</td> </tr> </tbody> </table> <p>6. After the reaction, perform data analysis.</p>	Step	Temperature	Time	Cycles	Pre-denaturation	95°C	10 min	1 cycle	Denaturation	95°C	5 sec		Annealing	58°C	25 sec	40 cycles	Extension	72°C	30 sec	Detection		Scan		Final extension	65°C	5 min	1 cycle	Melting	65-95°C	1 sec	-
Step	Temperature	Time	Cycles																														
Pre-denaturation	95°C	10 min	1 cycle																														
Denaturation	95°C	5 sec																															
Annealing	58°C	25 sec	40 cycles																														
Extension	72°C	30 sec																															
Detection		Scan																															
Final extension	65°C	5 min	1 cycle																														
Melting	65-95°C	1 sec	-																														
	 <b>Control primers</b>	<p><b>1) Non-Template Control (NTC)</b></p> <ul style="list-style-type: none"> <li>- NTC is a negative control for checking on random or on reagent contamination.</li> <li>- Just add the pre-mixture containing nuclease-free water but excluding the template into NTC well.</li> <li>- If the value of Ct<sup>NTC</sup> is less than 35, there is overall DNA contamination in your PCR system. In this case, clean up the equipment and replace all the reagents to new ones.</li> </ul> <p><b>2) Genomic DNA Control (GDC)</b></p> <ul style="list-style-type: none"> <li>- GDC primer is for the detection of non-transcribed genomic DNA contamination.</li> <li>- In GDC well, primers which target genomic DNA are coated.</li> <li>- Add pre-mixture (your template, 2X Master Mix and nuclease free water) into the GDC well.</li> <li>- If the value of Ct<sup>GDC</sup> is less than 35, gDNA contamination might have occurred in your RNA samples. In this case, you ought to conduct an additional DNase treatment to clean up your samples.</li> </ul> <p><b>3) Positive PCR Control (PPC)</b></p> <ul style="list-style-type: none"> <li>- PPC primer is for the PCR test.</li> <li>- The PPC well contains positive template and primers, so just add 2X Master Mix and nuclease-free water into the PPC well.</li> <li>- The value of Ct<sup>PPC</sup> should be referred to the quick manual provided together.</li> </ul>																															