

Bioneer Corporation
8-11 Munpyeongseo-ro, Daedeok-gu, Daejeon, 34302
Republic of Korea
Tel: +82-42-930-8777 (Korea : 1588-9788)
Fax: +82-42-930-8688 E-mail: sales@bioneer.com

Bioneer Inc.
155 Filbert St. Suite 216
Oakland, CA 94607, USA
Toll Free: +1-877-264-4300 Fax: +1-510-865-0350
E-mail: order.usa@bioneer.us.com

Bioneer R&D Center
Korea Bio Park BLDG #B-702, 700 Daewangpangyo-ro
Bundang-gu, Seongnam-si, Gyeonggi-do, 13488
Republic of Korea
Tel: +82-31-628-0500 Fax: +82-31-628-0555

Synthetic gene manual

1. Use it when you want to increase

- a. Add 20ul of DW or TE buffer into delivered DNA. (Final 100~250ng /ul)
- b. Prepare competent cell. (Commercial competent *E. coli* or self-made competent *E. coli*, DH5 α strain is recommended)
- c. Carefully transfer 1ul of plasmid DNA into the tube prepared in Step b.
- d. Place the tube on ice for 20 minutes.
- e. Heat-shock the cells for 60~90 seconds in a water bath at exactly 42°C.
- f. Immediately place the competent cell + DNA mixture back to ice for 3 minutes.
- g. Put the mixture into 10ml of LB+Ampicillin (50ug/ml~100ug/ml) culture medium.
- h. Incubate for overnight at 37°C with shaking (200rpm, about 16hr) and isolate plasmid DNA using *Accuprep*® Plasmid extraction kit.
- i. Restriction enzyme cut or PCR depending on the purpose and use after quantitative measurement of DNA.

1. Use it when you want to experiment

- a. Add 20ul of DW or TE buffer into delivered DNA. (Final 100~250ng /ul)
 - b. After quantitative measurement of DNA, restriction enzyme cut or PCR depending on the purpose and use immediately.
- ▶ Delivery form is 5ug of lyophilized plasmid DNA cloned in cloning site of pUC vector. (Can be store at room temperature)
 - ▶ After dissolving DNA in DW or TE buffer, it is recommended that place DNA at 4°C for 10 minutes.
 - ▶ Store at -20°C after adding DW or TE buffer.
 - ▶ Refrigeration of DNA can be a risk of degradation.