

## MagListo™ 5M Genomic DNA Extraction Kit (K-3603G)

### I Before You Begin

- 1) Completely dissolve one vial of **Proteinase K** in **1,250 µl** of nuclease free water. For short term storage, dissolved Proteinase K should be stored at 4°C. For long term storage, it is recommended to aliquot the enzyme into separate tubes and store at -20°C.
- 2) Completely dissolve one vial of **RNase A** in **600 µl** of nuclease-free water. For short term storage, dissolved RNase A should be stored at 4°C. For long term storage, it is recommended to aliquot the enzyme into separate tubes and store at -20°C.
- 3) If there is any precipitate in TL Buffer (Lysis), incubate on a heating block at 60°C.
- 4) Add correct amount of **absolute ethanol (not provided)** to **WM1 Buffer** and **WB2 Buffer**, respectively (see bottle label).
- 5) Pre-heat EA Buffer at 60°C before use.

### II DNA Extraction from Whole Blood and Buffy Coat

- 1) Add **20 µl** of **Proteinase K** to 1.5 ml tube.
- 2) Apply **200 µl** of **whole blood** and **buffy coat** to the tube containing proteinase K.  
(Note) If the sample volume is less than 200 µl, make the total volume 200 µl by adding 1X PBS (not provided).
- 3) Add **200 µl** of **GB Buffer** to the tube, mix thoroughly by vortexing, and incubate at 60°C for 10 min.
- 4) Add **400 µl** of **absolute ethanol (not provided)** and mix thoroughly by vortexing.
- 5) Add **100 µl** of **Magnetic NanoBead Solution** and mix thoroughly by vortexing.
- 6) Place the tube in **MagListo™-2 Magnetic Separation Rack** with the magnet plate attached and invert the rack gently, 3-4 times.
- 7) Without removing the tube from rack, remove the supernatant.
- 8) Detach the magnet plate from **MagListo™** stand. Add **500 µl** of **WM1 Buffer**.  
Mix thoroughly by vortexing.
- 9) Attach the magnet plate to **MagListo™** stand and invert the rack gently 3-4 times until the beads bind tightly to the magnet.
- 10) Without removing the tube from **MagListo™** rack, remove the supernatant.
- 11) Repeat the above step 8-10.
- 12) Detach the magnet plate from **MagListo™** stand. Add **700 µl** of **WB2 Buffer**.  
Mix thoroughly by vortexing, and repeat the above step 9-10.
- 13) Without removing the tube from **MagListo™**, add **700 µl** of **WE Buffer** to “the opposite side of bead”. Close the cap and gently invert the tube twice.
- 14) Discard the supernatant and completely remove the remaining supernatant by blotting.
- 15) Detach the magnet plate from **MagListo™** stand. Add **100 µl** of **EA Buffer**.  
Mix thoroughly by vortexing, and incubate at 60°C for 1 min.
- 16) Attach the magnet plate to **MagListo™** stand and invert the rack gently 3-4 times until the beads bind tightly to the magnet.
- 17) Without removing the tube from **MagListo™** rack, transfer the supernatant containing DNA carefully to a sterile microcentrifuge tube.

## MagListo™ 5M Genomic DNA Extraction Kit (K-3603G)

### III DNA Extraction from Cultured Cell

- 1) Centrifuge the cultured cells (~1x10<sup>6</sup> cells) for 10 min at 300 x g.
- 2) Resuspend the pellet in 200 µl of 1X PBS and place them to 1.5 ml tube.
- 3) Add 20 µl of **Proteinase K** to the tube.
- 4) Add 10 µl of **RNase A**. Incubate for 2 min at RT.
- 5) Add 200 µl of **GB Buffer**, mix thoroughly by vortexing, and incubate at 60°C for 10 min.
- 6) Add 400 µl of **absolute ethanol** and mix thoroughly by vortexing.
- 7) Add 100 µl of **Magnetic NanoBead Solution** and mix thoroughly by vortexing.
- 8) Place the tube in **MagListo™-2 Magnetic Separation Rack** with the magnet plate attached and invert the rack gently, 3-4 times.
- 9) Without removing the tube from **MagListo™** rack, remove the supernatant.
- 10) Detach the magnet plate from **MagListo™** stand. Add 700 µl of **WM1 Buffer**. Mix thoroughly by vortexing.
- 11) Attach the magnet plate to **MagListo™** stand and invert the rack gently 3-4 times until the beads bind tightly to the magnet.
- 12) Without removing the tube from **MagListo™** rack, remove the supernatant.
- 13) Go to **step 12** of "II. DNA Extraction from Whole Blood and Buffy Coat" in page 1 and continue the extraction process.

### IV DNA Extraction from Animal Tissue

- 1) Disrupt (or homogenize) the sample (~25 mg) then place them into 1.5 ml tube.
- 2) Add 180 µl of **TL Buffer** to the tube.
- 3) Add 20 µl of **Proteinase K** and mix thoroughly by vortexing.
- 4) Add 10 µl of **RNase A** and incubate for 2 min at RT.
- 5) Incubate at 60°C until the tissue is completely lysed.  
(Note) To increase the DNA purity, centrifuge the lysate at 13,000 rpm for 5 min and transfer the supernatant to 1.5 ml tube.
- 6) Add 200 µl of **GB Buffer** and mix thoroughly by vortexing.
- 7) Add 400 µl of **absolute ethanol (not provided)** and mix thoroughly by vortexing.
- 8) Add 100 µl of **Magnetic NanoBead Solution** and mix thoroughly by vortexing.
- 9) Place the tube in **MagListo™-2 Magnetic Separation Rack** with the magnet plate attached and invert the rack gently, 3-4 times.
- 10) Without removing the tube from rack, remove the supernatant.
- 11) Detach the magnet plate from **MagListo™** stand. Add 700 µl of **WM1 Buffer**. Mix thoroughly by vortexing.
- 12) Attach the magnet plate to **MagListo™** stand and invert the rack gently 3-4 times until the beads bind tightly to the magnet.
- 13) Without removing the tube from **MagListo™** rack, remove the supernatant.
- 14) Go to **step 12** of "II. DNA Extraction from Whole Blood and Buffy Coat" in page 1 and continue the extraction process.

※ For more information, please visit [www.bioneer.com](http://www.bioneer.com) and refer to the User Guide of this kit.

## MagListo™ 5M Genomic DNA Extraction Kit (K-3603G)

### V DNA Extraction from Bacterial Cells (Gram-Negative Bacteria)












- 1) Harvest up to  $1 \times 10^9$  bacterial cells by centrifugation at  $6000 \times g$  for 10 minutes to pellet cells. Discard the supernatant carefully with a pipette.
- 2) Resuspend the cell pellet in 180  $\mu$ l of TL Buffer by vortexing or pipetting. Transfer the cell suspension into the clean tube.
- 3) Go to **step 3** of "IV. DNA Extraction from Animal Tissue" in page 2 and continue the extraction process.

### VI DNA Extraction from Bacterial Cells (Gram-Positive Bacteria)

- 1) Harvest up to  $1 \times 10^9$  bacterial cells by centrifugation at  $6000 \times g$  for 10 minutes to pellet cells. Discard the supernatant carefully with a pipette.
- 2) Resuspend the cell pellet in 180  $\mu$ l of lysis buffer (not provided) by vortexing or pipetting. Transfer the cell suspension into clean tube.  
(Note) Lysis buffer for gram-positive bacteria can be prepared by using this formulation:  
20 mM Tris-HCl (pH 8.0), 2 mM sodium EDTA, and 1.2% Triton® X-100.
- 3) Add 20  $\mu$ l of lysozyme (100 mg/ml, not provided) and mix well by vortexing.
- 4) If RNA-free genomic DNA is required, add up to 10  $\mu$ l of RNase A to the sample and gently mix.
- 5) Incubate at 37°C for 30 minutes.
- 6) Add 20  $\mu$ l of Proteinase K and mix well by vortexing.
- 7) Add 200  $\mu$ l of GB Buffer and mix well by vortexing.
- 8) Incubate at 60°C for 30 minutes or until bacterial cells are completely lysed.
- 9) Go to **step 7** of "IV. DNA Extraction from Animal Tissue" in page 2 and continue the extraction process.

※ For more information, please visit [www.bioneer.com](http://www.bioneer.com) and refer to the User Guide of this kit.

# MagListo™ 5M Genomic DNA Extraction Kit (K-3603G)

Step	Image	Description
Lysis		Add Reagents according to sample type <b>1. Whole blood and Buffy coat</b> - Sample 200 µl - Proteinase K 20 µl, GB Buffer 200 µl
		<b>2. Cultured cell</b> - 1X PBS 200 µl + Cell pellet (~1x10 <sup>6</sup> ) - Proteinase K 20 µl, RNase A 10 µl, GB Buffer 200 µl
		<b>3. Animal tissue</b> - Disrupt sample ~25 mg - Proteinase K 20 µl, TL Buffer 180 µl, RNase A 10 µl
		<b>4. Bacterial cells (Gram-Negative, Positive Bacteria)</b> - ~1x10 <sup>9</sup> cells - Gram-Negative: TL Buffer 180 µl, Proteinase K 20 µl, RNase A 10 µl - Gram-Positive: Lysis Buffer 180 µl, Lysozyme 20 µl, RNase A 10 µl, Proteinase K 20 µl, GB Buffer 200 µl
Incubation		60°C Heating block for 10 min (For Animal tissue and Bacterial cells, incubate at, 60°C Heating block until it is completely lysed)
Precipitation		Add absolute ethanol 400 µl and mix (For Animal tissue and Gram-Negative Bacteria, add GB Buffer 200 µl and absolute ethanol 400 µl and mix)
Binding		Add Magnetic NanoBead 100 µl and mix until the beads are fully resuspended
		Follow these 3 Magnet Attach/Detach Step 1. Attach magnet 2. Discard the supernatant 3. Detach magnet
1 <sup>st</sup> Wash		Add WM1 Buffer 700 µl and mix until the beads are fully resuspended (For Whole blood, 1 <sup>st</sup> washing with WM1 Buffer 500 µl twice)
		Repeat the above magnet attach/detach step (step 1, 2 and 3)
2 <sup>nd</sup> Wash		Add WB2 Buffer 700 µl and mix until the beads are fully resuspended
		Repeat the above magnet attaching step (step 1 and 2)
3 <sup>rd</sup> Wash		Add WE Buffer 700 µl to "the opposite side of the bead pellet" Close the cap and gently invert the rack twice
		Repeat the above magnet detaching step (step 2 and 3)
Elution		Add EA Buffer 100 µl and mix
		60°C Heating block for 1 min