




Experimental Procedures

Steps		Procedure Details																											
Resuspension Protocol																													
1	 Resuspension of primer	<ol style="list-style-type: none"> Briefly centrifuge tubes containing primers to ensure that the pellet is located at the bottom of the tube. Dissolve primers in nuclease-free water or 0.1X TE buffer (1 mM Tris, 0.1 mM EDTA) as follows: for 100 reaction products, add 250 µl, and for 200 reaction products, add 500 µl, respectively. Mix the primers by pipetting or vortexing briefly and spin down. Store at -20°C in small aliquots and avoid repeated freeze and thaw cycles. 																											
Example of Real-time PCR Protocol																													
<ul style="list-style-type: none"> AccuTarget™ Real-Time PCR Primer Library is verified with the following protocol using BIONEER's AccuPower® GreenStar™ qPCR PreMix. 																													
2	 Preparation of reaction mixture	<ol style="list-style-type: none"> Add template DNA, primers, nuclease-free water, and 50X ROX dye (optional, not provided) into AccuPower® GreenStar™ qPCR PreMix (Cat. No. K-6200, not provided) tubes to make a total volume of 50 µl. Do not calculate the dried pellet. <ul style="list-style-type: none"> Preparation of reaction mixture <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="text-align: center;">Components</th> <th style="text-align: center;">50 µl reaction</th> </tr> </thead> <tbody> <tr> <td>Template DNA (40-200 ng)</td> <td style="text-align: center;">Variable</td> </tr> <tr> <td>Forward primer</td> <td style="text-align: center;">2.5 µl</td> </tr> <tr> <td>Reverse primer</td> <td style="text-align: center;">2.5 µl</td> </tr> <tr> <td>(Optional) 50X ROX dye*</td> <td style="text-align: center;">1 µl</td> </tr> <tr> <td>Nuclease-free water</td> <td style="text-align: center;">Variable</td> </tr> <tr> <td>Total volume</td> <td style="text-align: center;">50 µl</td> </tr> </tbody> </table> <p>* Note: ROX dye is used for normalization of intensity by background subtraction. The use of ROX dye is recommended for Applied Biosystems 7500 Real-Time PCR System, but not required for BIONEER Exicycler™ 96 Real-Time PCR System.</p> Seal real-time PCR tubes with an optical adhesive film (Cat. No. 3111-4110). Mix the reaction mixture by vortexing and briefly spin down. 	Components	50 µl reaction	Template DNA (40-200 ng)	Variable	Forward primer	2.5 µl	Reverse primer	2.5 µl	(Optional) 50X ROX dye*	1 µl	Nuclease-free water	Variable	Total volume	50 µl													
Components	50 µl reaction																												
Template DNA (40-200 ng)	Variable																												
Forward primer	2.5 µl																												
Reverse primer	2.5 µl																												
(Optional) 50X ROX dye*	1 µl																												
Nuclease-free water	Variable																												
Total volume	50 µl																												
3	 Real-time PCR	<ol style="list-style-type: none"> Perform the reaction under the following conditions. <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="text-align: center;">Step</th> <th style="text-align: center;">Temperature</th> <th style="text-align: center;">Time</th> <th style="text-align: center;">Cycles</th> </tr> </thead> <tbody> <tr> <td>Pre-denaturation</td> <td style="text-align: center;">94°C</td> <td style="text-align: center;">5 min</td> <td style="text-align: center;">1 cycle</td> </tr> <tr> <td>Denaturation</td> <td style="text-align: center;">95°C</td> <td style="text-align: center;">10 sec</td> <td></td> </tr> <tr> <td>Annealing</td> <td style="text-align: center;">58°C</td> <td style="text-align: center;">25 sec</td> <td rowspan="2" style="text-align: center;">40 cycles</td> </tr> <tr> <td>Extension</td> <td style="text-align: center;">72°C</td> <td style="text-align: center;">30 sec</td> </tr> <tr> <td>Detection (Scan)</td> <td></td> <td></td> <td></td> </tr> <tr> <td>Melting</td> <td style="text-align: center;">-</td> <td></td> <td style="text-align: center;">1 cycle</td> </tr> </tbody> </table> After the reaction is completed, analyze the results. 	Step	Temperature	Time	Cycles	Pre-denaturation	94°C	5 min	1 cycle	Denaturation	95°C	10 sec		Annealing	58°C	25 sec	40 cycles	Extension	72°C	30 sec	Detection (Scan)				Melting	-		1 cycle
Step	Temperature	Time	Cycles																										
Pre-denaturation	94°C	5 min	1 cycle																										
Denaturation	95°C	10 sec																											
Annealing	58°C	25 sec	40 cycles																										
Extension	72°C	30 sec																											
Detection (Scan)																													
Melting	-		1 cycle																										