




Experimental Procedures

Steps		Procedure Details																								
Resuspension Protocol																										
1	 <p>Resuspension of primer</p>	<ol style="list-style-type: none"> Briefly centrifuge tubes containing primers to ensure that the pellet is located at the bottom of the tube. Dissolve primers to a convenient stock concentration in nuclease-free water or 0.1X TE buffer (1 mM Tris, 0.1 mM EDTA). Mix the primers by pipetting or vortexing briefly and spin down. Store at -20°C in small aliquots and avoid repeated freeze and thaw cycles. 																								
Real-time PCR Protocol																										
<p>• This recommended protocol can be modified to get optimal results, according to real-time PCR instrument and target DNA sequences.</p>																										
1	 <p>Preparation of reaction mixture</p>	<ol style="list-style-type: none"> Add template DNA, primers, and nuclease-free water into <i>AccuPower</i>® GreenStar™ qPCR PreMix (K-6210, not provided) tubes to a total volume of 20 µl or 50 µl. Do not calculate the dried pellet. <ul style="list-style-type: none"> Preparation of reaction mixture <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="text-align: center;">Components</th> <th style="text-align: center;">20 µl reaction</th> <th style="text-align: center;">50 µl reaction</th> </tr> </thead> <tbody> <tr> <td>Template DNA (40-200 ng)</td> <td style="text-align: center;">Variable</td> <td style="text-align: center;">Variable</td> </tr> <tr> <td>Forward primer (10 pmol/µl)</td> <td style="text-align: center;">1-2 µl</td> <td style="text-align: center;">2-5 µl</td> </tr> <tr> <td>Reverse primer (10 pmol/µl)</td> <td style="text-align: center;">1-2 µl</td> <td style="text-align: center;">2-5 µl</td> </tr> <tr> <td>Nuclease-free water</td> <td style="text-align: center;">Variable</td> <td style="text-align: center;">Variable</td> </tr> <tr> <td>Total volume</td> <td style="text-align: center;">20 µl</td> <td style="text-align: center;">50 µl</td> </tr> </tbody> </table> Seal real-time PCR tubes with an optical adhesive film (3111-4110, not provided). Mix the reaction mixture by vortexing and briefly spin down. 	Components	20 µl reaction	50 µl reaction	Template DNA (40-200 ng)	Variable	Variable	Forward primer (10 pmol/µl)	1-2 µl	2-5 µl	Reverse primer (10 pmol/µl)	1-2 µl	2-5 µl	Nuclease-free water	Variable	Variable	Total volume	20 µl	50 µl						
Components	20 µl reaction	50 µl reaction																								
Template DNA (40-200 ng)	Variable	Variable																								
Forward primer (10 pmol/µl)	1-2 µl	2-5 µl																								
Reverse primer (10 pmol/µl)	1-2 µl	2-5 µl																								
Nuclease-free water	Variable	Variable																								
Total volume	20 µl	50 µl																								
2	 <p>Real-time PCR</p>	<ol style="list-style-type: none"> Perform the reaction under the following conditions. <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="text-align: center;">Step</th> <th style="text-align: center;">Temperature</th> <th style="text-align: center;">Time</th> <th style="text-align: center;">Cycles</th> </tr> </thead> <tbody> <tr> <td>Pre-denaturation</td> <td style="text-align: center;">95°C</td> <td style="text-align: center;">1-5 min</td> <td style="text-align: center;">1 cycle</td> </tr> <tr> <td>Denaturation</td> <td style="text-align: center;">95°C</td> <td style="text-align: center;">5-20 sec</td> <td></td> </tr> <tr> <td>Annealing/Extension</td> <td style="text-align: center;">55-60°C</td> <td style="text-align: center;">40-45 sec</td> <td style="text-align: center;">40-45 cycles</td> </tr> <tr> <td>Detection (Scan)</td> <td></td> <td></td> <td></td> </tr> <tr> <td>Melting</td> <td style="text-align: center;">-</td> <td></td> <td style="text-align: center;">1 cycle</td> </tr> </tbody> </table> <p>* Note: For maximum yield and specificity, temperatures and cycling times should be optimized for each new template DNA or primers.</p> 	Step	Temperature	Time	Cycles	Pre-denaturation	95°C	1-5 min	1 cycle	Denaturation	95°C	5-20 sec		Annealing/Extension	55-60°C	40-45 sec	40-45 cycles	Detection (Scan)				Melting	-		1 cycle
Step	Temperature	Time	Cycles																							
Pre-denaturation	95°C	1-5 min	1 cycle																							
Denaturation	95°C	5-20 sec																								
Annealing/Extension	55-60°C	40-45 sec	40-45 cycles																							
Detection (Scan)																										
Melting	-		1 cycle																							