




Experimental Procedures

Steps		Procedure Details																																							
1	 Dilute and mix oligonucleotides	<ol style="list-style-type: none"> Dilute oligonucleotides with Tris buffer (10 mM Tris, 0.1 mM EDTA, 50 mM NaCl, pH 8.0) or phosphate buffer (100 mM sodium phosphate, 150 mM NaCl, 0.1 mM EDTA, pH 7.5 or 8.0). * Note: Some salts are required for the hybridization of oligonucleotides. Mix the two sequences together in equal molar amounts. * Note: If different amounts are used, there will always be single-stranded sequences left over. 																																							
2	 Option 1. Anneal with a water bath or heating block	<ol style="list-style-type: none"> Incubate oligonucleotides at 95°C for 5 min. Turn off the hotplate of water bath (or heating block) and allow oligonucleotides to slowly cool to room temperature. Store at -20°C in small aliquots. 																																							
	 Option 2. Anneal with a thermal cycler	<ol style="list-style-type: none"> Perform the reaction under the following conditions. <ul style="list-style-type: none"> Simple protocol <table border="1" style="width: 100%; border-collapse: collapse; text-align: center;"> <thead> <tr> <th>Step</th> <th>Temperature</th> <th>Time</th> <th>Cycles</th> </tr> </thead> <tbody> <tr> <td>Step 1</td> <td>95°C</td> <td>5 min</td> <td>1 cycle</td> </tr> <tr> <td>Step 2</td> <td>95°C (-1°C/cycle)</td> <td>1 min/cycle</td> <td>70 cycles</td> </tr> <tr> <td>Step 3</td> <td>4°C</td> <td></td> <td>Hold</td> </tr> </tbody> </table> Advanced protocol (Example in which the oligonucleotide pair has a T_m of 55°C) <table border="1" style="width: 100%; border-collapse: collapse; text-align: center;"> <thead> <tr> <th>Step</th> <th>Temperature</th> <th>Time</th> <th>Cycles</th> </tr> </thead> <tbody> <tr> <td>Step 1</td> <td>95°C</td> <td>5 min</td> <td>1 cycle</td> </tr> <tr> <td>Step 2</td> <td>95°C (-1°C/cycle)</td> <td>1 min/cycle</td> <td>40 cycles*</td> </tr> <tr> <td>Step 3</td> <td>55°C</td> <td>30 min</td> <td>1 cycle</td> </tr> <tr> <td>Step 4</td> <td>55°C (-1°C/cycle)</td> <td></td> <td>20 cycles*</td> </tr> <tr> <td>Step 5</td> <td>4°C</td> <td></td> <td>Hold</td> </tr> </tbody> </table> * The number of cycles in step 2 and 4 can be according to T_m value of oligonucleotides to be annealed. Store at -20°C in small aliquots. 	Step	Temperature	Time	Cycles	Step 1	95°C	5 min	1 cycle	Step 2	95°C (-1°C/cycle)	1 min/cycle	70 cycles	Step 3	4°C		Hold	Step	Temperature	Time	Cycles	Step 1	95°C	5 min	1 cycle	Step 2	95°C (-1°C/cycle)	1 min/cycle	40 cycles*	Step 3	55°C	30 min	1 cycle	Step 4	55°C (-1°C/cycle)		20 cycles*	Step 5	4°C	
Step	Temperature	Time	Cycles																																						
Step 1	95°C	5 min	1 cycle																																						
Step 2	95°C (-1°C/cycle)	1 min/cycle	70 cycles																																						
Step 3	4°C		Hold																																						
Step	Temperature	Time	Cycles																																						
Step 1	95°C	5 min	1 cycle																																						
Step 2	95°C (-1°C/cycle)	1 min/cycle	40 cycles*																																						
Step 3	55°C	30 min	1 cycle																																						
Step 4	55°C (-1°C/cycle)		20 cycles*																																						
Step 5	4°C		Hold																																						