

MagListo™ 5M Forensic Sample DNA Extraction Kit (K-3615)

I Before You Begin

- 1) Completely dissolve one vial of **Proteinase K** in **1,250 µl** of nuclease free water. For short term storage, dissolved Proteinase K should be stored at **4°C**. For long term storage, it is recommended to aliquot the enzyme into separate tubes and store at **-20°C**.
- 2) Completely dissolve **Poly(A)** in **1,000 µl** of nuclease free water. Divide the solution into aliquots of appropriate amount and store the aliquots at **-20°C**. Repeating of freeze-thawing of the aliquots more than 3 times will result in the degradation of RNA.

II DNA Extraction from Forensic Sample

※ Handling samples

Dried body fluid spot or fingerprint (FTA card, paper, cloth, etc.): Punch out ~7 mm or cut out ~2 cm²

Hair: Cut the hair 1 cm length from the hair root

Bone & teeth: Homogenize the bone or teeth ~100 mg

Chewing gum: Cut out ~30 mg

Cigarette butts: Cut out ~2 cm² from the end of the cigarette butt

Busccal swab: Cut the swab from its stick by hand or scissors

- 1) Apply the forensic sample to a 1.5 ml tube.
- 2) Add **300 µl** of **TL Buffer** and **10 µl** of **Proteinase K**, and mix thoroughly by vortexing.
(Optional) Add 20 µl of 1 M DTT and mix thoroughly by vortexing. If the sample is hair, hair clipping or semen stains, this step is necessary to increase yield.
- 3) Incubate at 60°C for at least 1 hour.
- 4) Add **300 µl** of **GB Buffer**, and mix thoroughly by vortexing and incubate at 60°C for 20 min.
(Optional) If the sample contains low number of cell (< 1×10⁶) or has a small amount of DNA (<100 ng), we recommend to add **2 µl** of dissolved **Poly(A)** solution to each sample after addition of GB Buffer.
- 5) Centrifuge at 13,000 rpm for 1 min and transfer the supernatant only into a new 1.5 ml tube.
- 6) Add **600 µl** of **absolute ethanol** to the tube containing the supernatant, and mix thoroughly by vortexing.
- 7) Add **100 µl** of **Magnetic NanoBead solution**, and mix thoroughly by vortexing.
(Note) Please vortex Magnetic NanoBead Solution well before use.
- 8) Place the tube in **MagListo™-2 Magnetic Separation Rack** with the magnet plate attached and invert the rack gently 3~4 times.
- 9) Without removing the tube from **MagListo™** rack, remove the supernatant.
- 10) Detach the magnet plate from **MagListo™** stand. Add **700 µl** of **WM1 Buffer**.
Mix thoroughly by vortexing.
- 11) Attach the magnet plate to **MagListo™** stand and invert the rack gently 3~4 times until the beads bind tightly to the magnet.

MagListo™ 5M Forensic Sample DNA Extraction Kit (K-3615)

- Without removing the tube from *MagListo*™ rack, remove the supernatant.
- Detach the magnet plate from *MagListo*™ stand. Add **700 µl** of **W2 Buffer**.
Mix thoroughly by vortexing, and repeat the above step 11~12.
- Without removing the tube from *MagListo*™ rack, add **700 µl** of **WE Buffer** to “the opposite side of bead”. Close the cap and gently invert the tube twice.
- Discard the supernatant and completely remove the remaining supernatant by blotting.
- Detach the magnet plate from *MagListo*™ stand. Add **30~80 µl** of **EA Buffer**.
Mix thoroughly by vortexing and incubate at 60°C for 1 min.
- Vortex the tubes for 15 sec, then place them on the *MagListo*™ Magnetic Separation Rack. Attach the magnet plate to *MagListo*™ stand and invert the rack gently 3~4 times until the beads bind tightly to the magnet.
- Without removing the tube from *MagListo*™ rack, transfer the supernatant containing DNA carefully to a sterile microcentrifuge tube.

III DNA Extraction from Small Amount of Tissue

- Disrupt (or homogenize) **the sample (~10 mg)** then place them into a 1.5 ml tube.
- Add **90 µl** of **TL Buffer** and **10 µl** of **Proteinase K solution**, mix thoroughly by vortexing and incubate at 60°C until the sample is lysed.
(Optional) If RNA-free DNA is required, add up to 200 µg of RNase A (KB-0101, not provided) and incubate for 2 min at RT.
- Add **100 µl** of **GB Buffer**, and mix thoroughly by vortexing.
- Add **200 µl** of **absolute ethanol** and mix thoroughly by vortexing.
- Go to **step 7** of “II. DNA Extraction from Forensic Sample” in page 1 and continue extraction process.

IV DNA Extraction from Small Volumes of Blood and Saliva

- Add **10 µl** of **Proteinase K solution** to a 1.5 ml tube.
- Apply **1~100 µl** of **blood or saliva**.
(Note) If the sample volume is less than 100 µl, make the total volume 100 µl by adding 1X PBS (not provided).
- Add **100 µl** of **GB Buffer**, mix thoroughly by vortexing and incubate at 60°C for 10 min.
- Add **200 µl** of **absolute ethanol** and mix thoroughly by vortexing.
- Go to **step 7** of “II. DNA Extraction from Forensic Sample” in page 1 and continue extraction process.

※ For more information, please visit www.bioneer.com and refer to the User's Guide of this kit.

MagListo™ 5M Forensic Sample DNA Extraction Kit (K-3615)

Before You Begin

- 1) Completely dissolve **Proteinase K** in **1,250 µl** of nuclease free water.
- 2) Completely dissolve **Poly(A)** in **1,000 µl** of nuclease free water.

Step	Image	Description
Lysis		Add Reagents according to sample type and resuspend the sample completely 1. Forensic samples - FTA card, Paper, Cloth etc. ~2 cm ² , Hair 1 cm, Bone & Teeth ~100 mg, Chewing gum ~30 mg, Cigarette butts ~2 cm ² , Busccal swab - Proteinase K 10 µl, TL Buffer 300 µl 2. Small amount of tissue - Sample ~10 mg - Proteinase K 10 µl, TL Buffer 90 µl 3. Small volumes of blood & saliva - Sample 100 µl - Proteinase K 10 µl, GB Buffer 100 µl
		1. For forensic sample , 60°C for at least 1 hour 2. For tissue , 60°C until the sample is completely lysed 3. For blood & saliva , 60°C for 10 min
		1. For forensic sample , add GB Buffer 300 µl , mix and 60°C for 20 min 2. For tissue , add GB Buffer 100 µl and mix
Precipitation		1. For forensic sample , centrifuge at 13,000 rpm for 1 min and transfer the supernatant into a new tube. Add absolute ethanol 600 µl and mix 2. For tissue , add absolute ethanol 200 µl and mix 3. For blood & saliva , add absolute ethanol 200 µl and mix
Binding		Add Magnetic NanoBead 100 µl and mix until the beads are fully resuspended
		Follow these 3 Magnet Attach/Detach Step 1. Attach magnet 2. Discard the supernatant 3. Detach magnet
1 st Wash		Add WM1 Buffer 700 µl and mix until the beads are fully resuspended Repeat the above magnet attach/detach step (step 1, 2 and 3)
2 nd Wash		Add W2 Buffer 700 µl and mix until the beads are fully resuspended Repeat the above magnet attach step (step 1 and 2)
3 rd Wash		Add WE Buffer 700 µl to "the opposite side of the bead pellet" Close the cap and gently invert the rack twice
		Repeat the above magnet detach step (step 2 and 3)
Elution		Add EA Buffer 30~80 µl and mix
		60°C Heating block for 1 min