

MagListo™ cfDNA Extraction Kit (K-3619)

I Before You Begin

- 1) Completely dissolve one vial of **Proteinase K** in **1,250 µl** of nuclease free water. For short term storage, dissolved Proteinase K should be stored at **4°C**. For long term storage, it is recommended to aliquot the enzyme into separate tubes and store at **-20°C**.
- 2) CL Buffer, CB buffer and CW1 buffer contain chaotropic salt. You should take the appropriate laboratory safety precautions and wear gloves when handling.
- 3) Prepare isopropanol.

II Pretreatment of Samples

- a. Plasma from whole blood
 - 1) Centrifuge the blood samples at **2,000 x g** for 10 min.
 - 2) Transfer the supernatant to a new 50 ml tube.
 - 3) Repeat the above step 1.
- b. Serum from whole blood
 - 1) Allow the blood samples to clot for 30 min ~ 1 hour at room temperature before centrifugation.
 - 2) Centrifuge the samples at **2,000 x g** for 10 min.
 - 3) Transfer the supernatant to a new 50 ml tube.
- c. Urine and Saliva
 - 1) Centrifuge the urine samples at **3,000 x g** for 30 min.
 - 2) Transfer the supernatant to a 50 ml tube.

III cfDNA Extraction from Plasma/Serum/Urine/Saliva

- 1) Transfer the collected samples (~5 ml) to a 50 ml tube.
- 2) Add **40 µl** of **Proteinase K** solution per 1 ml of sample and mix well.
- 3) Add **same sample volume** of **CL Buffer** to the collected sample tube.
- 4) Incubate at **60°C** for 10 min.
- 5) Add **375 µl** of **isopropanol** per 1 ml of sample to the tube and mix well.
- 6) Add **200 µl** of **Magnetic NanoBeads solution** to the tube and mix thoroughly by vortex mixer until the beads are fully resuspended.
(Note) Please vortex Magnetic NanoBeads solution well before use.
- 7) After incubation plate the tube for 5 min in the ice or 10 min at room temperature.
- 8) Place the tube in **MagListo™-50 Magnetic Separation Rack** with the magnet plate attached and invert the rack gently 3~4 times.
- 9) Without removing the tubes from rack, carefully pour the supernatant out and completely remove the remaining supernatant using paper towel by blotting.

MagListo™ cfDNA Extraction Kit (K-3619)










- 10) Detach the magnet plate from *MagListo*™ stand. Add **1 ml** of **CB Buffer** to the tube.
Close the cap and mix by vortex mixer until the beads are fully resuspended.
- 11) Take the CB buffer mixed with magnetic nano beads and transfer into the new 1.5 ml or 2 ml tube.
- 12) Place the tube in *MagListo*™-2 **Magnetic Separation Rack** with the magnet plate attached and invert the rack 3~4 times gently until the beads bind tightly to magnet.
- 13) Without removing the tubes from rack, discard supernatant carefully and completely remove the remaining supernatant onto a paper towel by blotting down.
- 14) Detach the magnet plate from *MagListo*™ stand. Add **1 ml** of **CW1 Buffer**, close the cap and mix thoroughly by vortexing. Then repeat the above step 12-13.
- 15) Detach the magnet plate from *MagListo*™ stand. Add **1 ml** of **TA Buffer**, close the cap and mix thoroughly by vortexing. Then repeat the above step 12-13.
- 16) Detach the magnet plate from *MagListo*™ stand. Add **30 µl (1 ml) / 50 µl (3 ml) / 70 µl (5 ml)** of **EA Buffer** to the tube and resuspend by pipetting.
- 17) Incubate the tubes at 60°C for 1 min and vortex the tubes for 15 sec.
- 18) Centrifuge at 1,000 rpm for 1 min.
- 19) Carefully transfer the supernatant containing cfDNA to a sterile tube.

※ For more information, please visit www.bioneer.com and refer to the User's Guide of this kit.

MagListo™ cfDNA Extraction Kit (K-3619)

Before You Begin

- 1) Completely dissolve one vial of **Proteinase K** in **1,250 µl** of nuclease free water.
- 2) Prepare Isopropanol.

| Step | Image | Description |
|----------------------|---|--|
| Lysis |  | Add Reagents according to sample type 1. Serum, Saliva, Plasma, Urine - Collected sample ~5 ml - Proteinase K 40 µl per 1 ml of sample - CL Buffer with equal volume to sample |
| | | Incubate at 60°C for 10 min |
| Precipitation |  | Add Isopropanol 375 µl per 1 ml of sample and mix |
| Binding |  | Add Magnetic NanoBead 200 µl and mix until the beads are fully resuspended Incubate at RT for 10 min or on ice for 5 min |
| |  | Follow these 3 steps 1. Attach magnet 2. Discard the supernatant 3. Detach magnet |
| 1 st Wash |  | Add CB Buffer 1 ml and mix until the beads are fully resuspended |
| Transfer |  | Transfer into new 1.5 ml or 2 ml tube |
| | | Repeat the above magnet attach/detach step (Step 1,2 and 3) |
| 2 nd Wash |  | Add CW1 Buffer 1 ml and mix until the beads are fully resuspended |
| | | Repeat the above magnet attach/detach step (Step 1,2 and 3) |
| 3 rd Wash |  | Add TA Buffer 1 ml and mix until the beads are fully resuspended |
| | | Repeat the above magnet attach/detach step (Step 1,2 and 3) |
| Elution |  | Add EA Buffer 30 µl (1 ml) / 50 µl (3 ml) / 70 µl (5 ml) and mix |