

MagListo™ Genomic DNA Extraction Kit for ExiPrep™96 Lite

I Before You Begin

- 1) Completely dissolve one vial of **Proteinase K** in **1,250 µl** of nuclease-free water. For short term storage, dissolved Proteinase K should be stored at **4°C**. For long term storage, it is recommended to aliquot the enzyme into separate tubes and store at **-20°C**.
- 2) Completely dissolve **RNase A** in **600 µl** of nuclease-free water. For short term storage, dissolved RNase A should be stored at **4°C**.

II Sample Preparation

A. DNA Extraction from Cultured Cells

- 1) Harvest cultured cells (10^4 - 10^6 cells) by centrifugation at **300 x g** for 5 min to pellet cells. Discard the supernatant carefully without disturbing a cell pellet.
- 2) Completely resuspend the cell pellet in **200 µl** of **1x PBS buffer**.
- 3) Add **20 µl** of **Proteinase K** to the sample.
- 4) Add **10 µl** of **RNase A** to the sample, mix thoroughly, and incubate for 2 min at room temperature.
- 5) Add **200 µl** of **GB Buffer** to the sample and mix immediately by vortexing.
- 6) Incubate at **60°C** for 10 min.
- 7) Transfer the lysate to a new 96-well dome plate.

B. DNA Extraction from Whole Blood and Buffy Coat

- 1) Apply **200 µl** of whole blood or buffy coat sample to a clean 1.5 ml tube.
- 2) Add **20 µl** of **Proteinase K** to the sample.
- 3) Add **200 µl** of **GB Buffer** to the sample and mix immediately by vortexing.
- 4) Incubate at **60°C** for 10 min.
- 5) Transfer the lysate to a new 96-well dome plate.

C. DNA Extraction from Animal Tissue

- 1) Homogenize **25 – 50 mg** of fresh or thawed animal tissue sample with a homogenizer and place them into a clean 1.5 ml tube.
- 2) Add **200 µl** of **TL Buffer** and **20 µl** of **Proteinase K** and **RNase A** to the sample and mix thoroughly by vortexing.
- 3) Incubate at **60°C** until the sample has been completely lysis (1 hour).
- 4) Add **200 µl** of **GB Buffer** to the lysate and mix immediately by vortexing.
- 5) Transfer the lysate to a new 96-well dome plate.

D. DNA Extraction from Bacterial Cells (Gram-Negative Bacteria)

- 1) Harvest up to 1×10^9 bacterial cells by centrifugation at **8,000 rpm** for 5 min to pellet cells. Discard the supernatant carefully with a pipette.
- 2) Completely resuspend the cell pellet in **180 µl** of **TL Buffer** by vortexing or pipetting.
- 3) Add **20 µl** of **Proteinase K** and **RNase A** and mix thoroughly by vortexing.

MagListo™ Genomic DNA Extraction Kit for ExiPrep™96 Lite

- 4) Incubate at 60°C for 1 hour or until bacterial cells are completely lysed.
- 5) Add **200 µl** of **GB Buffer** to the sample and mix immediately by vortexing.
- 6) Transfer the lysate to a new 96-well dome plate.

E. DNA Extraction from Bacterial Cells (Gram-Positive Bacteria)

- 1) Harvest up to 1×10^9 bacterial cells by centrifugation at 8,000 rpm for 5 min to pellet cells. Discard the supernatant carefully with a pipette.
- 2) Completely resuspend the cell pellet in **180 µl** of **Lysis Buffer** (for Gram-Positive bacteria, not provided) by vortexing or pipetting.
(Note) Lysis Buffer for Gram-Positive bacteria: 20 mM Tris-HCl (pH 8.0), 2 mM EDTA, 1.2% Triton® X-100
- 3) Add **20 µl** of **lysozyme (100 mg/ml)**, not provided) and **10 µl** of **RNase A** and mix thoroughly by vortexing.
- 4) Incubate at 37°C for 30 min.
- 5) Add **20 µl** of **Proteinase K** and mix immediately by vortexing.
- 6) Add **200 µl** of **GB Buffer** and mix immediately by vortexing.
- 7) Incubate at 60°C for 30 min or until bacterial cells are completely lysed.
- 8) Transfer the lysate to a new 96-well dome plate.

III Loading the Kit to the Instrument

- 1) Add **400 µl** of **100 % Ethanol** to the 96-well dome plate containing lysate using multichannel pipettes.
- 2) Aliquot the solution from *MagListo™* Genomic DNA Extraction Kit to each of the new 96-well dome plate using multichannel pipette.

Cartridge No.	Solution	Volume
①	Lysate + Ethanol	Up to 830 µl
②	Magnetic Nano Bead solution	100 µl
③	WM1 Buffer	700 µl
④	WM1 Buffer	700 µl
⑤	W2 Buffer	700 µl
⑥	80 % Ethanol	700 µl
⑦	EA Buffer	100 µl

- 3) Press the 'Plate' Button on the instrument.
- 4) Place the Magnetic Rod Cover to the Magnetic Rod.
- 5) Place the plate onto the proper position of the base plate.
- 6) Press the 'Standard Protocol' Button and select '**K-3603/K-3615 Genomic DNA (V1.0)**'.
- 7) Press the 'Run' Button to start the selected protocol.