

Bioneer
Healthier Future for
Humanity with
Genomic Technology

Agaro-Power™

Convenient and Powerful Agarose Gel Electrophoresis Device



- The well size, gel length, gel thickness can be adjusted.
- Additional sealing tape is not needed.
- The tank lid is also specially designed.
- 96 loading adapters that fit multi-channeled pipette spacing are provided.
- The progress of the electrophoresis can be monitored due to transparent lid.

Agaro-Power™ is a convenient agarose gel electrophoresis device, which can load large amounts of samples in a short time. Agaro-Power™ as a molded unit has higher durability than those of the assembly type of the conventional electrophoresis devices. Furthermore, a power supply is also offered with no extra charge.

Technical Specifications

Distance between electrodes	228 mm
Buffer volume	700 ml
Gel Tray size	163(L) x 149(W) mm
Comb(No. of teeth)	25 well or 34 well/line x 1 ea 21 well or 34 well/line x 1 ea 15 well or 34 well/line x 1 ea 13 well or 34 well/line x 1 ea
Caster size	170(L) x 156(W) mm
Input Voltage	110 Vac 50/60 Hz, 220Vac 50/60 Hz
Output Voltage	High: 150V ±10% / Low: 75V ±10% (selectable)

Notice to Purchaser

Some products may not be available in some countries. For details, please contact your local distributor.

Ordering Information

Cat. No.	Product Description	Size
A-7020	Agaro-Power™ System	1 system
A-7020-1	Power Supply	1 ea
A-7020-2	Agaro Tank/Lid	1 ea
A-7020-3	Comb (25-well / 34-well) / (21-well / 34-well) / (15-well / 34-well) / (13-well / 34-well)	1 set
A-7020-4	Gel Caster/Tray	1 ea

Agaro-Power™ is a trademark of Bioneer Corporation. Copyright 2016. Bioneer Corporation. All Rights Reserved.