

ExiPrep™ 16 Plus Kit Series

Kit Description with Protocol Number and Runtime

2022.06.15
V.2 (EN)

Cat. No.	Instrument	Product Description		
A-5030	ExiPrep™ 16 Plus	Fully-automated instrument for nucleic acid extraction		
Cat. No.	Related Kit Description	Sample Sources	Protocol. No.	Run Time ± 2 min
For Extraction of Genomic DNA				
K-4211	ExiPrep™ Plus Blood Genomic DNA Kit	Whole blood	101, Genomic DNA: Whole Blood	1hr 50min
K-4214	ExiPrep™ Plus Bacteria Genomic DNA Kit	Gram (+) bacteria (~1x10 ⁹ cells)	108, Genomic DNA: Gram (+) bacteria	1hr 52min
		Gram (-) bacteria (~1x10 ⁹ cells)	109, Genomic DNA: Gram (-) bacteria	
		Sputum	115, Genomic DNA: Sputum	
		BAL	116, Genomic DNA: BAL	
		Saliva	117, Genomic DNA: Saliva	
		Swabs	118, Genomic DNA: Swab	
		Urine	119, Genomic DNA: Urine	
		Respiratory sample	125, Genomic DNA: Respiratory sample	
K-4215	ExiPrep™ Plus Plant Genomic DNA Kit	Plant tissue, seed (100 mg)	110, Genomic DNA: Yeast	1hr 30min
			111, Genomic DNA: Yeast	1hr 30min
K-4217	ExiPrep™ Plus Seed Genomic DNA Kit	Plant seed, tissue, root	105, Genomic DNA: Plant seed	1hr 17min
For Extraction of Total RNA				
K-4244	ExiPrep™ Plus Plant Total RNA Kit	Plant tissue	204, Total RNA: Plant tissue	1hr 20min
		Seed (100 mg)	205, Total RNA: Plant seed	
For Extraction of Viral DNA/RNA				
K-4271	ExiPrep™ Plus Viral DNA/RNA Kit	Plasma	612, Viral DNA/RNA: Plasma	1hr 43min
		Serum	613, Viral DNA/RNA: Serum	
		Swab	618, Viral DNA/RNA: Swab	
		Urine	619, Viral DNA/RNA: Urine	
		Cell free body fluid	621, Viral DNA/RNA: Cell free body fluid	
		Respiratory sample	625, Viral DNA/RNA: Respiratory sample	

Contact Us

Bioneer Corporation

8-11 Munpyeongseo-ro, Daedeok-gu
Daejeon, 34302, Republic of Korea
Tel: +82-42-930-8777 (Korea: 1588-9788)
Fax: +82-42-930-8688
E-mail: sales@bioneer.com

Bioneer Corporation

Bioneer Global Center, 71, Techono 2-ro,
Yuseong-gu, Daejeon, 34013, Republic of Korea
Tel: +82-42-939-6333
Fax: +82-42-939-6444
E-mail: sales@bioneer.com

Bioneer Inc.

155 Filbert St. Suite 216
Oakland, CA 94607, USA
Toll Free: +1-877-264-4300
Fax: +1-510-865-0350
E-mail: order.usa@bioneer.us.com

Bioneer R&D Center

Korea Bio Park BLDG #B-702
700 Daewangpangyo-ro, Bundang-gu, Seongnam-si
Gyeonggi-do, 13488, Republic of Korea
Tel: +82-31-628-0500
Fax: +82-31-628-0555